


Pg. 2 full body exerciser



Pg. 15 hand therapy



Pg. 61 electrotherapy



Pg. 23 / 25 measurement



Pg. 51 paraffin wax treatment



Pg. 32 sensory evaluation



Pg. 71 traction



Pg. 41 splinting materials



Pg. 53 hot and cold therapy



CANDo® Band

Preferred Colors® exercise band

the economic alternative to TheraBand® when TheraBand® colors are requested



- Made with Tear-Stop® technology for reduced puncture and tear “snap” breakage
- Powderless band leaves no residue
- Available in the Preferred Colors® sequence
- Textured for easy gripping

Perf 100® 100-yard dispenser with perforations

		Latex	Latex-Free
yellow	x-light	10-5191	10-5691
red	light	10-5192	10-5692
green	medium	10-5193	10-5693
blue	heavy	10-5194	10-5694
black	x-heavy	10-5195	10-5695
yel red grn blu blk	set of 5 (1 ea)	10-5198	10-5698
	set of 5 w/ Dispens-a-Band® rack	10-5199	10-5699

Twin-Pak® 100-yard box

		Latex	Latex-Free
tan	xx-light	10-5490	10-5650
yellow	x-light	10-5491	10-5651
red	light	10-5492	10-5652
green	medium	10-5493	10-5653
blue	heavy	10-5494	10-5654
black	x-heavy	10-5495	10-5655
silver	xx-heavy	10-5496	10-5656
gold	xxx-heavy	10-5497	10-5657
yel red grn blu blk	set of 5 (1 ea)	10-5498	10-5658

50-yard dispenser box

		Latex	Latex-Free
tan	xx-light	10-5220	10-5620
yellow	x-light	10-5221	10-5621
red	light	10-5222	10-5622
green	medium	10-5223	10-5623
blue	heavy	10-5224	10-5624
black	x-heavy	10-5225	10-5625
silver	xx-heavy	10-5226	10-5626
gold	xxx-heavy	10-5227	10-5627
yel red grn blu blk	set of 5 (1 ea)	10-5228	10-5628
	set of 5 w/ Dispens-a-Band® rack	10-5229	10-5629

25-yard dispenser box

		Latex	Latex-Free
tan	xx-light	10-5270	10-5630
yellow	x-light	10-5271	10-5631
red	light	10-5272	10-5632
green	medium	10-5273	10-5633
blue	heavy	10-5274	10-5634
black	x-heavy	10-5275	10-5635
silver	xx-heavy	10-5276	10-5636
gold	xxx-heavy	10-5277	10-5637
yel red grn blu blk	set of 5 (1 ea)	10-5278	10-5638
	set of 5 w/ Dispens-a-Band® rack	10-5279	10-5639

PERF 100

the better way to dispense bands

No scissors needed

Pull and separate at perforation to dispense

60 easy-to-dispense 5-foot exercisers in each box (same length each time)



Twin-Pak®
2 50-yard boxes



50-yards

Twin-Pak® offers 15% savings vs. 50-yard boxes



25-yards



6-yards

6-yard dispenser box

		Latex	Latex-Free
tan	xx-light	10-5210	10-5610
yellow	x-light	10-5211	10-5611
red	light	10-5212	10-5612
green	medium	10-5213	10-5613
blue	heavy	10-5214	10-5614
black	x-heavy	10-5215	10-5615
silver	xx-heavy	10-5216	10-5616
gold	xxx-heavy	10-5217	10-5617
yel red grn blu blk	set of 5 (1 ea)	10-5218	10-5618

CANDO[®] Singles[®] Preferred Colors[®] exercise band

polybagged for easy distribution

- Polybagged 5-foot Singles[®] pre-cuts are quick and easy to distribute
- Purchase in dispenser box or individual polybagged Singles[®]

(30) 5-foot Singles[®] in dispenser box

		Latex	Latex-Free
tan	xx-light	10-5180	10-5740
yellow	x-light	10-5181	10-5741
red	light	10-5182	10-5742
green	medium	10-5183	10-5743
blue	heavy	10-5184	10-5744
black	x-heavy	10-5185	10-5745
silver	xx-heavy	10-5186	10-5746
gold	xxx-heavy	10-5187	10-5747
yel red grn blu blk	set of 5 (1 ea)	10-5188	10-5748

5-foot Singles[®]

		Latex	Latex-Free
tan	xx-light	10-6450	10-5750
yellow	x-light	10-6451	10-5751
red	light	10-6452	10-5752
green	medium	10-6453	10-5753
blue	heavy	10-6454	10-5754
black	x-heavy	10-6455	10-5755
silver	xx-heavy	10-6456	10-5756
gold	xxx-heavy	10-6457	10-5757
yel red grn blu blk	set of 5 (1 ea)	10-6458	10-5758

(30) 5-foot Singles[®]
in dispenser box



PEP pack[®]

complete progressive rehab program

- 3 progressive 5-foot Singles[®] exercise bands
- Includes exercise chart
- Ideal as patient give-away or for retail sale



PEP[®] packs

PEP® packs				Latex	Latex-Free	
yellow	red	green	easy	set of 3 (1 ea)	10-5280	10-5680
green	blue	black	moderate	set of 3 (1 ea)	10-5282	10-5682
black	silver	gold	difficult	set of 3 (1 ea)	10-5284	10-5684

CANDO[®] tubing Preferred Colors[®] exercise tubing

100-foot dispenser box

		Latex	Latex-Free
tan	xx-light	10-5520	10-5720
yellow	x-light	10-5521	10-5721
red	light	10-5522	10-5722
green	medium	10-5523	10-5723
blue	heavy	10-5524	10-5724
black	x-heavy	10-5525	10-5725
silver	xx-heavy	10-5526	10-5726
gold	xxx-heavy	10-5527	10-5727
yel red grn blu blk	set of 5 (1 ea)	10-5529	10-5729

25-foot bag

		Latex	Latex-Free
tan	xx-light	10-5510	10-5710
yellow	x-light	10-5511	10-5711
red	light	10-5512	10-5712
green	medium	10-5513	10-5713
blue	heavy	10-5514	10-5714
black	x-heavy	10-5515	10-5715
silver	xx-heavy	10-5516	10-5716
gold	xxx-heavy	10-5517	10-5717
yel red grn blu blk	set of 5 (1 ea)	10-5519	10-5719

- Can be used in place of exercise band
- Use for hand therapy and athletic training
- The economic alternative when TheraBand[®] colors are requested
- Use alone or with handles and anchors



25-foot bag



100-foot dispenser box



- PEP[®] packs contain 3 progressive 6-foot exercise tubes
- Exercise chart included
- Ideal for patient give-aways or retail sale



PEP[®] pack

PEP[®] packs

			Latex	Latex-Free
yellow red green	easy	set of 3 (1 ea)	10-5380	10-5687
green blue black	moderate	set of 3 (1 ea)	10-5382	10-5688
black silver gold	difficult	set of 3 (1 ea)	10-5384	10-5689

Sup-R Band® and Sup-R Tubing® latex-free exercisers



Preferred Colors® latex-free exercise band

Performs like latex without latex allergenic properties

- Excels compared to other latex and latex-free bands
- Same force to stretch and contract as TheraBand® latex bands
- Available in the Preferred Colors® sequence
- Powderless for less mess; Textured for easy gripping
- Eliminates the need to stock both latex and latex-free bands
 - o Avoid risk of confusing latex and latex-free stock
 - o The only band your clinic will need
 - o Save money by stocking only Sup-R Band®



latex-free

Twin-Pak® 100-yard box

	xx-light	10-6330
	x-light	10-6331
	light	10-6332
	medium	10-6333
	heavy	10-6334
	x-heavy	10-6335
	xx-heavy	10-6336
	xxx-heavy	10-6337
	set of 5 (1 ea)	10-6338

50-yard dispenser box

	xx-light	10-6320
	x-light	10-6321
	light	10-6322
	medium	10-6323
	heavy	10-6324
	x-heavy	10-6325
	xx-heavy	10-6326
	xxx-heavy	10-6327
	set of 5 (1 ea)	10-6328
	set of 5 w/ Dispens-a-Band® rack	10-6329

6-yard dispenser box

	xx-light	10-6310
	x-light	10-6311
	light	10-6312
	medium	10-6313
	heavy	10-6314
	x-heavy	10-6315
	xx-heavy	10-6316
	xxx-heavy	10-6317
	set of 5 (1 ea)	10-6318

(30) 5-foot Singles® in dispenser box

Latex-Free

	x-light	10-6371
	light	10-6372
	medium	10-6373
	heavy	10-6374
	x-heavy	10-6375
	set of 5 (1 ea)	10-6378



(30) 5-foot Singles®
in dispenser box



Twin-Pak®
2 50-yard boxes

5-foot Singles®

Latex-Free

	x-light	10-6301
	light	10-6302
	medium	10-6303
	heavy	10-6304
	x-heavy	10-6305
	set of 5 (1 ea)	10-6308



50-yard roll

5-Foot Singles®



6-yard roll

PEP® packs

Latex-Free

			easy	set of 3 (1 ea)	10-6380
			moderate	set of 3 (1 ea)	10-6382
			difficult	set of 3 (1 ea)	10-6384



PEP® Pack



Preferred Colors® latex-free exercise tubing

- Performs like latex with the non-allergenic properties of latex-free
- Available in TheraBand® colors
- Can't confuse latex and latex-free tubing



100-foot dispenser box



latex-free



PEP® Pack



25-foot bag

PEP® packs

			easy	set of 3 (1 ea)	10-5880
			moderate	set of 3 (1 ea)	10-5882
			difficult	set of 3 (1 ea)	10-5884

100-foot dispenser box

Latex-Free

	xx-light	10-5860
	x-light	10-5861
	light	10-5862
	medium	10-5863
	heavy	10-5864
	x-heavy	10-5865
	xx-heavy	10-5866
	xxx-heavy	10-5867
	set of 5 (1 ea)	10-5869

25-foot bag


Latex-Free

	xx-light	10-5870
	x-light	10-5871
	light	10-5872
	medium	10-5873
	heavy	10-5874
	x-heavy	10-5875
	xx-heavy	10-5876
	xxx-heavy	10-5877
	set of 5 (1 ea)	10-5879



Val-u-Band® Berry Colors® exercise band

from the makers of CanDo® Band



- Save money by purchasing Val-u-Band® Berry Colors® band
- Save even more by purchasing Twin-Pak®
- Same great quality as CanDo® and TheraBand®
- Use with all dispenser systems
- Available in latex and latex-free formulations



Clinic Rolls



Twin-Pak®
(2) 50-yard boxes



50-yards



6-yards



50-yard set with
Dispens-a-Band® rack

50-yard dispenser box

			Latex	Latex-Free
pear	0	xx-light	10-6220	10-6120
peach	1	x-light	10-6221	10-6121
orange	2	light	10-6222	10-6122
lime	3	medium	10-6223	10-6123
blueberry	4	heavy	10-6224	10-6124
plum	5	x-heavy	10-6225	10-6125
mushroom	6	xx-heavy	10-6226	10-6126
pineapple*	7	xxx-heavy	10-6227	10-6127
1 2 3 4 5	set of 5 (1 ea)		10-6228	10-6128
	set w/ Dispens-a-Band® rack		10-6229	10-6129

*2 25-yard boxes (pineapple only)

**Twin-Pak® offers 15% savings
vs. 50-yard boxes**

Twin-Pak® 100-yard box

			Latex	Latex-Free
pear	0	xx-light	10-6230	10-6130
peach	1	x-light	10-6231	10-6131
orange	2	light	10-6232	10-6132
lime	3	medium	10-6233	10-6133
blueberry	4	heavy	10-6234	10-6134
plum	5	x-heavy	10-6235	10-6135
mushroom	6	xx-heavy	10-6236	10-6136
pineapple*	7	xxx-heavy	10-6237	10-6137
1 2 3 4 5	set of 5 (1 ea)		10-6238	10-6138

6-yard dispenser box

			Latex	Latex-Free
pear	0	xx-light	10-6210	10-6110
peach	1	x-light	10-6211	10-6111
orange	2	light	10-6212	10-6112
lime	3	medium	10-6213	10-6113
blueberry	4	heavy	10-6214	10-6114
plum	5	x-heavy	10-6215	10-6115
mushroom	6	xx-heavy	10-6216	10-6116
pineapple	7	xxx-heavy	10-6217	10-6117
1 2 3 4 5	set of 5 (1 ea)		10-6218	10-6118

Val-u-Band® Singles® Berry Colors® exercise band

from the makers of CanDo® Band

- Polybagged 5-foot Singles® pre-cuts are quick and easy to distribute
- Purchase polybagged Singles® in a dispenser box



(30) 5-foot Singles®
in dispenser box



5-foot Singles®

(30) 5-foot Singles® in dispenser box

			Latex	Latex-Free
pear	0	xx-light	10-6260	10-6160
peach	1	x-light	10-6261	10-6161
orange	2	light	10-6262	10-6162
lime	3	medium	10-6263	10-6163
blueberry	4	heavy	10-6264	10-6164
plum	5	x-heavy	10-6265	10-6165
mushroom	6	xx-heavy	10-6266	10-6166
pineapple	7	xxx-heavy	10-6267	10-6167
1 2 3 4 5	set of 5 (1 ea)		10-6268	10-6168


			Latex	Latex-Free
pear	0	xx-light	10-6270	10-6170
peach	1	x-light	10-6271	10-6171
orange	2	light	10-6272	10-6172
lime	3	medium	10-6273	10-6173
blueberry	4	heavy	10-6274	10-6174
plum	5	x-heavy	10-6275	10-6175
mushroom	6	xx-heavy	10-6276	10-6176
pineapple	7	xxx-heavy	10-6277	10-6177
1 2 3 4 5	set of 5 (1 ea)		10-6278	10-6178

REP Band® Berry Colors® latex-free exercise band

- REP Band® is only available in the latex-free formulation
- Only available in 4" width



latex-free



Twin-Pak®
(2) 50-yard boxes



6-yards



50-yards

Twin-Pak® 100-yard box

	1	2	3	4	5	set of 5 (1 ea)	10-1915
peach	1	x-light					10-1910
orange	2	light					10-1911
lime	3	medium					10-1912
blueberry	4	heavy					10-1913
plum	5	x-heavy					10-1914

Latex-Free

50-yard dispenser box

	1	2	3	4	5	set of 5 (1 ea)	10-1095
peach	1	x-light					10-1089
orange	2	light					10-1090
lime	3	medium					10-1091
blueberry	4	heavy					10-1092
plum	5	x-heavy					10-1093

Latex-Free

set of 5 (1 ea) w/ Dispens-a-Band® rack 10-5488

6-yard dispenser box

	1	2	3	4	5	set of 5 (1 ea)	10-1079
peach	1	x-light					10-1074
orange	2	light					10-1075
lime	3	medium					10-1076
blueberry	4	heavy					10-1077
plum	5	x-heavy					10-1078

Latex-Free

CanDo® exercise loops

- High quality latex band pre-formed into a loop
- Available in 3 lengths to perform upper, lower and full body exercises
- **Full Body Set** contains 10", 15" and 30" exercisers for total body workout




3 sizes

resistance of 10" loop is equivalent to 20" of flat band folded and tied together



upper body



lower body



whole body

			each			set
			10" lower	15" upper	30" whole body	full body set
yellow	x-light		10-5251	10-5261	10-5291	10-5841
red	light		10-5252	10-5262	10-5292	10-5842
green	medium		10-5253	10-5263	10-5293	10-5843
blue	heavy		10-5254	10-5264	10-5294	10-5844
black	x-heavy		10-5255	10-5265	10-5295	10-5845
silver	xx-heavy		10-5256	10-5266	10-5296	10-5846
gold	xxx-heavy		10-5257	10-5267	10-5297	10-5847
yel red grn blu blk	set of 5 (1 ea)		10-5259	10-5269	10-5299	10-5849

CanDo® tubing-with-handles

- 3 sizes to accommodate upper, lower and full body workouts
- Two durable foam-covered handles for comfort and firm grip
- **Full Body Set** contains 18", 36" and 48" exercisers for total body workout



*length of tubing from end to end when unit is laid flat


			each			set
			18 inch*	36 inch*	48 inch*	full body set
tan	xx-light		10-5540	10-5550	10-5560	10-5820
yellow	x-light		10-5541	10-5551	10-5561	10-5821
red	light		10-5542	10-5552	10-5562	10-5822
green	medium		10-5543	10-5553	10-5563	10-5823
blue	heavy		10-5544	10-5554	10-5564	10-5824
black	x-heavy		10-5545	10-5555	10-5565	10-5825
silver	xx-heavy		10-5546	10-5556	10-5566	10-5826
yel red grn blu blk	set of 5 (1 ea)		10-5549	10-5559	10-5569	10-5829

CanDo® tubing with cuffs exerciser

- Use on arms or legs for strength training, conditioning and rehabilitation
- Continuous tubing loop with two padded ankle/wrist cuffs
- Exerciser has 36" circumference



36" circumference



lower body



upper body

			18 inch*
tan	xx-light		10-5760
yellow	x-light		10-5761
red	light		10-5762
green	medium		10-5763
blue	heavy		10-5764
black	x-heavy		10-5765
silver	xx-heavy		10-5766
gold	xxx-heavy		10-5767
yel red grn blu blk	set of 5 (1 ea)		10-5769

*length of tubing from end to end when unit is laid flat


NOTE: Wear eye protection during the use of any band or tubing exercisers.

CanDo® accessory kit



- For all brands of exercise band and tubing, the kit includes accessories to enhance your workout

**Accessory kit includes:
A, B & C**



accessory kit

10-5334 accessory kit

standard foam handles with cinch

- Makes exercising with band and tubing easier
- Thread band / tubing through cinch
- Cinch makes adjusting band or tubing length easy
- 6" long foam-coated bar with heavy-duty webbing



10-5330	1 pair
10-5330-10	10 pair
10-5330-50	50 pair

double loop stirrup anchor

- Loop stirrup has two webbing loops: one large, one small
- Place band or tubing through small loop and put other loop around any stationary object such as your foot, a door anchor, or a wall anchor

10-5310	1 each
10-5310-10	10 each
10-5310-50	50 each



sports handle

- 12" foam-covered rigid handle
- 3 band / tubing cinch attachment positions
- Thread band / tubing through cinch to secure and adjust
- Use sports handle with band/tubing anchored to object or foot

10-5336	1 each
10-5336-10	10 each
10-5336-50	50 each



door jamb anchors

- Place flexible disc or nub anchor in door jamb at any height, even on top of door
- Close door to secure anchor
- Attach band or tubing to end of anchor
- Thread band / tubing through cinch to secure and adjust



nub anchor with loop

10-5320	1 each
10-5320-10	10 each
10-5320-50	50 each

disc anchor with cinch

10-5325	1 each
10-5325-10	10 each
10-5325-50	50 each

CanDo® Dispens-a-Band® rack

- Wood rack holds five 100-yard Perf 100® dispenser boxes and any brand 50-yard dispenser box
- Dispense directly from package



Dispens-a-Band® rack only

10-5111 rack only (holds 5 dispenser boxes)

Dispens-a-Band® rack with 5 band set

10-5229	50-yd CanDo® rolls (yellow-black)
10-5629	50-yd CanDo® latex-free rolls (yellow-black)
10-5199	100-yd Perf 100® rolls (yellow-black)
10-5699	100-yd Perf 100® latex-free rolls (yellow-black)
10-6329	50-yd Sup-R Band® latex-free rolls (yellow-black)

universal band racks

- Lightweight, stackable wood rack holds five 50-yard rolls of exercise band rolls
- Dispenses both core and coreless bands



universal band rack only

10-5104	single rack (holds 5 rolls)
10-5105	duplex rack (holds 10 rolls)
10-5106	scissors

universal band rack with 5 band set

10-5482	50-yd CanDo® rolls (yellow-black)
10-5483	50-yd CanDo® latex-free rolls (yellow-black)
10-5480	100-yd Perf 100® rolls (yellow-black)
10-5481	100-yd Perf 100® latex-free rolls (yellow-black)
10-5485	50-yd Sup-R Band® latex-free rolls (yellow-black)

CanDo® MGE™ exerciser

- Easy-to-use grips allow for unlimited exercise options
- Each exerciser contains both large (6") and small (2") grips. Use small grips as anchor or for finger exercises
- Cut exerciser from roll into any length (number of grips) desired
- Wrap around hands and feet to perform high impact upper and lower body exercises
- **Cloth exerciser contains no latex, no scent, and no powder coating**
- Can be machine washed and dried; use with any antiseptic spray to disinfect exerciser
- Easily change resistance by moving 1 grip up or down the exerciser, or by using a different resistance
- Progressively exercise with 8 different resistance levels
- Available in 9-grip (6 ft. exercisers) and 135-grip (30 yd. rolls) 15 exercisers



		9 grips, 6 foot exerciser	9 grips, 6 foot exerciser (24 ea)	135 grips, 30 yard roll, 15 exercisers
tan	xx-light	10-6600	10-6600-24	10-6610
yellow	x-light	10-6601	10-6601-24	10-6611
red	light	10-6602	10-6602-24	10-6612
green	medium	10-6603	10-6603-24	10-6613
blue	heavy	10-6604	10-6604-24	10-6614
black	x-heavy	10-6605	10-6605-24	10-6615
silver	xx-heavy	10-6606	10-6606-24	10-6616
gold	xxx-heavy	10-6607	10-6607-24	10-6617



10-1171 Protex™ disinfectant 32 oz trigger spray

9 grip (6 foot) exercisers

FabStretch™ 4-level incline board

- 4 incline levels: 5°, 15°, 25° and 35°
- Stretch calves and upper and lower back at the same time
- Made of composite plastic with 500 lb capacity
- Folds flat (3") to store easily
- Non-slip top is 14" x 14"



10-1179 incline board

CanDo® overdoor shoulder exerciser, single and double pulley

- Increase and maintain ROM on all planes of shoulder movement
- Rope adjusts for any height user; stopper prevents over-extension
- Foam-covered handles provide a comfortable grip
- **Visualizer®**: color-coded rope intervals reference ROM progress
- **Double pulley** has 2 pulleys and door bracket



Single pulley
w/ nub anchor (ea)
w/ nub anchor (25 ea)
w/ door bracket (ea)
w/ door bracket (25 ea)
w/ disc anchor (ea)
w/ disc anchor (25 ea)
Double pulley
double pulley (ea)
double pulley (25 ea)




Standard
50-1023
50-1023-25
50-1024
50-1024-25
50-1025
50-1025-25
Standard
50-1010
50-1010-25



Visualizer®
50-1035
50-1035-25
50-1036
50-1036-25
50-1037
50-1037-25
Visualizer®
50-1034
50-1034-25

CanDo® stretch strap

- Dynamic stretch strap has 8 elastic pockets to allow the user to ease into a stretch
- Designed for "facilitated stretching", an active form of stretching that uses isometric contractions to achieve greater flexibility gains than from static stretching



10-1385 single stretch strap
10-1385-25 25 stretch straps

CanDo® back stretcher



- Regular use can improve posture and relieve discomforts caused by poor posture, muscle spasms, headaches, and muscle tension

10-7900 back stretcher

CanDo® leg stretcher



- Use to reduce the risk of injury before and after exercise
- Increases flexibility, range of motion, and performance
- Gives support for controlled stretching

10-1178 single leg stretch

CanDo® Magneciser® pedal exerciser

- Quiet, magnetic, bi-directional operation with variable resistance (turn knob to change resistance)
- 5-function digital display

01-8030 pedal exerciser



CanDo® deluxe pedal exerciser

- Bi-directional exerciser can be used for upper- or lower-body pedaling
- Resistance can be adjusted with one easy-to-use dial
- Fully automated electronic monitor with LED display that shows: speed, distance (per session), timer, odometer (cumulative distance), calories burned, scan mode

10-0717 deluxe pedal exerciser



CanDo® passive pedal exerciser

- Bi-directional exerciser offers active and passive (assisted) motion for upper and lower body exercises
- Resistance / speed can be adjusted with wireless remote to assist patients
- 10 pre-set assisted speed / resistance programs
- Large LCD screen

10-2890 passive pedal exerciser



Magneciser® shoulder wheel

- Magnetic adjustable resistance for smooth and quiet operation
- Increase range-of-motion, muscle strength, mobility and coordination
- Bi-directional arm rotation
- Arm length adjusts from 13" (33 cm) to 22" (56 cm); equivalent to 26" - 44" diameter shoulder wheel
- Height-adjustable wall mount for standing, sitting or for wheelchair access
- 1¾" away from wall

01-8025 magnetic shoulder exerciser



CanDo® pedal exercisers

- Portable bi-directional upper- and lower-body pedaling
- Wide leg spread and rubber non-marring feet keep exerciser stable
- Use wraparound adjustable straps to secure feet or arms
- Adjust resistance by turning tension knob

10-0710 pedal exerciser

10-0711 knocked down, 3 piece



CanDo® deluxe folding pedal exercisers

- Folds for easy storage and transportation
- Use for upper- and lower-body pedaling exercises
- Adjustable resistance
- Durable metal frame
- Bi-directional
- Non-skid feet
- Adjustable hand/foot straps
- **Digital model** keeps track of exercise performed: time, calories burned, distance traveled and scan mode

10-0718 folding pedal exerciser

10-0712 folding pedal exerciser w/ digital readout



digital folding
pedal exerciser



folded

digital display

PhysioTrainer bi-directional upper body exerciser

- Heavy duty arm and leg trainer is designed for clinic or home use
- Use on table top or floor
- Easy-to-read, self-powered display monitors progress
- Adjust resistance with tension control knob
- Use with available manual or motorized table to reach optimal height for the user



tension control



shown with
motorized
table

69-0131 PhysioTrainer Bi-Directional UBE

69-0134 PhysioTrainer w/ table

69-0135 PhysioTrainer w/ motorized table

69-0132 PhysioTrainer Adjustable UBE Table

69-0133 Motorized Adjustable UBE Table



display

the Cuff® weights

- Flexible weight contours to snugly fit the ankle, wrist... and even thigh
- Color-coded weight levels
- Marked with both lbs and kgs
- Sold individually
- the Cuff® deluxe weight** has an extra-long loop strap with extra hook tabs
- Mobile rack** holds up to 500 lbs



10-0150 mobile weight rack

the Cuff® original			
¼ lb	white	10-0200	
½ lb	walnut	10-0201	
¾ lb	orange	10-0202	
1 lb	blue	10-0203	
1½ lb	olive	10-0204	
2 lb	white	10-0205	
2½ lb	red	10-0206	
3 lb	gold	10-0207	
4 lb	turquoise	10-0208	
5 lb	black	10-0209	
6 lb	beige	10-0210	
7 lb	lemon	10-0211	
7½ lb	orange	10-0212	
8 lb	red	10-0213	
9 lb	parchment	10-0214	
10 lb	brown	10-0215	
12½ lb	olive	10-0216	
15 lb	tan	10-0217	
20 lb	blue	10-0218	
25 lb	green	10-0219	

the Cuff® deluxe

10-2500
10-2501
10-2502
10-2503
10-2504
10-2505
10-2506
10-2507
10-2508
10-2509
10-2510
10-2511
10-2512
10-2513
10-2514
10-2515
10-2516
10-2517
10-2518
10-2519

the Cuff® original

Functional set (7 pieces)

1 ea: 1, 2, 3, 4, 5, 7½, 10 lb		
functional set	10-0250	10-2550
set with rack	10-0251	10-2551

Economy set (20 pieces)

2 ea: ¼, ½, ¾, 1, 1½, 2, 2½, 3, 4, 5 lb		
economy set	10-0252	10-2552
set with rack	10-0253	10-2553

Standard set (24 pieces)

2 ea: ¼, ½, ¾, 1, 1½, 2, 2½, 3, 4, 5, 7½, 10 lb		
standard set	10-0254	10-2554
set with rack	10-0255	10-2555

Deluxe set (32 pieces)

2 ea: ¼, ½, ¾, 1, 1½, 2, 2½, 3, 4, 5, 6, 7, 7½, 8, 9, 10 lb		
deluxe set	10-0256	10-2556
set with rack	10-0257	10-2557

Knee-Hip-Leg set (8 pieces)

2 ea: 10, 12½, 15, 20 lb		
knee-hip-leg set	10-0258	10-2558

the Cuff® variable weights

- One weight serves the purpose of many
- Add or remove weight inserts to change weight
- Available for use on ankle or wrist; for adults and children



		each weight	inserts	pairs
adult	wrist	2 lb	10 x 0.2 lb	10-3340-2
	wrist	4 lb	20 x 0.2 lb	10-3341-2
	ankle	5 lb	10 x 0.5 lb	10-3331-2
	ankle	10 lb	20 x 0.5 lb	10-3332-2
	ankle	20 lb	20 x 1.0 lb	10-3333-2
children	wrist	2 lb	12 x 0.17 lb	10-3345-2
	ankle	2 lb	12 x 0.17 lb	10-3335-2

CanDo® weight straps

- Soft and comfortable padded weights won't slip off your wrist or ankle
- Hook/loop strap with D-ring closure for snug fit
- Sold in pairs



	weight (each)		weight (pair)	
10-3380-2	½ lb	(0.21 kg)	2 ea	1 lb (0.45 kg)
10-3381-2	1 lb	(0.45 kg)	2 ea	2 lb (0.91 kg)
10-3382-2	1½ lb	(0.68 kg)	2 ea	3 lb (1.36 kg)
10-3383-2	2 lb	(0.91 kg)	2 ea	4 lb (1.81 kg)
10-3384-2	2½ lb	(1.13 kg)	2 ea	5 lb (2.27 kg)

CanDo® SoftGrip® hand weights

- Not intimidating to use
- Easy to grasp
- Comfortable and soft, but not floppy



	weight (lb / kg)	each	pair
tan	0.5 / 0.23	10-0350-1	10-0350-2
yellow	1 / 0.45	10-0351-1	10-0351-2
red	1.5 / 0.68	10-0352-1	10-0352-2
green	2 / 0.91	10-0353-1	10-0353-2
blue	2.5 / 1.13	10-0354-1	10-0354-2
black	3 / 1.36	10-0355-1	10-0355-2
silver	4 / 1.81	10-0356-1	10-0356-2
gold	5 / 2.27	10-0357-1	10-0357-2
tan	6 / 2.72	10-0358-1	10-0358-2
yellow	7 / 3.18	10-0359-1	10-0359-2
red	7.5 / 3.40	10-0360-1	10-0360-2
green	8 / 3.63	10-0361-1	10-0361-2
blue	9 / 4.08	10-0362-1	10-0362-2
black	10 / 4.54	10-0363-1	10-0363-2

CanDo® WaTE™ hand balls

- Small and flexible enough to grasp with one hand, yet big enough to hold with two hands
- Each ball is 5" (12.7 cm) diameter



5" diameter	weight (lb / kg)	each
tan	1.1 / 0.5	10-3160
yellow	2.2 / 1.0	10-3161
red	3.3 / 1.5	10-3162
green	4.4 / 2.0	10-3163
blue	5.5 / 2.5	10-3164
black	6.6 / 3.0	10-3165
tan, yellow, red, green, blue, black	set of 6 (1 ea)	10-3166
tan, yellow, red, green, blue, black	set of 6 w/ 1-tier rack	10-3189
tan, yellow, red, green, blue, black	set of 6 w/ 2-tier rack	10-3188

CanDo® vinyl-coated iron dumbbells

- Vinyl-coated, color-coded dumbbells feel comfortable and will not scratch floor
- economy set (10 pieces)**
2 ea: 1 through 5 lb
- 10-0563 economy
10-0564 w/ wall rack
10-0578 w/ floor rack
10-0575 wall rack only



	lbs	kgs	each	pair
pink	1	0.45	10-0550-1	10-0550-2
purple	2	0.91	10-0551-1	10-0551-2
green	3	1.36	10-0552-1	10-0552-2
light blue	4	1.81	10-0553-1	10-0553-2
blue	5	2.27	10-0554-1	10-0554-2
red	6	2.72	10-0555-1	10-0555-2
purple	7	3.18	10-0556-1	10-0556-2
black	8	3.63	10-0557-1	10-0557-2
yellow	9	4.08	10-0558-1	10-0558-2
orange	10	4.54	10-0559-1	10-0559-2
silver	15	6.80	10-0560-1	10-0560-2
brown	20	9.07	10-0561-1	10-0561-2



wall rack

CanDo® vinyl-coated iron kettlebells


- For rehabilitation, strength training and muscle toning
- Color-coded for easy weight identification
- Vinyl-coating gives a comfortable feel and reduces potential floor scratching



	weight (lb / kg)	diameter (inch / cm)	each
yellow	5 / 2.3	8 / 20.3	10-3191
red	7.5 / 3.4	8 / 20.3	10-3192
green	10 / 4.5	9 / 22.9	10-3193
blue	15 / 6.8	9 / 22.9	10-3194
black	20 / 9.1	10 / 26.4	10-3195
silver	25 / 11.4	11 / 27.9	10-3196
gold	30 / 13.6	11 / 27.9	10-3197
yel red grn blu blk	set of 5, 1 ea (5, 7, 10, 15, 20 lb)		10-3199

CanDo® WaTE™ exercise bars

- Easy-to-grip and store color-coded exercise WaTE™ bars from 1 to 25 lb; exercise instructions included
- Jumbo® WaTE™ bars are sleek, foam-covered with color-coded end caps



Jumbo®
WaTE™
bar




Slim®
WaTE™ bars

Slim® WaTE™ bar (1.5" diameter)

	lb	kg	L (in)	
tan	1	0.45	36"	10-1600
tan stripe	1½	0.68	36"	10-1601
yellow	2	0.91	36"	10-1602
yellow stripe	2½	1.14	36"	10-1603
red	3	1.36	36"	10-1604
red stripe	3½	1.59	36"	10-1605
green	4	1.82	36"	10-1606
green stripe	4½	2.05	36"	10-1607
blue	5	2.27	36"	10-1608
blue stripe	6	2.73	36"	10-1609
black	7	3.18	36"	10-1610
black stripe	7½	3.41	36"	10-1623
silver	8	3.64	36"	10-1611
silver stripe	9	4.09	36"	10-1612
gold	10	4.55	36"	10-1613

Jumbo® WaTE™ bar (2" diameter)

	lb	kg	L (in)	
tan	2½	1.14	48"	10-1636
yellow	5	2.27	48"	10-1637
red	7½	3.41	48"	10-1638
lt green	10	4.55	48"	10-1639
dk green	12½	5.68	48"	10-1614
lt blue	15	6.82	48"	10-1615
dk blue	17½	7.95	48"	10-1616
black	20	9.09	48"	10-1617
silver	22½	10.23	55"	10-1618
gold	25	11.36	55"	10-1619



Jumbo® WaTE™ bars

Slim® Introductory set

10-1624 set of 3,
(1 ea: 1, 2, 3 lb)

Slim® LiTE® set

10-1620 set of 5,
(1 ea: 1, 2, 3, 4, 5 lb)

Mini™ WaTE™ exercise bar – pairs

- Easy-to-grip, easy-to-store, color-coded weighted bars
- Easier to use than traditional weight bars
- 14" long, 1" diameter



	pair
yellow	10-1650 (2) 1 lb
red	10-1651 (2) 1½ lb
green	10-1652 (2) 2 lb
blue	10-1653 (2) 2½ lb
black	10-1654 (2) 3 lb

CanDo® Handy Grip™ balls

- Ergonomic shape, texture and strap make gripping easier
- Enables a variety of strength training activities and functional movements
- Secure ball to hand with adjustable strap for firm grasp



	weight (lb / kg)	each
 tan	1 / 0.45	10-0440
 yellow	2 / 0.91	10-0441
 red	3 / 1.36	10-0442
 green	4 / 1.81	10-0443
 blue	5 / 2.27	10-0444
 set of 5 (1 ea)		10-0445



CanDo® rubber medicine balls

- Textured surface for easy grip
- Color-coded for easy identification
- Constructed with rubber to bounce



	weight		diameter		
	lbs	kgs	in	cm	
tan	1	0.45	8	20.3	10-3140
yellow	2	0.90	8	20.3	10-3141
red	4	1.80	8	20.3	10-3142
green	7	3.20	9	22.9	10-3143
blue	11	5.00	9	22.9	10-3144
black	15	6.80	10	25.4	10-3145
silver	20	9.10	11	27.9	10-3147
gold	30	13.60	27.9	11.0	10-3148
yel red grn blu blk	set (2, 4, 7, 11, 15 lb)				10-3146

CanDo® P.T. pliable medicine balls

- Textured tactile surface is easy to grip and catch
- Increase/decrease air in ball to enhance grasp
- Use for upper body rehabilitation and fitness exercises



	weight		diameter		
	lbs	kgs	in	cm	
tan	1	0.45	5	12.7	10-3170
yellow	2	0.90	5	12.7	10-3171
red	4	1.80	5	12.7	10-3172
green	7	3.20	7	17.8	10-3173
blue	11	5.00	7	17.8	10-3174
black	15	6.80	9	22.9	10-3175
silver	20	9.10	9	22.9	10-3177
gold	30	13.60	9	22.9	10-3178
yel red grn blu blk	set (2, 4, 7, 11, 15 lb)				10-3176

CanDo® round rebounder

- Returns ball with speed, force and angle consistent with throw
- Adjustable incline angle: 0, 28, 33, 39 & 42 degrees
- Built-in floor rack to store balls
- Minimal assembly required
- Can ship by UPS or FedEx
- 40 lbs; 38½" L x 36" W x 42" H
- 27" diameter surface



rebounder

10-3112 round rebounder

rebounder sets

- 10-3136 with 5 P.T. balls: 2, 4, 7, 11, 15 lb
- 10-3135 with balls and horizontal plastic rack (10-3168)
- 10-3137 with balls and vertical metal rack (10-3120)

CanDo® square rebounder

- Returns ball with speed, force and angle consistent with throw
- Adjustable incline angle adjust: 0, 28, 33, 39 & 42 degrees
- Built-in floor rack to store balls
- Minimal assembly required
- Can ship by UPS or FedEx
- 59 lbs; 39½" L x 39½" W x 52½" H
- 39" diagonal surface (27.5" x 27.5")



rebounder with
5 PT balls

rebounder

10-3110 square rebounder

rebounder sets

- 10-3132 with 5 P.T. balls: 2, 4, 7, 11, 15 lb
- 10-3131 with balls and horizontal plastic rack (10-3168)
- 10-3133 with balls and vertical metal rack (10-3120)

CanDo® weight & mirror rack

- Easily store up to 6 band or tubing rolls, 22 dumbbells and 48 the Cuff® weights
- Mirror features ANSI safety backing

Accessories include:

- 16 the Cuff® weights (2 ea: 1, 1.5, 2, 2.5, 3, 4 & 5 lbs; 1 ea: 7.5, 10 lbs)
- 20 dumbbells (2 ea: 1 - 10 lbs)
- 5 - 50 yard exercise band rolls (yellow - black)

- 15-4255 with accessories
- 15-4234 without accessories



CanDo® Rickshaw™ rehab exerciser

- Designed for people who use a wheelchair
- Use to strengthen arms and shoulder muscles
- Use facing either forwards or backwards
- Arms can be used individually or in tandem by adding crosspiece
- Heavy-duty steel construction
- Free weights sold separately



10-0690 CanDo® Rickshaw™ Rehab Exerciser

CanDo® graded pinch exerciser

- Squeeze pinchpin to improve pinch strength and coordination
- Place pinchpin on one of 3 different diameter horizontal bars or the vertical bar
- Complete set includes 5 sets (10 ea) of color-coded graded resistance pinchpins, three horizontal rods, and a vertical rod for shoulder ROM exercise



complete set

10-0840 pinch exerciser with 50 pins

replacement pins

		force to pinch
yellow	10-0841-10	10 ea 1 lb
red	10-0842-10	10 ea 2 lb
green	10-0843-10	10 ea 4 lb
blue	10-0844-10	10 ea 6 lb
black	10-0845-10	10 ea 8 lb
set of 50 (10 ea)	10-0846-50	

CanDo® multi-function hand & wrist exerciser

- Improve flexion/extension, pronation/supination, and ulnar/radial deviation
- Mix and match different sized semi-spheres and handles
- Set includes:
3 semi-spheres and 3 handles



use with handle



use semi-sphere alone



set

10-2995 3 semi-spheres, 3 handles

individual parts

10-2900	ball only, x-small, black
10-2902	ball only, medium, green
10-2904	ball only, x-large, yellow
10-2905	handle only, x-small, black
10-2907	handle only, medium, green
10-2909	handle only, x-large, yellow

CanDo® wrist wheel



supination



pronation

- Use for pronation/supination exercises
- Rubberized outside surface provides a secure, non-slip surface during exercise
- Hand bar is large enough to accommodate any hand size
- Off center hand bar gives two travel distances

10-2960 CanDo wrist wheel

CanDo® web - flexion, extension and wrist exerciser

- Perform flexion, extension, opposition and supination exercises
- Modify exercises by adjusting the hand position or depth of finger insertion
- Modify resistance by switching to a different resistance color web
- Available in 14" and 7" dia. sizes; latex and latex-free



tan	xx-light	10-0850
yellow	x-light	10-0851
red	light	10-0852
green	medium	10-0853
blue	heavy	10-0854
black	x-heavy	10-0855
tan yel red grn blu blk	set of 6 (1 ea)	10-0856

large 14" diameter

latex	latex-free
10-0850	10-0870
10-0851	10-0871
10-0852	10-0872
10-0853	10-0873
10-0854	10-0874
10-0855	10-0875
10-0856	10-0876

small 7" diameter

latex	latex-free
10-0860	10-0880
10-0861	10-0881
10-0862	10-0882
10-0863	10-0883
10-0864	10-0884
10-0865	10-0885
10-0866	10-0886

skate and powder board

- Skates may be used on floor, table or powder board
- Powder board with laminated walnut top and folding legs



10-1130	6" x 6" skate
10-1131	6" x 12" skate
10-1120	powder board

finger / shoulder ladder


- Offers progressive mobility of arm at shoulder, elbow and wrist joints
- Solid wood units have 36 steps at 1¼" intervals





10-1160 finger / shoulder ladder

shoulder wheel

- Range-of-motion arc can be varied from 10" to 39" by moving handle
- Adjustable tension
- Adjustable height unit accommodates seated and standing therapy
- 37½" diameter



10-1150 shoulder wheel



TheraPutty®

RESISTIVE HAND EXERCISE MATERIAL

- Each color-coded putty has a different firmness ranging from xx-soft to x-firm
- TheraPutty® is gluten, casein, and latex-free; non-toxic if swallowed
- TheraPutty® material does not dry out
- TheraPutty® Plus** contains non-toxic anti-microbial agents to kill micro-organisms
- Puff® Putty** has a light and airy feel vs. traditional TheraPutty® material

Extra Containers with lids



	qty	
2/3 oz	25	10-0940
4/6 oz	25	10-0941
1 lb	10	10-0943
5 lb	10	10-0944

TheraPutty® hand exercise material

		tan	yellow	red	green	blue	black
		xx-soft	x-soft	soft	medium	firm	x-firm
patient	2 oz	10-0956	10-0900	10-0901	10-0902	10-0903	10-1466
	3 oz	10-0957	10-0967	10-0968	10-0969	10-0970	10-1467
	4 oz	10-0958	10-0905	10-0906	10-0907	10-0908	10-1468
	6 oz	10-0959	10-0910	10-0911	10-0912	10-0913	10-1469
clinic	1 lb	10-0994	10-0918	10-0919	10-0920	10-0921	10-1470
	5 lb	10-0995	10-0923	10-0924	10-0925	10-0926	10-1471
	50 lb	10-1460	10-1461	10-1462	10-1463	10-1464	10-1473

set of 6 (1 ea)

		tan	yel	red	grn	blu	blk
2 oz							
3 oz							
4 oz							
6 oz							
1 lb							
5 lb							

TheraPutty® Plus Anti-microbial hand exercise material

		tan	yellow	red	green	blue	black
		xx-soft	x-soft	soft	medium	firm	x-firm
patient	2 oz	10-2600	10-2601	10-2602	10-2603	10-2604	10-2605
	3 oz	10-2610	10-2611	10-2612	10-2613	10-2614	10-2615
	4 oz	10-2620	10-2621	10-2622	10-2623	10-2624	10-2625
	6 oz	10-2630	10-2631	10-2632	10-2633	10-2634	10-2635
clinic	1 lb	10-2640	10-2641	10-2642	10-2643	10-2644	10-2645
	5 lb	10-2650	10-2651	10-2652	10-2653	10-2654	10-2655
	50 lb	10-2660	10-2661	10-2662	10-2663	10-2664	10-2665

set of 6 (1 ea)

		tan	yel	red	grn	blu	blk
2 oz							
3 oz							
4 oz							
6 oz							
1 lb							
5 lb							

Puff® Putty hand exercise material

		tan	yellow	red	green	blue	black
		xx-soft	x-soft	soft	medium	firm	x-firm
cc	2 oz	10-1400	10-1401	10-1402	10-1403	10-1404	10-1405
	3 oz	10-1410	10-1411	10-1412	10-1413	10-1414	10-1415
	4 oz	10-1420	10-1421	10-1422	10-1423	10-1424	10-1425
	13 oz	10-1430	10-1431	10-1432	10-1433	10-1434	10-1435
	3.5 lb	10-1440	10-1441	10-1442	10-1443	10-1444	10-1445

set of 6 (1 ea)

		tan	yel	red	grn	blu	blk
2 oz							
3 oz							
4 oz							
13 oz							
3.5 lb							

puttycise® exercise tools

- Tools used to simulate functional activity
- Push, pull, or turn tools through putty for specific exercises
- Change putty resistance to increase or decrease difficulty
- Use with 6 oz or more of TheraPutty® material
- Teaches joint protection techniques

Puttycise® exercise kits

- TheraPutty® exercise material containers (1 each of 4 resistances as indicated below)
- Puttycise® tools (Set of 5, 1 each)
- Carrying bag and instructions


	easy	moderate	difficult
size	easy	moderate	difficult
6 oz	10-2822	10-2828	10-2834
1 lb	10-2824	10-2830	10-2836
5 lb	10-2826	10-2832	10-2838

Puttycise® tools

10-2821	tool set w/ carrying bag and instructions
10-2820	5-piece tool set
10-2810	knob turn tool
10-2811	L-bar tool
10-2812	peg turn tool
10-2813	key turn tool
10-2814	cap turn tool
10-2819	carrying bag
10-2818	instructions
10-2840	putty work board



CanDo® variable strength putty



- One putty with variable resistance
- Work pliable pellets into base putty to increase resistance as your grip strength increases
- Pellet pack contains 4 pellets

sets

10-2670	2½ oz base putty, 1 pack
10-2671	1 lb base putty, 6 packs
10-2672	5 lb base putty, 32 packs

base putty only

10-2680	2½ oz base putty
10-2681	1 lb base putty
10-2682	5 lb base putty

pellet packs (4 pellets per pack)

10-2690	1 pellet pack
---------	---------------

Digi-Flex Multi® hand exerciser

Customize hand therapy like never before. Adjust each finger's resistance to fit the patient's rehab program.

Clinic sets are available with 5 bases or 5 pre-built exercisers, 20 or 32 piece resistive button sets, storage case and rack.

More exercise bases, exercisers and resistive buttons are available.



pinch release tabs and pull up to remove button



insert button until a "snap" is heard



use your customized hand exerciser



small clinic set with 5 bases, deluxe (10-3840)



large clinic set with 5 pre-built exercisers, standard (10-3846)

patented

buttons only with storage box

10-3759 32 buttons 4 ea tan yel red grn blu blk sil gld

10-3758 20 buttons 4 ea yel red grn blu blk

small clinic set with 5 bases

10-3840 deluxe, 32 buttons 4 ea tan yel red grn blu blk sil gld

10-3841 standard, 20 buttons 4 ea yel red grn blu blk

large clinic set with 5 pre-built exercisers

10-3845 deluxe, 32 buttons 4 ea tan yel red grn blu blk sil gld

10-3846 standard, 20 buttons 4 ea yel red grn blu blk

individual finger lock buttons

Insert the white finger lock button when finger movement or exertion is counter-indicated.

10-3738-4 4 each with case



additional bases


10-3830 1 base

10-3832 5 bases



Digi-Flex® hand exercisers

- Use to rehabilitate injuries or conditions affecting the hand and upper extremities
- Develops isolated finger strength, flexibility, and coordination
- Develop whole hand strength
- **Digi-Flex®** has upper and lower spring to allow for more hand motion flexibility. Use for individual finger or hand
- **Digi-Flex LiTE®** is an inexpensive option for patient giveaway or sale. It has individual finger locks and a single spring design
- **Digi-Flex Thumb®** allows for thumb involvement. The thumb piece can attach at three positions to accommodate all thumb sizes and swivels for all thumb positions



exercise isolated finger by placing the palm bar over the continuous base

Digi-Flex®

Digi-Flex LiTE®

Digi-Flex Thumb®

patented

patented


force to compress (lb)
finger hand

Digi-Flex®

Digi-Flex LiTE®

Digi-Flex Thumb®

		finger	hand	
	tan	xx-light	0.75	2.5
	yellow	x-light	1.50	5.0
	red	light	3.00	10.0
	green	medium	5.00	16.0
	blue	heavy	7.00	23.0
	black	x-heavy	9.00	31.0
	silver	xx-heavy	11.00	36.0
	gold	xxx-heavy	13.00	43.0
yel red grn blu blk	set of 5 (1 each)			10-0745
yel red grn blu blk	set with plastic stand			10-0746
yel red grn blu blk	set with metal stand			10-0749
tan yel red grn blu blk sil gld	set of 8 (1 ea)			10-0758
tan yel red grn blu blk sil gld	set of 8 (1 ea)			10-0759
w/ 2 metal stands				



metal stand

Digi-Flex® Stands

- Holds Digi-Flex® units
- 10-0750 metal stand
10-0747 plastic stand

CanDo® gel hand exercise balls

- Hand exerciser for grip strength, dexterity, mobility and fine and gross motor skills
- Heat or chill for hot or cold therapy



standard (round) 2" dia.



large ergonomic (egg-shaped)
2.15" dia. on long dimension
1.75" dia. on short dimension



	xx-light	10-1490
	x-light	10-1491
	light	10-1492
	medium	10-1493
	heavy	10-1494
	x-heavy	10-1495
	set of 6 (1 ea)	10-1496

standard (round)

10-1490
10-1491
10-1492
10-1493
10-1494
10-1495
10-1496

large ergonomic (egg-shaped)

10-1890
10-1891
10-1892
10-1893
10-1894
10-1895
10-1896


CanDo® fixed resistance grips

- Resistance levels range from earliest rehabilitation after a stroke (yellow) to advanced athletic training (black)
- Force to squeeze grip doubles with each level
- Sold in pairs



		resistance (lbs)	
	x-easy	3	10-1804
	easy	6	10-1805
	moderate	12	10-1806
	difficult	24	10-1807
	x-difficult	48	10-1808
	set of 5 (1 ea)		10-1809

power hand grip



- Adjustable grippers
- Interchange springs or spring position to change tension

10-0815	pro (5-150 lb)
10-0813	deluxe (10-100 lb)
10-0816	adjustable (5-50 lb)

CanDo® rubber band hand exerciser

- Vary resistance by increasing or decreasing number of rubber bands, or by changing to a different color-coded set of bands
- Comes with 5 red (light) bands or with 5 sets of color-coded bands. Set includes 25 bands – 5 of each color: tan (xx-light), yellow (x-light), red (light), green (medium) and blue (heavy)



Exercisers

with 5 red (light) bands	
1 set	10-1860
10 sets	10-1860-10
50 sets	10-1860-50

with 25 bands (5 each color / resistance)	
1 set	10-1865
10 sets	10-1865-10
50 sets	10-1865-50



Additional Bands

25 additional resistance bands

	xx-light	10-1850
	x-light	10-1851
	light	10-1852
	medium	10-1853
	heavy	10-1854

Set of 25 (5 each color)

	10-1855
--	---------

CanDo® flexion gloves

- Increase flexion of MP, PIP and DIP joints
- Made of soft stretch nylon with eyelets at each finger tip
- Standard glove has one wrist hook, deluxe gloves have four Velcro wrist tabs
- Deluxe glove can be purchased with or without a thumb fixture
- Measure width at MP joint for sizing: small / medium 2.5" - 3.5"; large / x-large 3.5" - 4.5"



standard



deluxe




deluxe with thumb

		left	right
standard	small / medium	10-4000L	10-4000R
	large / x-large	10-4001L	10-4001R
deluxe	small / medium	10-4002L	10-4002R
	large / x-large	10-4003L	10-4003R
deluxe with thumb	small / medium	10-4004L	10-4004R
	large / x-large	10-4005L	10-4005R


CanDo® Twist-n' Bend® exercisers

- Change resistance by varying distance between hands and/or by using a different resistance bar
- Each resistance has the same bar diameter 1 3/4"
- Use the bar to Twist, Bend, Shake, or Oscillate



12" bend

36" shake



24" twist

	xx-light	
	x-light	
	light	
	medium	
	heavy	
	x-heavy	
	set of 6 (1 ea)	

Twist-n' Bend®

12" long

10-1510
10-1511
10-1512
10-1513
10-1514
10-1515
10-1516

Twist-Bend-Shake®

24" long

10-1520
10-1521
10-1522
10-1523
10-1524
10-1525
10-1526

36" long


10-1530
10-1531
10-1532
10-1533
10-1534
10-1535
10-1536



36" oscillate

Digi-Squeeze[®] hand exerciser

- Develops isolated finger strength, flexibility and coordination
- 3 ball sizes for small, medium and large hands
- 5 color-coded resistances
- Comes with exercise pamphlet



patented

		small	medium	large
yellow	x-light	10-1970	10-1980	10-1990
red	light	10-1971	10-1981	10-1991
green	medium	10-1972	10-1982	10-1992
blue	heavy	10-1973	10-1983	10-1993
black	x-heavy	10-1974	10-1984	10-1994
yel red grn blu blk	set of 5 (1 ea)	10-1975	10-1985	10-1995
yel red grn blu blk	set of 5 (1 ea) with stand	10-1976	10-1986	10-1996

Digi-Extend n' Squeeze[®] hand exerciser

- Provides a complete flexion and extension exercise program
- Combines the best features of Digi-Extend[®] and Digi-Squeeze[®]
- Ergonomic shape positions fingers and hand into ideal placement for best pressure distribution
- 5 color-coded balls with 4 color-coded extension bands to accommodate lightest to heaviest requirements
- Available in 3 ball sizes to fit all hands
- Each ball includes extension bands: 3 each tan (xx-light), yellow (x-light), red (light), green (med)



patented

		small	medium	large
yellow	x-light	10-2270	10-2280	10-2290
red	light	10-2271	10-2281	10-2291
green	medium	10-2272	10-2282	10-2292
blue	heavy	10-2273	10-2283	10-2293
black	x-heavy	10-2274	10-2284	10-2294
yel red grn blu blk	set of 5 (1 ea)	10-2275	10-2285	10-2295
yel red grn blu blk	set of 5 (1 ea) with stand	10-2276	10-2286	10-2296

Digi-Extend[®] hand exerciser

- Facilitates extension and flexion exercises
- Isolate extension of all fingers and individual finger joints
- Vary resistance by changing to a different color-coded resistance band
- 4 color-coded resistance band levels
- For rehabilitation or carpal tunnel, arthritis, stroke, fractures, tendon injury, nerve lacerations, tennis elbow, and more
- Exercise pamphlet included
- Clinic set contains 4 Digi-Extend[®] units with color-coded bands and metal stand



patented



clinic set with stand


- exerciser**
- 10-0775 Digi-Extend[®] exerciser set
- 10-0757 clinic set (4 piece) w/stand
- stand**
- 10-0756 stand only

replacement resistance bands - set of 25

tan	xx-light	10-1850
yellow	x-light	10-1851
red	light	10-1852
green	medium	10-1853
blue	heavy	10-1854
tan yel red grn blu	set of 25 (5 ea)	10-1855

Xtensor[®] hand strengthener

- Effective for hand therapy, tennis elbow, flexor tendon repair, post surgery rehabilitation
- Instructions included
- Comes with set of finger bands



patented

- 10-0960B Xtensor[®] exerciser
- 10-0962 foam cushion pad for Xtensor[®]
- 10-0963-8 additional finger bands (8 ea)
- 10-0964 additional thumb band
- 2 comfort pads are included with each Xtensor[®]

HandMaster Plus[™] exerciser

- Strengthens hand muscles
- Use for flexion and extension
- 3 ball firmnesses with elastic finger cord



- 10-0784 soft, purple with yellow cord
- 10-0785 medium, red with black cord
- 10-0786 hard, orange with gray cord
- 10-0787 set of 3 (1 each)

CanDo® inflatable exercise balls

- Use to improve balance, coordination, flexibility and strength, or just for fun
- Non-slip surface is ribbed for extra security

- **Standard Balls** support up to 290 lb (131kg)

Deluxe ABS Balls

- are made from a specially formulated slow deflate material. Supports a weight up to 450 lb (205 kg); ABS static weight limit 150 lb (68 kg)

- Inflator pump sold separately

standard ball - polybagged

	cm	in	color
30-1800	30	11.8	blue
30-1801	45	17.7	yellow
30-1802	55	21.6	orange
30-1803	65	25.6	green
30-1804	75	29.5	red
30-1805	85	33.5	blue
30-1806	95	37.4	red
30-1841	105	41.3	blue
30-1807	120	47.3	orange
30-1808	150	59.1	yellow

deluxe ABS ball - polybagged

	cm	in	color
30-1851	45	17.7	yellow
30-1852	55	21.6	orange
30-1853	65	25.6	green
30-1854	75	29.5	red
30-1855	85	33.5	blue
30-1857	105	42	red



CanDo® ball chairs

- Mobile ergonomic seating system combines the benefits of ball therapy with the convenience of a chair base
- Ball held securely in place by the base
- Removable back
- Light assembly required
- **Plastic Ball Chairs** are available in adult (55 cm) and child (35 cm) sizes (ball color may vary). Casters glide on carpet and hard surfaces. Locking casters available for added security (sold separately)
- **Metal Ball Chairs** are available with or without arms. They feature locking casters that glide on carpet or hard surfaces; Includes 55 cm ball (color may vary)



Plastic Ball Chair, child size, 35 cm



Metal Ball Chair, with arms, 55 cm

Plastic Ball Chair

- 30-1792 adult size w/ 55 cm ball
- 30-1795 child size w/ 35 cm ball

Metal Ball Chair

- 30-1790 without arms w/ 55 cm ball
- 30-1791 with arms w/ 55 cm ball

Accessories

- 30-1797 locking casters, 2 each
- 30-1798 replacement ball, adult size, 55 cm
- 30-1789 replacement ball, child size, 35 cm

TOGU® ABS® Powerball® exercise balls

- Anti-burst system (ABS) provides an added level of puncture-resistant safety
- ABS threshold is the dynamic weight limit at which the ball will slowly deflate without bursting
- "Double groove" design provides traction preventing the ball from rolling away
- Available in four sizes to meet each user's requirements



ABS® Powerball®

ABS limit: 200lb / 90kg

Weight Capacity: 1,100lb / 500kg

30-4000	17.7" / 45cm	yellow
30-4001	21.7" / 55cm	orange
30-4002	25.6" / 65cm	green
30-4003	29.5" / 75cm	red

ABS® Powerball® premium

ABS limit: 265lb / 120kg

Weight Capacity: 2,200lb / 1,000kg

30-4010	17.7" / 45cm	yellow
30-4011	21.7" / 55cm	orange
30-4012	25.6" / 65cm	green
30-4013	29.5" / 75cm	red

CanDo® vestibular discs

- Inflatable disc mimics movement and shape of an inflatable ball
- Available in small and large sizes
- Available in six colors: blue (B), black (BLK), red (R), green (G), yellow (Y), and silver (S). Silver available in 35 cm only.
- Colors available for personal preference; does not identify size or firmness



30-1870	small, 14" / 35 cm
30-1868	large, 24" / 60 cm

60 cm disc

CanDo® donut balls

- Rounded shape offers unique balance challenge for core strengthening
- Requires inflation (pump not included)



	cm	in
	dia x h	di x h
30-1951	yellow 45 x 25	17.7 x 9.8
30-1952	orange 55 x 30	21.7 x 11.8
30-1953	green 65 x 35	25.6 x 13.8
30-1954	red 75 x 40	29.5 x 15.7
30-1955	blue 85 x 45	33.5 x 17.7

CanDo® core trainer dome

- Inflatable dome on rigid base can be used flat-side up or down
- 21" diameter
- Complete with pump, resistance tubing and instructional DVD
- Economic alternative to Bosu® balance trainer



30-1902 CanDo® core trainer dome

Bosu® balance trainer

- Versatile balance trainer
- Use as a board or a ball
- Inflate/deflate ball to change instability
- 25" diameter



- 30-1900 exerciser with pump, plug, DVD, 120 page book and 20 page manual
- 30-1901 exerciser with DVD and 20 page manual

Bosu® NexGen™ Pro exerciser

- Unique textured markings enhance user grip
- Smooth, non-skid, non-marking base
- 26" diameter; 350 lb weight capacity



30-1905 Bosu® NexGen™ Pro

CanDo® vestibular wedges

- Can be used for balance training
- Inflatable wedge requires active participation from user for correct sitting posture
- Dynamic seating and static positioning
- Inflate or deflate wedge to vary difficulty



30-1925 large wedge (13½" x 13½")
30-1926 small wedge (10" x 10")

TOGU® Jumper® stability dome

- Functional balance trainer for jumping and swinging, strength, endurance, and conditioning
- **Pro** model has a solid plate for extra stability
- **Mini** model is ideal for one leg exercises and rehabilitation
- 440lb / 200kg weight capacity




- 30-4050 Regular 20.5" dia / 52 cm
- 30-4051 Pro 20.5" dia / 52 cm
- 30-4052 Mini 14.2" dia / 36 cm



CanDo® foam rollers

- **White PE Foam Roller** is the industry standard and provides moderate firmness
- **Black Composite Foam Roller** is made of molded foam for increased firmness
- **EVA Foam Roller** is made from extra firm EVA foam, yet it maintains a soft feel
- **Gray Plus PE Foam Roller** combines the comfortable firmness of a white roller with a gray color to conceal dirt and aging. Features rounded, soft edges
- **Full-Skin Extruded Rollers** have a full skin to prevent warping; not cut from round roll



	White PE Foam		Black Composite		EVA Foam Roller	
	Round	Half-Round	Round	Half-Round	Round	Half-Round
Dia x L						
6" x 36"	30-2100	30-2110	30-2280	30-2290	30-2200	30-2210
6" x 12"	30-2101	30-2111	30-2281	30-2291	30-2201	30-2211
	Gray Plus PE		full-skin extruded			
	Round			Half-Round		
Dia x L						
6" x 36"	30-2500			30-2340		
6" x 12"	30-2502			30-2341		



CanDo® MVP® balance system


All you need to have a full balance system in one product!

- MVP® system offers an all-in-one balance system (rocker AND wobble board)
- The board stability can be altered to various difficulty levels
- Perfect for balance, coordination, strengthening and rehabilitation of the ankles and knees, ankle range-of-motion and flexibility, core stability, improved posture, and proprioceptive training

Complete MVP® set includes:

- 10 color-coded instability balls (2 of each level from easiest to most difficult)
- 3 boards (16" hard, 20" medium and 30" hard diameter)
- Wall rack (for 10 balls)

10-1904 complete set



tie band or tubing in slots for dynamic exercise



select a board and 1 or 2 instability balls



screw the instability ball into the board



use your wobble (1 ball) or rocker (2 balls) board

Additional boards

- Circular boards have 3 ball placement options
- Rectangular board has 5 ball placement options

10-1776 16" diameter (hard)
10-1786 20" diameter (medium)
10-1796 30" diameter (easy)
10-1756 15"L x 18"W rectangular board

Additional instability balls and rack

	challenge	height	each
yellow	x-easy	1"	10-1760
red	easy	1.5"	10-1761
green	moderate	2"	10-1762
blue	hard	2.5"	10-1763
black	x-hard	3"	10-1764
10 ball set (2 ea)			10-1766
10 ball set (2 ea) with wall rack			10-1743
wall rack only -- holds 10 balls			10-1744

CanDo® Multi-Axial platform system

- Multi-Axial Platform System offers adjustable ROM (range-of-motion) and weight training to maximize lower body rehabilitation
- Reversible board (1 side for left foot and 1 side for right foot)
- Balls screw into bottom of board to provide vestibular challenge
- Place weight rod and weights on top of board to modify difficulty

Multi-Axial platform system sets include:

- Multi-Axial reversible board
- 5 color-coded instability balls (1 of each difficulty level from easiest to most difficult)
- 2 weight rods
- select sets (10-1731) and (10-1735) include 10 lbs of weight plates (4 - 2.5 lb weight plates)

10-1730 set with board, balls, rods
10-1731 set with board, balls, rods, weights, wall rack
10-1735 set with board, balls, rods, weights, storage tub



use with MVP instability balls



Shuttle® MiniPress®



- Smooth resistance for flexion and extension rehab or maintenance upper and lower body exercises
- Modify resistance by adding/removing resistance cords
- **Mini-Press Lite 2:** 4 position adjustable foot plate, 14" distance, measures 13" x 36"
- **Mini-Press:** 6 position adjustable foot plate, 20" distance, measures 12" x 42"

Mini-Press Lite 2


10-0638 1-50 lb, resistance

Mini-Press

10-0634 37-80 lb, resistance

AIREX® mats


- Airex® mats are lightweight, comfortable, and do not slide on the floor
- Made from easy-to-clean, closed-cell foam with an “integrally molded skin” for extra toughness
- Transverse ribbing pattern is molded into the mat to resist tearing and to give the mat enough flexibility to be “rolled-up” for carrying and storage
- Ideal for individual workouts, yet tough enough for use in fitness centers and therapy clinics
- Each mat is durable and easy to clean
- Corona 200 and Coronella 200 mats have an anti-bacterial treatment



size (D x L x W)	name/color	individually boxed
1 x 78 x 39"	Hercules, green	32-1230G
1 x 78 x 39"	Hercules, blue	32-1230B
5/8 x 78 x 48"	Atlas, green	32-1234G
5/8 x 78 x 48"	Atlas, red	32-1234R
5/8 x 72 x 39"	Corona, green	32-1236G
5/8 x 72 x 39"	Corona, red	32-1236R
5/8 x 72 x 39"	Corona, blue	32-1236B
5/8 x 78 x 39"	Corona 200, charcoal	32-1257C
5/8 x 78 x 39"	Corona 200, platinum	32-1257PL
5/8 x 78 x 39"	Corona 200, terra	32-1257TER
5/8 x 72 x 23"	Coronella, green	32-1238G
5/8 x 72 x 23"	Coronella, red	32-1238R
5/8 x 72 x 23"	Coronella, blue	32-1238B
5/8 x 78 x 23"	Coronella 200, charcoal	32-1258C
5/8 x 78 x 23"	Coronella 200, platinum	32-1258PL
5/8 x 78 x 23"	Coronella 200, terra	32-1258TER
2/5 x 72 x 23"	Fitline 180, aqua	32-1247A
2/5 x 72 x 23"	Fitline 180, pink	32-1247PNK
2/5 x 72 x 23"	Fitline 180, lime	32-1247LIM
2/5 x 72 x 23"	Fitline 180, charcoal	32-1247C
2/5 x 56 x 23"	Fitline 140, aqua	32-1248A
2/5 x 56 x 23"	Fitline 140, pink	32-1248PNK
2/5 x 56 x 23"	Fitline 140, charcoal	32-1248C
3/5 x 48 x 23"	Fitness 120, blue	32-1240B
3/10 x 75 x 23"	Pilates 190, purple	32-1232P
3/10 x 75 x 23"	Pilates 190, black	32-1232BLK

CanDo® balance pad

- Closed-cell foam provides excellent elasticity and energy absorption
- Coarse surface for traction on bottom and top surfaces
- Sanitized coating resists sweat and saliva and protects against microbes and odors
- 16"L x 20"W x 2.5"H, blue



32-1500B pad (ea)
32-1500B-10 pad (10 ea)

AIREX® balance pads

- Balance pads for vestibular training and exercises
- The difficulty level can be increased by stacking pads

Standard (20" x 16")

- Use for individual or group use
- 2¼" thick, blue

standard 30-1910



Solid (18" x 16")

- Extremely firm and durable
- 2" thick, blue

solid 30-1912



Elite (20" x 16")

- Waffle textures on top and bottom offer slip-resistance
- 2½" thick

blue 30-1915
kiwi 30-1915KIW
lava 30-1915LAV
pink 30-1915PNK



XL (40" x 16")

- More surface area
- 2¼" thick, blue

XL 30-1908



Mini (10" x 16")

- Ideal for traveling
- 2½" thick; lava color

mini 30-1909



CanDo® PER premium yoga mat

- Lightweight, yet cushioned and strong
- Non-skid mat surface
- 68"L x 24"W

thickness each
1/4" / 0.64 cm 30-2401B
1/6" / 0.17 cm 30-2400B



TOGU® Challenge Disc 2.0

- Use to measure and improve balancing and coordination skills
- Every move on the Challenge Disc is measured by an integrated motion sensor and results are displayed within the Challenge Disc APP using Bluetooth 4.0 technology
- Challenge Disc APP is a free download available on desktop, tablet or mobile devices (Windows, Mac, Android, and iOS)
- Use the APP with the disc to perform selective exercises systematically, efficiently and in a targeted manner
- Ideal for therapists, patients, trainers, athletes, and fitness enthusiasts of all performance levels and ages



30-4930 TOGU® Challenge Disc 2.0



Baseline® BIMS™ digital (Load-Cell) 5-position dynamometers

- Adjustable 5-position handle fits all size hands
- Accurate and repeatable grip and pinch strength measurements
- Purchase the dynamometer that meets your practice needs: **functional**, **clinic** or **deluxe**
- Subject/Test information can be inputted at time of test: dominant hand, handle position, exertion/rest times, units (lb/kg)
- 3 piece hand set includes grip and pinch dynamometers, plastic finger goniometer (12-1014) and carry case
- Bluetooth enabled to accept future APPs

2
year
warranty



Test protocols

- **Live Test:** Shows real time and maximum value. Test result is not stored
- **Quick Test:** Shows real time and maximum value. Maximum is stored
- **GST (General Strength Test):** 3 tests are performed at the same handle position. Peak, mean, SD, and COV are stored
- **RET (Rapid Exchange Test):** 6 tests (alternating hands) are performed at the same handle position
- **MMVE (Modified Maximum Voluntary Effort Test):** 5 tests are performed; 1 at each handle position
- **MVE (Maximum Voluntary Effort Test):** 15 tests are performed per hand; 3 at each handle position
- **Fatigue (Work Test):** 1 extended time test is performed. Results can be used to compare strength exerted (total work performed) over 2 or more user defined time periods

Functional model

- Shows strength readings in real time. Holds maximum value until reset. Includes **Live test**

Clinic model

- Ideal for the hands-on practitioner. Shows results in real time and stores subject/test information and test results for recall. Includes **Live**, **Quick**, and **GST** tests.

Deluxe model

- Ideal for Workers Comp, FCE (Functional Capacity Evaluation), research application and day-to-day testing. Shows results in real time and stores subject/test information and test results for recall. Includes **Live**, **Quick**, **GST**, **RET**, **MMVE**, **MVE**, and **Fatigue** tests

12-0072 grip (300 lb / 135 kg)
12-0082 pinch (100 lb / 45 kg)
12-0092 3-piece hand set

12-0070 grip (300 lb / 135 kg)
12-0080 pinch (100 lb / 45 kg)
12-0090 3-piece hand set

12-0071 grip (300 lb / 135 kg)
12-0081 pinch (100 lb / 45 kg)
12-0091 3-piece hand set

Baseline® 5 position hydraulic hand grip dynamometers

- The Baseline® hydraulic hand dynamometer gives accurate grip strength readings without the subject being able to “feel” the handle move
- Readings are repeatable
- Maximum reading remains until reset
- Readings in lb and kg
- CE certified
- Manufactured and serviced in the USA since 1978
- Recalibration available once/year based on warranty
- Five position handle and body contours assure results consistent with published Baseline® and Jamar® studies
- Internationally accepted design assures reliability, accuracy and repeatability
- Lightweight units include protective carry case, instructions for use and norms



shown in case with norms

200 lb capacity

- Inexpensive and reliable
- 200 lb (90 kg) capacity
- 2½” diameter dial
- Attractive red finish



LITE®

1 year warranty

12-0241 200 lb LITE®

- Our most popular style
- 200 lb (90 kg) capacity
- 2½” diameter dial
- Blue hammertone finish



Standard

2 year warranty

12-0240 200 lb Standard

- Maximum warranty
- 200 lb (90 kg) capacity
- 2½” diameter dial
- Black finish



HD®

5 year warranty

12-0221 200 lb HD®

Extended range 300 lb capacity

- Extended range 300 lb (136 kg) capacity
- Extra large 3½” diameter dial
- Blue hammertone finish



Hi-Res® ER™

2 year warranty

12-0246 300 lb ER™ Hi-Res®



- Extended range 300 lb (136 kg) capacity
- Electronic zero calibration system
- Max clear button
- Toggle maximum and current reading
- Toggle lb / kg
- 2½” diameter LCD gauge
- 2 AAA batteries included
- Blue hammertone finish

digital

2 year warranty

12-0247 300 lb digital

Baseline[®] Smedley spring dynamometer

- Adjustable handle comfortably accommodates any hand size
- Dial indicator remains at maximum until reset
- Readings in both lb and kg
- Includes protective carrying case



2
year
warranty

12-0281 adult (220 lb / 100 kg)
12-0282 child (110 lb / 50 kg)

Baseline[®] electronic Smedley dynamometer

2
year
warranty

- Digital, lightweight smedley style dynamometer offers accurate grip strength results
- 5-preset grip sizes
- Captures and displays maximum grip force
- 200 lb (90 kg) capacity
- Assesses results according to age and gender
- Saves/stores results for up to 19 users for fast retrieval
- Requires 2 AAA batteries included



12-0286 digital smedley (200 lb / 90 kg)

Baseline[®] pneumatic squeeze dynamometers

Single Bulb

1
year
warranty

- The inexpensive way to measure grip and pinch strength
- Available with and without a maximum force indicator (reset) that remains at the maximum reading until reset
- 2" bulb, fluid-filled
- Calibrated in PSI



with reset

12-0291 30 psi
12-0293 15 psi

without reset

12-0290 30 psi
12-0292 15 psi

Three-Bulb Set

2
year
warranty

- Three different sized bulbs (1", 1½" and 2") can be fitted to the gauge
- Maximum reading remains until reset
- Small bulb for use with children
- Unit is calibrated in PSI (0-30) and kPa (0-200)



12-0297 3-bulb dynamometer

MicroFET[®] HandGRIP[®] dynamometer

- LCD display
- For routine grip strength screening and initial/ongoing evaluations
- Durable, lightweight body
- Grip tests: maximum test, 5-position test, rapid exchange test, and hand fatigue test



2
year
warranty

12-0277W MF HandGRIP, wireless
12-0277WC with clinical software
12-0277WD with data collection software
12-0277WCD with both software packages



DynX[®] dynamometer

- For grip testing and grip therapy
- Perform maximum, endurance and rapid exchange tests
- Perform fixed effort or stepped effort grip therapy
- Store data and calculate statistics (average, standard deviation and coefficient of variation)
- Connect USB wire to data port to download archived data and statistics to a computer
- Use grip adjuster to vary grip size
- 200 lb (91 kg) capacity
- 9V battery included

2
year
warranty



12-0455 DynX[®] dynamometer

Jamar[®] dynamometers

Plus +

- Easy-to-read LCD display displays lb or kg
- 200 lb (90 kg) capacity
- Allows up to 5 trials for left and right hand
- Auto calculates average standard variation and coefficient of variation
- Rapid exchange test with audible signal

Smart


- wirelessly connects to the Jamar APP to record results



12-0604 Jamar Plus +[®] digital
12-0608 Jamar[®] Smart digital

Baseline[®] MMT (manual muscle test) kit

Use any hand dynamometer for manual muscle testing!



- MMT handle converts ANY hydraulic grip strength (hand) dynamometer to an MMT (Manual Muscle Test) dynamometer
- MMT handle comes with 3 push pad attachments: straight, curved and small circular

MMT handle with 3 push pads

12-0580 fits all hydraulic dynamometers

Baseline[®] 5-position hydraulic pinch gauge

- Easy-to-adjust paddle accommodates any hand size - standard pinch gauges only have one pinch width
- Pinch width without paddle is the same as standard pinch gauges (mechanical and hydraulic)
- 5 pinch positions allow for tip, palmar, key, rapid exchange, and 5-position pinch tests
- Users can't see results during testing
- Provides results for same test protocols as 5-position hand grip dynamometers
- Maximum reading remains until reset
- 50 lb (23 kg) capacity
- Reads in both lbs and kgs
- Results consistent with published studies




12-0480 5-position hydraulic pinch gauge
12-0480P additional paddle for pinch gauge

2
year
warranty

Baseline[®] mechanical pinch gauge

- Measure tip, key and palmar pinch strength in lb and kg
- Results are consistent with published Mathiowetz studies
- Indicator remains at the maximum reading until reset
- Includes protective case

2
year
warranty



for orthopedic and sports medicine

12-0200 30 lb (13.6 kg), blue
12-0201 60 lb (27 kg), red

for weak or damaged hands

12-0203 10 lb (4.5 kg), silver
12-0202 2 lb (0.9 kg), gold

Baseline[®] hydraulic pinch gauge

- Measure tip, key and palmar pinch strength in lb and kg
- Includes protective case
- Hi-Res models have a larger 3.5" diameter dial
- Digital model clearly shows results on 2.5" LCD screen



12-0226 50 lb LITE[®] 1 year warranty
12-0235 50 lb standard 2 year warranty
12-0222 50 lb HD[®] 5 year warranty

extended range (100 lb / 45 kg capacity)

12-0228 100 lb ER[™] Hi-Res[®] 2 year warranty
12-0237 100 lb digital 2 year warranty

Baseline[®] electronic pinch gauge

- LCD display shows grip force in lb or kg
- Perform tip, key, and palmar tests
- Stores 99 readings
- 60 lb (27 kg) capacity
- 3V battery included

2
year
warranty

12-0475 Baseline[®] electronic pinch gauge





standard
7-piece
set

Choose your warranty

LITE®
1 year
warranty

**Standard,
Digital, Hi-Res®**
2 year
warranty

HD®
5 year
warranty



standard
3-piece
set

Baseline® 7-piece hand sets

- Protective carrying case stores instruments
- Set includes choice of hydraulic grip dynamometer and mechanical pinch gauge
- Also includes stainless steel goniometer (12-1010), 2-point discriminator with 3rd point (12-1481), Wartenberg pinwheel (12-1450), finger circumference gauge (12-1222) and functional finger motion gauge (12-1060)

grip / pinch			warranty
12-0100	200 lb (12-0240) / 30 lb (12-0200)	standard	2 year
12-0102	200 lb (12-0241) / 30 lb (12-0200)	LITE®	1 year
12-0117	200 lb (12-0221) / 30 lb (12-0200)	HD®	5 year
12-0111	300 lb (12-0246) / 60 lb (12-0201)	ER™ Hi-Res®	2 year
12-0113	300 lb (12-0247) / 60 lb (12-0201)	ER™ digital	2 year

Baseline® 3-piece hand sets

- Protective carrying case stores instruments
- Includes choice of hydraulic grip dynamometer, hydraulic pinch gauge and 6" stainless steel goniometer (12-1010)

grip / pinch			warranty
12-0103	200 lb (12-0240) / 50 lb (12-0235)	standard	2 year
12-0101	200 lb (12-0241) / 50 lb (12-0226)	LITE®	1 year
12-0116	200 lb (12-0221) / 50 lb (12-0222)	HD®	5 year
12-0110	300 lb (12-0246) / 100 lb (12-0228)	ER™ Hi-Res®	2 year
12-0112	300 lb (12-0247) / 100 lb (12-0237)	ER™ digital	2 year

Baseline® 3-piece 5-position hand set

- Pinch gauge faces away from subject (toward practitioner) to block readings from subject
- Set includes: 5-position pinch gauge, 200 lb standard grip dynamometer and black plastic digit goniometer
- 2 year warranty



12-0160 set w/ case 2 year warranty

set contents

12-0240	standard grip dynamometer (200 lb)
12-0480	5-position pinch gauge (50 lb)
12-1014	black plastic digit goniometer

Baseline® 8-piece hand sets

8-piece portable carrying case contains:

- 200 lb (90 kg) hydraulic grip dynamometer
- 50 lb (23 kg) hydraulic pinch gauge
- Finger goniometer (12-1010)
- Finger circumference gauge (12-1222)
- 180° 6" plastic goniometer (12-1005HR)
- Measuring tape (12-1210)
- Discrim-a-Gon® (12-1492)
- 5-pc tactile™ monofilament set (12-1662)



grip / pinch			warranty
12-0128	200 lb (12-0240) / 50 lb (12-0235)	standard	2 year
12-0126	200 lb (12-0241) / 50 lb (12-0226)	LITE®	1 year
12-0127	200 lb (12-0221) / 50 lb (12-0222)	HD®	5 year

Baseline® 3-piece wrist dynamometer sets

Wrist strength can now be measured

- Simply and accurately measure the strength of the wrist muscles during flexion, extension, abduction and adduction and the forearm muscles during supination and pronation
- Maximum reading remains until reset
- Measurement in lb and kg (500 lb / 226 kg)
- Includes protective carrying case
- Manufactured and serviced in the USA
- CE certified
- Set includes door knob grip, table mount and protective carry case



12-0250 analog (dial) readout
12-0251 digital (LCD) readout

2 year warranty

optional measurement accessories

- Measurement accessories easily snap on and off



door knob grip
radial / ulnar deviation
(comes standard with dynamometer set)
12-0261

rod / pipe grip
flexion / extension
12-0264



shovel handle
supination / pronation
12-0262



lever handle
abduction / adduction
12-0263

Baseline® push-pull dynamometers

- Lightweight, small and portable; easy to grasp with or without handles
- Use for muscle strength management, job task analysis, or functional capacity evaluation
- Available in hydraulic, electronic or mechanical spring mechanism
- Comes with 3 push pads and 2 pull hooks

Hydraulic

analog
12-0392 50 lb / 23 kg
12-0393 100 lb / 45 kg
12-0394 250 lb / 113 kg
12-0388 500 lb / 226 kg

digital
12-0397 50 lb / 23 kg
12-0398 100 lb / 45 kg
12-0399 250 lb / 113 kg
12-0387 500 lb / 226 kg



2 year warranty

Electronic

- "Load-cell" technology assures accuracy, repeatability and reliability

12-0340 50 lb / 23 kg
12-0341 100 lb / 45 kg
12-0342 250 lb / 113 kg
12-0343 500 lb / 226 kg



2 year warranty

Mechanical spring

12-0300 2.2 lb / 1 kg
12-0301 5.5 lb / 2.5 kg
12-0302 11 lb / 5 kg
12-0304 22 lb / 10 kg
12-0303 66 lb / 30 kg




1 year warranty

ErgoFET 500™ push-pull dynamometer



2 year warranty

- Measure peak force, average or 3-second average force
- Push and pull attachments allow for a wide variety of strength measurement tests
- Use removable handles when measuring higher forces
- 500 lb push capacity, 200 lb pull capacity
- Includes: 3 push attachments, 2 pull hooks, D-ring strap, extension post, test forms, handles, wall chart, case, 3.7V rechargeable lithium battery and AC adapter



12-0460W EF500, wireless
12-0460WD EF500, with data collection software

MicroFET[®] 2 MMT (manual muscle tester)

- Displays peak force and time duration
- Versatile push attachments ensure proper muscle isolation for increased testing validity
- Includes: 3 push attachments, Test forms, Wall chart, Case, two 3.7LV lithium batteries, AC adapter, Battery charger
- 300 lb (136 kg) capacity
- Connects via bluetooth with plug 'n play software



2
year
warranty

12-0381W	MF2, wireless
12-0381WC	MF2, with clinical software
12-0381WD	MF2, with data collection software
12-0381WCD	MF2, with both software packages

ErgoFET[®] MMT with pull function

- Displays peak force and time duration
- Push attachments ensure proper muscle isolation for increased testing validity
- Pull attachments enable complete test capabilities
- Includes: 3 push attachments, 2 pull hooks, D-ring strap, Extension post, Test forms, AC adapter, 2 3.7V rechargeable lithium batteries
- 300 lb (136 kg) capacity



2
year
warranty

12-0461W	EF, wireless
12-0461WD	EF, with data collection software
12-0461WC	EF, with clinical software
12-0461WCD	EF, with both software packages

Lafayette MMT (manual muscle tester)

- Measures peak force and time to achieve it, total test time, time within selectable ranges and average force
- Comes with 3 push pads
- Optional data storage for up to 150 tests
- 300 lb (136 kg) capacity
- Optional data transfer software allows bluetooth communication between MMT and computer
- Rechargeable battery with auto shut-off
- Easy-to-read graphical LCD display reads in lb, kg, and newtons



2
year
warranty

12-0380	Lafayette MMT
12-0380S	Software/wireless data

Baseline[®] back-leg-chest dynamometers

- Measure strength of back, leg and chest muscles
- Base provides secure footing
- Chain length is adjustable to accommodate for height differences or to vary the point of force application
- Comes fully assembled with 5-foot chain
- Oversize base 24" x 24"; regular base 18" x 18"
- Pointer remains at maximum until reset
- Reads lb and kg



12-0400	adult with standard base	660 lb (300 kg)
12-0403	adult with oversize base	660 lb (300 kg)
12-0401	adolescent with standard base	330 lb (150 kg)
12-0402	child with standard base	165 lb (75 kg)

Baseline[®] fingertip pulse oximeters

Standard Adult

- Measures pulse and SpO₂
- Compact and lightweight
- Includes lanyard and 2 AAA batteries



12-1926	fingertip pulse oximeter
12-1926-25	25 each

Deluxe Adult

- 2-color screen displays SpO₂, PR, waveform and pulse bar
- 6 display modes with adjustable brightness
- Auto shut-off
- 2 AAA batteries and lanyard included



12-1927	deluxe pulse oximeter
12-1927-25	25 each

Pediatric

- Small design is perfect for pediatric use
- Measures pulse and SpO₂
- 6 display modes with adjustable brightness
- Battery low indicator



12-1935	pediatric pulse oximeter
12-1935-25	25 each

Baseline® volumetric edema gauges

- Objectively measure body part volume by using the fluid displacement method
- Body part is placed in a water-filled tank; the water overflow is collected and measured
- Amount of edema can be measured and quantified

Each set includes:

- Measurement tank
- Overflow beaker
- Graduated cylinder




upper extremity (L x W x H)

- 12-3500 hand set 5" x 3" x 9"
- 12-3504 forearm set 6" x 6" x 24"
- 12-3501 arm set 7" x 7" x 30"

lower extremity (L x W x H)

- 12-3502 foot set 13" x 5" x 9"
- 12-3505 partial leg set 13" x 6" x 24"
- 12-3508 large partial leg set 13" x 8" x 16"

ISOM measurement chart




- 36" x 30" wall chart illustrates measurement and recording rules from the American Academy of Orthopedic Surgeons
- Goniometry, x-ray evaluation & anthropometric measurement

12-1095 ISOM chart

Baseline® finger flexion gauge

- Clear plastic gauge simplifies and speeds measurement of composite finger flexion and finger-thumb opposition
- Scale measures up to 10 cm



12-1060 finger flexion gauge

Baseline® circumference tape

- Measure the circumference of a body part in both cm and inches



12-1205 tape
12-1205-25 25 each

Baseline® measuring tape



- Reads in cm and inches
- Locks in place until reset button is pressed
- 3 lengths available

12-1210 5-foot tape (ea)
12-1210-25 25 each
12-1211 6-foot tape (ea)
12-1211-25 25 each
12-1212 10-foot tape (ea)
12-1212-25 25 each

Baseline® Gulick measurement tape



- Spring exerts constant force on tape for accurate and repeatable results
- Push-button retractor
- Reads in cm and inches

12-1201 5-foot (ea)
12-1201-25 25 each
12-1203 6-foot tape (ea)
12-1203-25 25 each
12-1204 10-foot tape (ea)
12-1204-25 25 each

Baseline® finger circumference gauge



- Wrap webbing around finger or digit and read circumference in both cm and inches

12-1220 gauge (in)
12-1221 gauge (cm)
12-1222 gauge (in/cm)
12-1222-25 25 each

Baseline® finger (digit) goniometers

- Measures 110° flexion through 40° hyperextension in 2° increments on plastic model (5° increments on metal model)
- Arm has 10 cm/4" linear scale
- Use one finger to move the measurement head through flexion, extension and hyperextension
- Readings can be taken using only one hand; freeing the therapist's other hand to stabilize the patient during measurement or to record the measurement; standard finger goniometers require two hands to take a finger ROM measurement

Use only one hand to measure your patient!



plastic



stainless steel



plastic digit goniometer

- 12-1014 plastic digit goniometer, each
12-1014-25 plastic digit goniometer, 25 each

stainless steel digit goniometer

- 12-1016 stainless steel digit goniometer, each
12-1016-25 stainless steel digit goniometer, 25 each

180° 8" goniometer



- Two opposing scales marked in 1° increments
- Tension between the arms is controlled by the thumb knob
- Stainless steel

12-1040 8", 180°

180° 14" goniometer



- Two opposing scales marked in 1° increments
- Extra-long arms for measuring all joints
- Stainless steel

12-1041 14", 180°

360° 14" goniometer



- Two opposing 180° scales and one 360° scale, all marked in 1° increments
- Arm tension is controlled by thumb knob
- Stainless steel

12-1050 14", 180°

180° 6" goniometer



- Robinson pocket-sized goniometer has a 180° scale in 5° increments and a 5" (12½ cm) linear scale
- Stainless steel

12-1042 6", 180°
12-1042-25 25 each

Baseline® 180° finger (small joint) goniometer



- Stainless steel goniometer measures ROM of finger joints (metacarpophalangeal and interphalangeal) and other small joints including toes
- Head has opposing 150° scales in 5° increments
- Linear scale in cm and inches
- Deluxe model has larger head and full 180° range


	each	25 each
standard finger, 6"	12-1010	12-1010-25
deluxe finger, 6"	12-1011	12-1011-25
small finger, 3½"	12-1015	12-1015-25

6-piece goniometer set



- 6 goniometers in a padded carrying case
 - 6" 180° Robinson pocket
 - 14" 180° stainless steel
 - 14" 360° stainless steel
 - 6" stainless steel finger
 - 8" 180° x-ray
 - 8" 180° stainless steel
- 12-1043 goniometer set

180° X-ray Radiopaque goniometer




- Document angles on X-ray negatives
- 180° head has two opposing scales marked in 2½° increments
- Black aluminum with easy-to-read white radiopaque markings


12-1030 14"
12-1031 8"
12-1032 6"

Baseline® plastic goniometers

- Scale reads in 1 degree increments; linear scale reads in centimeters and inches
- Hi-Res® option has white background behind the gradations and numbers for enhanced readability




Absolute+Axis®



- Take readings from absolute vertical or horizontal

Absolute+Axis® attachment slides onto any 12" goniometer



Patient

clinic - plastic ISOM

12" ISOM 360°	each	12-1000
12" ISOM Hi-Res® 360°		12-1000HR
8" ISOM 360°		12-1001
8" ISOM Hi-Res® 360°		12-1001HR
6" ISOM 360°		12-1002
6" ISOM Hi-Res® 360°		12-1002HR

patient - plastic pocket

6" pocket 180°		12-1005
6" pocket Hi-Res® 180°		12-1005HR
6" rulometer 360°		12-1006
6" rulometer Hi-Res® 360°		12-1006HR

finger - plastic

flexion/hyper-extension		12-1012
flexion/hyper-extension Hi-Res®		12-1012HR

plastic sets

standard set - 6 pieces		12-1028
Hi-Res® set - 6 pieces		12-1028HR

absolute axis - plastic

12" Absolute+Axis®		12-1025
12" Absolute+Axis® Hi-Res®		12-1025HR
AA attachment only		12-1026

Finger



finger flexion/hyper-extension

25 each

12-1000-25
12-1000HR-25
12-1001-25
12-1001HR-25
12-1002-25
12-1002HR-25
12-1005-25
12-1005HR-25
12-1006-25
12-1006HR-25
12-1012-25
12-1012HR-25
12-1028-25
12-1028HR-25
12-1025-25
12-1025HR-25
12-1026-25


Baseline® XTender® extendable goniometer



- Arms extend to measure large and small joints
- Extendable arms allow the goniometer to lie over anatomical bony landmarks
- Head is the same as a traditional goniometer for easy use
- Metal legs extend from 10" to 36"

12-1034HR extendable goniometer

Lafayette Gollehon goniometer



- For measuring large and small joints
- 4x magnifying lens to read the dual 180° high-resolution scale and two stainless steel arms extending 9½" to 26"
- Smooth precision glide

12-1035 Gollehon

Baseline® digital goniometers

- Digital plastic goniometer reads 0-360° on a LCD screen; 10" linear scale
- Absolute+Axis® goniometer has vertical and horizontal levels incorporated into its digital goniometer arm. It reads 0-180° on a LCD screen




	each
10" plastic 360°	12-1049
9" Absolute+Axis®	12-1027

25 each

12-1049-25
12-1027-25

Baseline® large joint (arthrodial) protractor

- Measure ROM for all major articulations, cervical rotation, lateral flexion of the head and anterior-posterior cervical flexion



12-1076 arthrodial protractor

Baseline® clinical monofilament sensory evaluators

- Same trusted product in new color-coded system
- Measure sensitivity, nerve damage or impairment, changes in neurological status
- Color code reference for recording results: green, normal; blue, diminished light touch; purple, diminished protective sensation; red, loss of protective sensation
- **Fold-Up™ Monofilaments:** top folds to protect monofilament without an additional case
- **Tactile™ Monofilaments:** come with protective case

size	force		color	Tactile™ w/ case	Fold-Up™	
1.65	(0.008 gm)	1 ea	normal	12-1640	12-1720	
2.36	(0.02 gm)	1 ea		12-1641	12-1721	
2.44	(0.04 gm)	1 ea		12-1642	12-1722	
2.83	(0.07 gm)	1 ea		12-1643	12-1723	
3.22	(0.16 gm)	1 ea	diminished light touch	12-1644	12-1724	
3.61	(0.4 gm)	1 ea		12-1645	12-1725	
3.84	(0.6 gm)	1 ea	diminished protective sensation	12-1646	12-1726	
4.08	(1.0 gm)	1 ea		12-1647	12-1727	
4.17	(1.4 gm)	1 ea		12-1648	12-1728	
4.31	(2.0 gm)	1 ea		12-1649	12-1729	
4.56	(4.0 gm)	1 ea	loss of protective sensation	12-1650	12-1730	
4.74	(6.0 gm)	1 ea		12-1651	12-1731	
4.93	(8.0 gm)	1 ea		12-1652	12-1732	
5.07	(10.0 gm)	1 ea		12-1653	12-1733	
5.07	(10.0 gm)	100 ea		12-1653-100	12-1733-100	
5.18	(15.0 gm)	1 ea		12-1654	12-1734	
5.46	(26.0 gm)	1 ea		12-1655	12-1735	
5.88	(60.0 gm)	1 ea		12-1656	12-1736	
6.10	(100.0 gm)	1 ea		12-1657	12-1737	
6.45	(180.0 gm)	1 ea		12-1658	12-1738	
6.65	(300.0 gm)	1 ea		12-1659	12-1739	
2-piece, normal, 2.83, 3.61				12-1660	12-1745	
2-piece, protective sensation 4.56, 5.07				12-1668	12-1746	
5-piece, hand, 2.83, 3.61, 4.31, 4.56, 6.65				12-1662	12-1740	
6-piece, foot, 2.83, 3.61, 4.31, 4.56, 5.07, 6.65				12-1664	12-1741	
20-piece, complete				12-1666	12-1743	



monofilament accessories

- Use to keep track of evaluation results;
1 pad with 100 sheets
- 12-1360 hand form, 1 pad
12-1361 foot form, 1 pad
12-1362 5 pencils (R, G, B, P, BLK)

Baseline® Fold-Up™ monofilaments



Baseline® Tactile™ monofilaments



Baseline® disposable monofilaments

- Measure and detect skin sensitivity, nerve damage or impairment, and changes in neurological status
- Used for screening peripheral nerve impairment
- Disposable 5.07 / 10 gm monofilaments



12-1670	LEAP	1 ea
12-1670-20	LEAP	20 ea
12-1670-40	LEAP	40 ea
12-1671	ADA	1 ea
12-1671-20	ADA	20 ea
12-1671-40	ADA	40 ea
TAA compliant for US government		
12-1672S-1	LEAP / ADA / LOPS	1 ea
12-1672-20	LEAP / ADA / LOPS	20 ea
12-1672-40	LEAP / ADA / LOPS	40 ea

Baseline® folding home monofilaments

- Inexpensive home monofilaments detect changes in neurological status
- Screen for peripheral nerve impairment and diabetic foot
- Monofilament folds into base



	size	force
12-1395	4.08	1 gm
12-1396	5.07	10 gm
12-1397	6.10	100 gm
12-1398	3 piece set	

Baseline® retracting monofilaments

- Filament retracts into handle
- Filament detects changes in neurological status
- Screen for peripheral nerve impairment and diabetic foot



	size	force
12-1390	4.08	1 gm
12-1391	5.07	10 gm
12-1392	6.10	100 gm
12-1393	3 piece set	

Baseline® 4 monofilament pen

- Choose 1 of 4 sensory evaluators (0.07, 0.4, 4.0, 10.0 gm)
- Monofilaments safely retracts into pen barrel when not in use
- Can be carried in shirt pocket



12-1370 4 monofilament pen

WEST™ sensory evaluators


- WEST (Weinstein Enhanced Sensory Test) uses five monofilaments to measure skin sensitivity and peripheral nerve involvement
- Tips are round and textured
- Includes case (4½" by 6½"), test form and manual
- WEST-D includes range from face to foot



12-1406	WEST hand	0.07, 0.2, 2, 4, 200gm
12-1407	WEST foot	0.5, 2, 10, 50, 200 gm
12-1408	WEST-D	0.02, 0.07, 0.5, 2, 10 gm

Baseline® dolorimeter


- Evaluate pain threshold
- Standardized (1.52 cm²) flat circular probe is pushed against subject's skin until pain threshold is reached
- Reads lb and kg




12-1440	2.2 lb / 1 kg
12-1441	5.5 lb / 2.5 kg
12-1444	11 lb / 5 kg
12-1442	22 lb / 10 kg
12-1443	66 lb / 30 kg

Baseline® Discrim-A-Gon® 2-point discriminator

- Easy-to-use, lightweight plastic wheel is the perfect sensory evaluation tool to test static and dynamic 1 and 2-point discrimination
- 2 separate 2-point discrimination octagons (D1 and D2) come in set
- Each octagon measures a different range of 8 labeled, fixed 2-point intervals ranging from 1 to 25 mm



set includes
2 separate
discrimination
octagons



12-1492	Discrim-A -Gon® set (2 discs)
12-1492-25	25 sets

Touch-Test® 2-point discriminator

- Tests static, moving 2-point discrimination
- Consists of two sturdy, rotating, plastic disks that are joined together
- Includes instructions and carrying case



12-1495 Touch Test® 2-point

Disk-Criminator™ 2-point discriminator

- For testing static and dynamic 1 and 2-point discrimination
- Therapeutic aid for desensitization and home programs following nerve repair, graft and innervated tissue damage
- Includes 2 disks that quantify innervation density from 1 to 25 mm



12-1490 2-point discriminator (2 discs)

Baseline® 2-point discriminators

- Evaluate cutaneous sensitivity and touch threshold
- Third point permits alternating between single and double point stimulation without changing setting
- Plastic points minimize influence of temperature
- Measures up to 14 cm




A plastic two point



B plastic with third point

A	12-1480	2-point discriminator
B	12-1481	with third point

Baseline[®] tuning forks



- Use to evaluate audio and vibratory sensation

tuning fork sets

6 piece set with case

128, 256, 512, 1024, 2048, 4096 cps 12-1460

2 piece sensory set

30 and 256 cps 12-1461

weighted tuning forks

	each	25 each
30 cps weighted	12-1465	12-1465-25
64 cps weighted	12-1464	12-1464-25
128 cps weighted	12-1466	12-1466-25
256 cps weighted	12-1467	12-1467-25

unweighted tuning forks

512 cps	12-1468	12-1468-25
1024 cps	12-1469	12-1469-25
2048 cps	12-1470	12-1470-25
4096 cps	12-1471	12-1471-25

student grade tuning forks


128 cps weighted	12-1475	12-1475-25
256 cps	12-1476	12-1476-25
512 cps	12-1477	12-1477-25

special tuning forks

		each
A Rydel-Seiffer tuning fork (64 and 128 cps)		12-1472
B variable frequency tuning fork (128 to 240 cps)		12-1462

Baseline[®] Taylor hammers

- Use to elicit responses during examination



Taylor Hammer


	each	25 each
C Taylor, red hammer	12-1500	12-1500-25
D with 128 cps tuning fork	12-1501	12-1501-25
E with 256 cps tuning fork	12-1502	12-1502-25

Taylor Hammer, Latex-Free

F black	12-1574	12-1574-25
G blue	12-1573	12-1573-25
H green	12-1572	12-1572-25
I yellow	12-1570	12-1570-25
J red	12-1571	12-1571-25

Baseline[®] Buck neurological hammers


- Two heads, point/brush for cutaneous sensation evaluation
- Buck hammer has a pinwheel
- Buck with tuning fork also has pinwheel



	each	25 each
K Buck hammer	12-1510	12-1510-25
L with 128 cps tuning fork	12-1511	12-1511-25
M with 256 cps tuning fork	12-1512	12-1512-25
N with pinwheel, brush, point	12-1516	12-1516-25
Buck, latex-free hammer (not shown)	12-1514	12-1514-25

Baseline[®] Babinski hammers

- Adjustable head can be oriented in two directions; use pointed bottom as a "gross" aesthesiometer




Babinski hammer

	each	25 each
standard, 8½"L	12-1520	12-1520-25
extra long body, 10"L	12-1521	12-1521-25
standard, latex-free, 8½"L	12-1524	12-1524-25
extra long body, latex-free 10"L	12-1526	12-1526-25

pinwheels

- Test for nerve damage or nerve dysfunction
- Used by orthopedists, neurologists, chiropractors and other healthcare practitioners to test a patient's ability to feel sharp contact with the skin
- Consists of a rotating pinwheel attached to a handle



Baseline[®] Cleanwheel[®]

	single use
12-1446	non-sterile Cleanwheel [®]
12-1446-25	25 each
12-1445	sterile Cleanwheel [®]
12-1445-25	25 each

Wartenberg



12-1450	Wartenberg
12-1450-25	25 each

Neuro-Aide



	single use
12-1455	Neuro-Aide
12-1455-25	25 each

Baseline® Bubble® inclinometer

- Place inclinometer near the joint to measure, turn dial until scale reads 0, take joint through its range, read range from dial
- Some measurements need the use of 2 inclinometers
- Standards codified in "AMA Guide to the Evaluation of Permanent Impairment, third edition"



- 12-1056 Bubble® inclinometer (each)
 12-1056-2 Bubble® inclinometer (set of 2)
 12-1056-25 Bubble® inclinometer (25 each)

Baseline® AcuAngle® inclinometers

- Place inclinometer near joint to be measured; turn dial to 0; take joint through its range; read ROM from dial
- AcuAngle® Inclinometer** has adjustable feet along scale (side-to-side) that adapts to body contours
- Pointer is dampened by fluid to assure accurate ROM measurements
- AcuAngle® Wrist Inclinometer** measures ROM of the wrist along all axes



- 12-1149 AcuAngle® inclinometer (each)
 12-1149-2 AcuAngle® inclinometer (set of 2)
 12-0502 AcuAngle® wrist inclinometer
 12-0503 AcuAngle® wrist handle only

Baseline® digital inclinometer

- Place inclinometer near the joint to be measured; press "zero" button; move joint through its range; press the "hold" button
- Read ROM directly after the joint has been taken through its range



- 12-1057 digital inclinometer
 12-1057-2 set of 2

Baseline® scoliometer


- Measure the degree of rotation of a deformity of the back
- Screen which people to refer for further medical evaluation
- 7½"L x 3½"W
- Includes storage pouch



- 12-1099 scoliometer

spinal assessment tools

- BROM Basic** - BROM rotation lateral unit plus two universal inclinometers
- CROM Basic** - CROM main frame, rotation arm and video
- CROM Deluxe** - Cervical Range-of-Motion instrument
- PALM** - PALpation Meter (skeletal alignment and leg length discrepancy instrument)
- BROM II** - Back Range-of-Motion instrument
- PAA** - Deluxe inclinometer
- Deluxe Universal Inclinometer**
- CROM 3**



- A 12-1151 BROM basic
 B 12-1155 CROM basic
 C 12-1156 CROM deluxe
 D 12-1180 PALM
 E 12-1150 BROM II
 F 12-1097 PAA universal inclinometer
 G 12-1066 Deluxe universal inclinometer
 H 12-1157 CROM 3



Baseline[®] Sit 'n Reach[®] inclinometers

- Perform sit 'n reach flexibility tests
- Measurements include: fitness program, Cooper Institute/YMCA, SHAPE America (formerly known as AAHPERD)
- Scales printed in cm and inches
- Sturdy metal construction with a built-in foot plate and a maximum stretch indicator that remains until reset
- **Deluxe model** offers adjustable start point; use to normalize scores to account for variation in the arm and leg length of the person being tested


12-1085 standard
12-1086 deluxe



Allen diagnostic tools

Occupational Therapy Treatment Goals Book

- Therapists use this book to interpret and document results of ADM
- Softcover, 360 pages



12-3161 occupational therapy treatment goals book

Module Instruction Manual

- Manual describes standardized procedures for 35 activities including making samples, instructions, prompts and rating behavior
- Includes score sheets
- Binder, 26 sections



12-3150 instruction manual, 2nd ed.

Cognitive Levels Guide

- Color-coded reference guide is the starting point to learn the concepts, terminology and criteria for the Allen Cognitive Level Modes of Function
- Spiral bound, 82 pages



12-3159 Allen cognitive levels guide

cognitive level screen tool

- Standardized ACLS screening tool used to obtain quick level of processing capacities, learning potential and performance abilities
- Manual is updated for administering ACLS and LACLS; 66 pages
- **ACLS kit** includes 72 waxed linen threads and 120" lace
- **LACLS kit** includes 1 shoelace and 120" lace



12-3154 screen tool - 3¾" x 4¾"
12-3155 large screen tool - 5¾" x 7"

Cognitive Level Manual

12-3152 for Screen-5 and large Screen-5
12-3153 for Screen-5 and large Screen-5 with tools

Replacement Kits

12-3156 for ACLS testing
12-3157 for LACLS testing

canvas placemet module

- Includes canvas placemats (10½" x 13") and precut felt hearts and circles
- For levels 3.0 to 4.6



12-3151 canvas placemats, 6 ea.

tile trivet module

- Includes assorted trivets (tiling area 4½" sq), presorted tiles and sandpaper
- Pattern designs vary
- For levels 3.0 to 5.8



12-3164 tile trivets, 6 ea.

leather key fob module


- Tests sequencing to evaluate learning
- Includes leather circles, crow beads, key rings, elastic cord and rawhide lacing
- For levels 4.0 to 5.0



12-3177 leather key fobs, 6 ea.

ribbon card module

- Includes precut (5" x 7") greeting cards with envelopes and precut ribbon for weaving
- For levels 4.0 to 5.0



12-3172 ribbon cards, 8 ea.

Purdue® pegboard

- Measure gross movements of hands, fingers and arms; and finger-tip dexterity
- Pins, collars, washers, examiner's manual and norms included



12-3000 Purdue®
12-3002 record blanks (25)
12-3001 extra parts:
25 collars, 45 washers, 55 pins

Lafayette® pegboard

- Measure a timed unilateral motor response
- Subject places 25 pegs in round holes as quickly as possible across the rows



12-3005 Lafayette

grooved pegboard


- Manipulative dexterity test consisting of 25 holes with randomly positioned slots
- Pegs with a "key" along one side must be rotated to match hole before inserted
- Requires more complex visual-motor coordination



12-3010 grooved
12-3011 extra pegs (32)


Roeder manipulation test

- Measures hand, arm and finger dexterity and speed
- Test sorting and assembling
- Complete with instructional manual, norms and record blanks



12-3080 Roeder
12-3081 record blanks (50)
12-3082 extra parts:
25 collars, 45 washers, 55 pins


Minnesota Manual Dexterity test



- Manual dexterity test measures simple eye-hand-finger movement
- Complete with carrying case, record blanks and norms

12-3050 MN Dexterity
12-3051 extra pegs, each
12-3052 record blanks (50)

complete Minnesota Dexterity test



- Measures eye-hand coordination, speed and dexterity
- Tests displacing, and one-hand and two-hand turning and placing
- 2 folding boards, 60 blocks, carrying case and manual

12-3060 complete MN
12-3051 extra pegs, each


work bench



- Improve strength and endurance of shoulders, elbows, wrists and hands; fine motor coordination
- Front allows bilateral work that duplicates a hidden, confined work area
- Toolbox with 6" and 8" adjustable wrenches, 2 screwdrivers and 4 hex keys

10-1163 work bench

hand-tool dexterity test



- Measures proficiency by using simple mechanic's tools
- Consists of tools and two uprights with bolts
- Disassemble all bolts from one upright and reassemble them on the other upright

10-1164 hand-tool test

O'Connor dexterity tests



- Finger test requires hand placement of three pins per hole
- Tweezers test requires the placement of a single pin in a 1/16" diameter hole

12-3020 O'Connor finger
12-3030 O'Connor tweezer
12-3031 extra pins (105)

Baseline® 9-hole pegboard



- Improve fine motor coordination and observe quality of movement and speed
- 14 pegs stored in two compartments within board
- Square board 4-7/8"; pegs are 1 1/2"L x 1/4" dia

12-3141 9-hole pegboard

Corbett Targeted Coin Test (CTCT)



- Report accuracy and speed of in hand manipulation
- Assess palm to finger tip in hand manipulation
- Assess fine motor manipulation with dart thrower's motion

12-3400 CTCT coin test

stopwatch

- Start-stop, timeout, split time, calendar modes



12-2100 stopwatch

timer

- Rings after pre-set period elapses



12-2004 60-minute timer

Baseline® pedometers

- Use on-the-go to easily monitor steps taken
- Select models track steps taken, distance traveled, calories burned and activity time



	each	25 each
LITE® steps taken	12-1954	12-1954-25
Standard steps, distance	12-1955	12-1955-25
Deluxe steps, distance, calories, time	12-1956	12-1956-25

HeartSmart on-site AED with slim carry case



- Portable defibrillator can be used on children and adults with sudden cardiac arrest

13-4501 on-site AED

stethoscopes

- Nurse's:** anodized aluminum chest piece, plastic ear tips
- Dual Head:** turn to convert to non-child ring; external spring binaural
- Sprague-Rappaport:** rotating drum stem and external spring binaural



12-2200	nurse's, each
12-2200-25	nurse's, 25 each
12-2210	dual-head, each
12-2210-25	dual-head, 25 each
12-2211	stainless steel, latex-free
12-2220	Sprague-Rappaport

aneroid sphygmomanometers

wall mount

- 90° swivel mount
- Latex-free cuff
- With adult inflation system



12-2261 wall mount

mobile floor

- Adjustable height caster base, double cuff container
- With adult inflation system



12-2262 mobile floor

portable

- High quality with adult inflation system and carry bag



12-2250 unit w/case
12-2250-25 25 each

cuffs only

- 2-tube cuff bladder system
- Use with handheld, wall mount or floor aneroid or mercurial sphygmomanometers



	W x L
12-2243 adult	5½ x 20"
12-2244 x-large	7 x 24"
12-2242 child	4 x 13½"
12-2241 infant	3 x 10 1/3"
12-2240 newborn	2 x 7¼"

Baseline® economy skinfold caliper



with floating measurement tips

- Plastic skinfold caliper is the inexpensive choice for skinfold analysis
- Spring-loaded arm assures reproducible results
- Floating tips assure accurate and reproducible results
- 50 mm scale
- With booklet and tables

12-1112 1 each
12-1112-25 25 each

Baseline® skinfold calipers

- 2-sided lightweight aluminum caliper is easy to use and read for left- or right-handed people; includes calibration block
- Floating tips to assure accurate and reproducible results
- Large, easy-to-read 70 mm scale
- Includes protective carrying case
- Calibration block verifies at 5-50 mm

skinfold calipers
12-1130 deluxe (dual sided) with calibration block
12-1110 skinfold caliper (single sided) without calibration block

accessory
12-1116 calibration block only

double sided
2-sided deluxe skinfold calipers viewable on front and back









Baseline® body fat scale



- Tempered glass scale measures weight, body fat %, water %, muscle %, bone density
- Features an LCD screen, auto on, 5 activity levels and memory for 12 users
- 330 lb/150 kg capacity
- 3V battery included

12-1190 body fat scale

anatomical charts



A hand and wrist
12-4609 P/L/S
D lymphatic system
12-4613 P/L

B shoulder and elbow
12-4619 P/L
E human musculature
12-4614 P/L/S

C nervous system
12-4628 P/L/S
F human skeleton
12-4620 P/L/S

SPECIFY: P (Paper), L (Laminated), or S (Sticky)

paper	(P)	20" x 26"
laminated	(L)	20" x 26"
sticky back	(S)	18" x 24"

muscle arm, 6-part models

Illustrates both the superficial and deeper muscles, five of which are removable. Tendons, blood vessels, nerves and bone components of the left arm and shoulder are shown.


Regular is $\frac{3}{4}$ life-size: 24" x 7" x 7", 4 lb
Deluxe is life-size: 31" x 10" x 10", 7 lb



12-4556 regular arm, 6-part
12-4558 deluxe arm, 6-part

hand and foot skeleton with ligaments and muscles

Features bones, muscles, tendons, ligaments, nerves, arteries and veins of the hand, lower forearm and foot.



12-4522 hand skeleton, 4-part, 13" x 5" x 5"; 1lb
12-4524 foot skeleton, 6-part, 9" x 7" x 12"; 2lb

hand and foot skeletons with ligaments


All the anatomically important ligaments and tendons are shown on these two detailed models. Hand includes the interosseous membrane between the radius and ulna and removable flexor retinaculum. Foot includes the Achilles and peroneus longus tendons of the ankle.



12-4521 hand skeleton, 12" x 6" x 4"; 1 lb
12-4523 foot skeleton, 9" x 7" x 12"; 1 lb

joint models

High-quality, life-size joint models with ligaments show the anatomy and all possible physiological movements



12-4513 deluxe shoulder*
12-4509 classic shoulder
12-4519 miniature shoulder*
12-4515 deluxe knee
12-4511 classic knee
12-4518 miniature knee*
*not shown

12-4516 deluxe elbow
12-4512 classic elbow
12-4520 miniature elbow*
12-4514 deluxe hip
12-4510 classic hip
12-4517 miniature hip*

sports therapy joint models




Shows the individual muscles and color-coded muscular origins (red) and insertions (blue). Cartilage is represented in light blue. By removing the individual muscles, movements of the muscles can be shown in detail.



12-4525 sports shoulder with removable muscles, 5-part with rotator cuffs. 7" x 7" x 9"; 2 lb
12-4528 elbow joint with removable muscles, 8-part 9" x 16" x 20"; 2 lb
12-4527 knee joint with removable muscles, 12 part 13" x 7" x 7"; 2 lb
12-4526 hip joint with removable muscles, 6-part 7" x 13" x 7"; 2 lb

Orfit® Classic sheets

- Non-coated material sticks easily to itself and to outrigger attachments
- Temporarily bonds during fabrication; permanent bonding possible with dry heat
- Turns transparent when heated for visibility
- Unlimited potential for fabrication of very large to very small orthoses
- Available in soft (high stretch) and stiff (low stretch)
- Soft version offers increased, accurate molding around anatomic corners; stiff version offers increases support for larger orthoses



perf	stretch	sheet	case qty
solid	soft	1/16" x 18" x 24"	24-5620-1 24-5620-4 4 ea
solid	soft	1/12" x 18" x 24"	24-5622-1 24-5622-4 4 ea
solid	soft	1/8" x 18" x 24"	24-5628-1 24-5628-4 4 ea
solid	soft	3/16" x 24" x 36"	24-5634-1 24-5634-4 4 ea
solid	stiff	1/8" x 18" x 24"	24-5636-1 24-5636-4 4 ea
solid	stiff	1/8" x 24" x 36"	24-5638-1 24-5638-4 4 ea
solid	stiff	3/16" x 24" x 36"	24-5640-1 24-5640-4 4 ea
25%	soft	1/12" x 18" x 24"	24-5625-1 24-5625-4 4 ea
25%	soft	1/8" x 18" x 24"	24-5629-1 24-5629-4 4 ea
3.5%	soft	1/8" x 18" x 24"	24-5630-1 24-5630-4 4 ea
3.5%	soft	3/16" x 24" x 36"	24-5635-1 24-5635-4 4 ea
3.5%	stiff	1/8" x 24" x 36"	24-5639-1 24-5639-4 4 ea
13%	soft	1/16" x 18" x 24"	24-5621-1 24-5621-4 4 ea
13%	soft	1/12" x 18" x 24"	24-5623-1 24-5623-4 4 ea

Orfit® Classic precuts

- Make splint creation quick and easy
- Precut pieces are anatomically shaped
- Easily remoldable
- Available in 3 sizes: small, medium and large



perf	thick	gauntlet thumb post	1 piece	case qty
	13%	1/16" small	24-5650-1	24-5650-2 2 ea
	13%	1/16" medium	24-5651-1	24-5651-2 2 ea
	13%	1/16" large	24-5652-1	24-5652-2 2 ea
	13%	1/12" small	24-5653-1	24-5653-2 2 ea
	13%	1/12" medium	24-5654-1	24-5654-2 2 ea
	13%	1/12" large	24-5655-1	24-5655-2 2 ea
	13%	1/12" small	24-5656-1	24-5656-2 2 ea
	13%	1/12" medium	24-5657-1	24-5657-2 2 ea
	13%	1/12" large	24-5658-1	24-5658-2 2 ea
	13%	1/12" small	24-5659-1	24-5659-2 2 ea
	13%	1/12" medium	24-5660-1	24-5660-2 2 ea
	13%	1/12" large	24-5661-1	24-5661-2 2 ea
	solid	1/8" small	24-5666-1	24-5666-2 2 ea
	solid	1/8" medium	24-5667-1	24-5667-2 2 ea
	solid	1/8" large	24-5668-1	24-5668-2 2 ea
	solid	1/8" small	24-5669-1	24-5669-2 2 ea
	solid	1/8" medium	24-5670-1	24-5670-2 2 ea
	solid	1/8" large	24-5671-1	24-5671-2 2 ea
	3.5%	1/8" small	24-5672-1	24-5672-2 2 ea
	3.5%	1/8" medium	24-5673-1	24-5673-2 2 ea
	3.5%	1/8" large	24-5674-1	24-5674-2 2 ea
	3.5%	1/8" small	24-5675-1	24-5675-2 2 ea
	3.5%	1/8" medium	24-5676-1	24-5676-2 2 ea
	3.5%	1/8" large	24-5677-1	24-5677-2 2 ea
	solid	1/8" left	24-5678L-1	
	solid	1/8" right	24-5678R-1	

Orfit® NS Soft sheets

- Coated material is non-stick, designed to prevent accidental bonding
- Allows for temporary bonding during splint creation
- Permanent bonding possible when coating removed and heat is applied
- Use as a base for different attachments after coating is removed
- Soft version is high stretch and allows for increased accurate molding around anatomic corners

perf	sheet	case qty
solid	1/16" x 18" x 24"	24-5680-1 24-5680-4 4 ea
solid	1/12" x 18" x 24"	24-5682-1 24-5682-4 4 ea
solid	1/8" x 18" x 24"	24-5688-1 24-5688-4 4 ea
25%	1/12" x 18" x 24"	24-5685-1 24-5685-4 4 ea
25%	1/8" x 18" x 24"	24-5689-1 24-5689-4 4 ea
3.5%	1/12" x 18" x 24"	24-5686-1 24-5686-4 4 ea
3.5%	1/8" x 18" x 24"	24-5690-1 24-5690-4 4 ea
13%	1/16" x 18" x 24"	24-5681-1 24-5681-4 4 ea
13%	1/12" x 18" x 24"	24-5683-1 24-5683-4 4 ea
13%	3/32" x 18" x 24"	24-5687-1 24-5687-4 4 ea



Orfit® Colors NS sheets

- Non-stick coating allows for temporary bonding during application
- Choice of unique colors
- Designed to appeal to patients of all ages
- Material has a high gloss and a shiny finish
- Available precuts make splint creation quick and easy



perf	sheet	case qty
solid sonic silver 1/12" x 18" x 24"	24-5770-1	24-5770-4 4 ea
solid sonic silver 1/8" x 18" x 24"	24-5772-1	24-5772-4 4 ea
solid atomic blue 1/12" x 18" x 24"	24-5774-1	24-5774-4 4 ea
solid atomic blue 1/8" x 18" x 24"	24-5776-1	24-5776-4 4 ea
solid gold 1/12" x 18" x 24"	24-5778-1	24-5778-4 4 ea
solid gold 1/8" x 18" x 24"	24-5780-1	24-5780-4 4 ea
solid violet 1/12" x 18" x 24"	24-5782-1	24-5782-4 4 ea
solid hot green 1/12" x 18" x 24"	24-5784-1	24-5784-4 4 ea
solid dominant black 1/12" x 18" x 24"	24-5786-1	24-5786-4 4 ea
solid dynamic red 1/12" x 18" x 24"	24-5788-1	24-5788-4 4 ea
solid ocean blue 1/12" x 18" x 24"	24-5790-1	24-5790-4 4 ea
3.5% sonic silver 1/8" x 18" x 24"	24-5773-1	24-5773-4 4 ea
3.5% atomic blue 1/8" x 18" x 24"	24-5777-1	24-5777-4 4 ea
3.5% gold 1/8" x 18" x 24"	24-5781-1	24-5781-4 4 ea
13% sonic silver 1/12" x 18" x 24"	24-5771-1	24-5771-4 4 ea
13% atomic blue 1/12" x 18" x 24"	24-5775-1	24-5775-4 4 ea
13% gold 1/12" x 18" x 24"	24-5779-1	24-5779-4 4 ea
13% violet 1/12" x 18" x 24"	24-5783-1	24-5783-4 4 ea
13% hot green 1/12" x 18" x 24"	24-5785-1	24-5785-4 4 ea
13% dominant black 1/12" x 18" x 24"	24-5787-1	24-5787-4 4 ea
13% dynamic red 1/12" x 18" x 24"	24-5789-1	24-5789-4 4 ea
13% ocean blue 1/12" x 18" x 24"	24-5791-1	24-5791-4 4 ea

Orfit® Colors NS precuts

perf thick

wrist + thumb



wrist + thumb



gauntlet immobilization

13% 1/12" atomic blue	24-5818-1
13% 1/12" ocean blue	24-5806-1
13% 1/12" sonic silver	24-5830-1
13% 1/12" dominant black	24-5842-1

gauntlet immobilization

13% 1/12" atomic blue	24-5821-1
13% 1/12" ocean blue	24-5809-1
13% 1/12" sonic silver	24-5833-1
13% 1/12" dominant black	24-5845-1

small

medium

large

24-5819-1	24-5820-1
24-5807-1	24-5808-1
24-5831-1	24-5832-1
24-5843-1	24-5844-1
24-5822-1	24-5823-1
24-5810-1	24-5811-1
24-5834-1	24-5835-1
24-5846-1	24-5847-1

Orfilight® sheets

- Non coated material is ideal for creating lightweight, supportive splints
- Easy to trim edges and soft foamy feel for patient comfort



perf	sheet	case qty
solid 1/8" x 18" x 24"	24-5722-1	24-5722-4 4 ea
3.5% 1/8" x 18" x 24"	24-5723-1	24-5723-4 4 ea
13% 1/16" x 18" x 24"	24-5720-1	24-5720-4 4 ea
13% 3/32" x 18" x 24"	24-5721-1	24-5721-4 4 ea

Orfilight® precuts



perf	thick	1 piece	case qty
rheumatoid arthritis resting			
3.5%	1/8" small	24-5940-1	24-5940-2 2 ea
3.5%	1/8" medium	24-5941-1	24-5941-2 2 ea
3.5%	1/8" large	24-5942-1	24-5942-2 2 ea

Orfilight® NS sheets

- Non-stick, coated material is ideal for creating lightweight, supportive splints
- Easy to trim edges and soft foamy feel for patient comfort
- Available in black and atomic blue colors



perf	sheet	case qty
solid black 1/8" x 18" x 24"	24-5742-1	24-5742-4 4 ea
3.5% black 1/8" x 18" x 24"	24-5743-1	24-5743-4 4 ea
13% black 1/16" x 18" x 24"	24-5740-1	24-5740-4 4 ea
13% black 3/32" x 18" x 24"	24-5741-1	24-5741-4 4 ea
solid atomic blue 1/8" x 18" x 24"	24-5762-1	24-5762-4 4 ea
3.5% atomic blue 1/8" x 18" x 24"	24-5763-1	24-5763-4 4 ea
13% atomic blue 1/16" x 18" x 24"	24-5760-1	24-5760-4 4 ea
13% atomic blue 3/32" x 18" x 24"	24-5761-1	24-5761-4 4 ea

Orficast™ thermoplastic tape

- Ideal for creating quick and easy immobilization splints
- Tape features thermoplastic threads knitted with cotton weave
- Feels soft, yet comfortable on client
- Offers full support and immobilization
- Product will adhere to self when warmed, even more so with dry heat
- **Orficast™ More** has more thermoplastic threads and is more rigid; allows for more orthotic fabrication options



Orficast™ (0.05" thick)

		roll
1/19" x 1" x 9.8 ft	blue	24-5600-1
1/19" x 1" x 9.8 ft	black	24-5602-1
1/19" x 2.5" x 9.8 ft	blue	24-5601-1
1/19" x 2.5" x 9.8 ft	black	24-5603-1

Orficast™ More (0.10" thick)

		roll
7/64" x 2.5" x 9.8 ft	blue	24-5610-1
7/64" x 2.5" x 9.8 ft	black	24-5613-1
7/64" x 5" x 9.8 ft	blue	24-5611-1
7/64" x 5" x 9.8 ft	black	24-5614-1
7/64" x 6" x 9.8 ft	blue	24-5612-1
7/64" x 6" x 9.8 ft	black	24-5615-1

case qty

24-5600-12	12 ea
24-5602-12	12 ea
24-5601-6	6 ea
24-5603-6	6 ea
24-5610-6	6 ea
24-5613-6	6 ea
24-5611-6	6 ea
24-5614-6	6 ea
24-5612-6	6 ea
24-5615-6	6 ea

Orfit® Natural NS Soft sheets

- Non-stick, anti-bacterial coating
- Ideal for hand splints
- Features a smooth white finish



perf		sheet	case qty
solid	1/16" x 18" x 24"	24-5710-1	24-5710-4 4 ea
solid	3/32" x 18" x 24"	24-5712-1	24-5712-4 4 ea
solid	1/8" x 18" x 24"	24-5714-1	24-5714-4 4 ea
25%	1/8" x 18" x 24"	24-5716-1	24-5716-4 4 ea
1%	1/8" x 18" x 24"	24-5715-1	24-5715-4 4 ea
13%	1/16" x 18" x 24"	24-5711-1	24-5711-4 4 ea
13%	3/32" x 18" x 24"	24-5713-1	24-5713-4 4 ea

Orfit® Flex NS sheets

- Non-stick, coated material
- Features high drapability and low resistance to stretch
- Anti-bacterial coating



perf		sheet	case qty
solid	1/8" x 18" x 24"	24-5902-1	24-5902-4 4 ea
solid	1/8" x 24" x 36"	24-5904-1	24-5904-4 4 ea
solid	3/32" x 18" x 24"	24-5900-1	24-5900-4 4 ea
1%	1/8" x 18" x 24"	24-5903-1	24-5903-4 4 ea
1%	1/8" x 24" x 36"	24-5905-1	24-5905-4 4 ea
1%	3/32" x 18" x 24"	24-5901-1	24-5901-4 4 ea

Orfit® Eco Black NS sheets

- Non-stick, coated black material is designed for large and rigid splints
- Made of high quality recycled materials
- Appealing smooth finish resists fingerprints
- Does not stick to bandages or dressings



perf		sheet	case qty
solid	3/32" x 18" x 24"	24-5870-1	24-5870-4 4 ea
solid	1/8" x 18" x 24"	24-5872-1	24-5872-4 4 ea
solid	1/8" x 24" x 36"	24-5874-1	24-5874-4 4 ea
3.5%	3/32" x 18" x 24"	24-5871-1	24-5871-4 4 ea
3.5%	1/8" x 18" x 24"	24-5873-1	24-5873-4 4 ea
3.5%	1/8" x 24" x 36"	24-5875-1	24-5875-4 4 ea

Orfit® Eco Black NS precuts

- Make splint creation quick and easy



perf		1 piece	case qty
Intrinsic Anti-Spastic Hand			
3.5%	1/8" small	24-5945-1	24-5945-2 2 ea
3.5%	1/8" medium	24-5946-1	24-5946-2 2 ea
3.5%	1/8" large	24-5947-1	24-5947-2 2 ea

Orfit® NS strips

- Strips are made of non-stick, coated material to prevent accidental bonding
- Use for corrective finger and hand orthoses
- Adhere to other orthoses to increase rigidity

1/12" x 18" x 1/5"
1/12" x 18" x 1/5"
1/12" x 18" x 1/5"
1/8" x 18" x 4/5"
1/8" x 18" x 4/5"
1/8" x 18" x 4/5"

sonic silver (10 pc) 24-5850-1
gold (10 pc) 24-5851-1
sonic silver (5 pc), gold (5 pc) 24-5852-1
gold (10 pc) 24-5853-1
atomic blue (10 pc) 24-5854-1
gold (5 pc), atomic blue (5 pc) 24-5855-1



Orfibrace® NS sheets

- Non-stick, coated material maintains perfect shape form and size for increased client comfort
- Anti-bacterial coating
- Cost effective when splints need to be remade



perf	sheet	case qty
solid 1/8" x 24" x 36"	24-5880-1	24-5880-4 4 ea
1% 1/8" x 24" x 36"	24-5881-1	24-5881-4 4 ea

Orfit® Allfit

- Small pellets of thermoplastic material
- Easily mold into any desired shape or pattern



bag	box of 6 bags
1 bag, 1.1 lbs (500 gm) 24-5910-1	24-5910-6 6 ea

Orfit® padding and liner Luxofoam

- Self-adhesive, soft PVC foam used for padding small surfaces
- Measures: 19.3" x 38.2" x 0.11"



24-5929 luxofoam, non perforated
24-5930 luxofoam, perforated

Orfilastic

- Terryloth lining with self-adhesive open cell foam backing
- Measures: 20" x 39.4" x 0.11"



24-5931 orfilastic
24-5932 orfilastic black

Thermoplastic and Elastic edging strip

- Use as an edging material on splints, external immobilization devices and rehabilitation aids
- Measures: 0.04" x 0.75" x 118"



24-5933 thermoplastic, natural color
24-5934 elastic soft velvet, black

Orfitube™ outrigger attachments

- Attach outriggers to splinting material to create mobilization orthoses (not shown)
- Secure to any type of base splint
- Precoiled wiring saves time in energy in the creation of more difficult splints

24-5912 orfitubes (10 pc), bending wires (2 pcs)
24-5914 elastic thread, max resistance 7 oz, 30 ft
24-5915 elastic thread, max resistance 14 oz, 30 ft
24-5916 elastic thread, max resistance 21 oz, 30ft
24-5917 orfitube ends (20 pc)
24-5918 single orfitube adjuster
24-5919 single orfitube adjuster (5 pc), allen wrench
24-5920 orfitube bending tool
24-5921 finger coil spring (10 pc)
24-5922 knuckle bender coil spring (10 pc)
24-5923 finger extension assist (10 pc)
24-5924 wrist extension coil spring (10 pc, 5 left & 5 right)
24-5928 stainless steel spring wire, 45 ft

- **Orfitube™ kit includes:**
10 orfitubes that are 330mm (13") long;
3 different types of elastic thread on a spool with different degrees of elastic strength;
1 spool of polyester monofilament for an optimal gliding effect; 2 bending wires

24-5911 orfitube kit

- **Radial Wrist Extension splint kit includes:**
2 Orfit® Classic Precuts (3.2mm (1/8") non perforated, stiff); and 2 lateral coil springs


24-5925 radial wrist extension splint kit, small
24-5926 radial wrist extension splint kit, medium
24-5927 radial wrist extension splint kit, large



Manosplint™ Ohio non-stick sheets

- Coated material has a shiny finish that resists fingerprints
- Will not adhere to self unless coating is removed
- 100% memory; maximum rigidity
- Moderate pressure needed for a conforming fit
- Material becomes translucent when heated for easy application
- Available in white and beige colors

PERFORATION PATTERNS



perf	sheet	case qty
solid 1/16" x 12" x 18" white	24-5210-1	24-5210-8 8 ea
solid 1/16" x 12" x 18" beige	24-5211-1	24-5211-8 8 ea
solid 1/16" x 18" x 24" white	24-5200-1	24-5200-4 4 ea
solid 1/16" x 18" x 24" beige	24-5201-1	24-5201-4 4 ea
solid 1/12" x 12" x 18" beige	24-5230-1	24-5230-8 8 ea
solid 1/12" x 18" x 24" beige	24-5220-1	24-5220-4 4 ea
solid 3/32" x 12" x 18" white	24-5250-1	24-5250-8 8 ea
solid 3/32" x 18" x 24" white	24-5240-1	24-5240-4 4 ea
solid 3/32" x 18" x 24" beige	24-5241-1	24-5241-4 4 ea
solid 1/8" x 12" x 18" beige	24-5280-1	24-5280-8 8 ea
solid 1/8" x 12" x 18" white	24-5281-1	24-5281-8 8 ea
solid 1/8" x 18" x 24" beige	24-5270-1	24-5270-4 4 ea
solid 1/8" x 18" x 24" white	24-5271-1	24-5271-4 4 ea
solid 1/8" x 24" x 36" white	24-5260-1	24-5260-2 2 ea
1% 1/8" x 12" x 18" white	24-5282-1	24-5282-8 8 ea
1% 1/8" x 18" x 24" white	24-5272-1	24-5272-4 4 ea
1% 1/8" x 24" x 36" white	24-5261-1	24-5261-2 2 ea
11% 1/16" x 12" x 18" white	24-5212-1	24-5212-8 8 ea
11% 1/16" x 12" x 18" beige	24-5213-1	24-5213-8 8 ea
11% 1/16" x 18" x 24" white	24-5202-1	24-5202-4 4 ea
11% 1/16" x 18" x 24" beige	24-5203-1	24-5203-4 4 ea
15% 1/12" x 12" x 18" beige	24-5231-1	24-5231-8 8 ea
15% 1/12" x 18" x 24" beige	24-5221-1	24-5221-4 4 ea
15% 3/32" x 12" x 18" beige	24-5251-1	24-5251-8 8 ea
15% 3/32" x 12" x 18" white	24-5252-1	24-5252-8 8 ea
15% 3/32" x 18" x 24" beige	24-5242-1	24-5242-4 4 ea
15% 3/32" x 18" x 24" white	24-5243-1	24-5243-4 4 ea
19% 1/8" x 12" x 18" white	24-5283-1	24-5283-8 8 ea
19% 1/8" x 18" x 24" white	24-5273-1	24-5273-4 4 ea
36% 3/32" x 12" x 18" white	24-5253-1	24-5253-8 8 ea
36% 3/32" x 18" x 24" white	24-5244-1	24-5244-4 4 ea
42% 1/8" x 12" x 18" white	24-5284-1	24-5284-8 8 ea
42% 1/8" x 18" x 24" white	24-5274-1	24-5274-4 4 ea

Manosplint™ Ohio S sticky sheets

- Non-coated material allows for self-bonding, keeping material in place during application
- 100% memory; maximum rigidity
- Moderate pressure needed for a conforming fit
- Material becomes translucent when heated for easy application



perf	sheet	case qty
solid 1/16" x 12" x 18" beige	24-5410-1	24-5410-8 8 ea
solid 1/16" x 18" x 24" beige	24-5400-1	24-5400-4 4 ea
solid 3/32" x 12" x 18" beige	24-5430-1	24-5430-8 8 ea
solid 3/32" x 18" x 24" beige	24-5420-1	24-5420-4 4 ea
solid 1/8" x 12" x 18" beige	24-5450-1	24-5450-8 8 ea
solid 1/8" x 18" x 24" beige	24-5440-1	24-5440-4 4 ea
11% 1/16" x 12" x 18" beige	24-5411-1	24-5411-8 8 ea
11% 1/16" x 18" x 24" beige	24-5401-1	24-5401-4 4 ea
15% 3/32" x 12" x 18" beige	24-5431-1	24-5431-8 8 ea
15% 3/32" x 18" x 24" beige	24-5421-1	24-5421-4 4 ea
19% 1/8" x 12" x 18" beige	24-5451-1	24-5451-8 8 ea
19% 1/8" x 18" x 24" beige	24-5441-1	24-5441-4 4 ea

Manosplint™ Ohio F fabric sheets

- Two-sided fabric-covered material makes fabrication efficient with comfort built-in
- Anti-microbial fabric reduces risk of bacteria and odor
- Unique hook-receptive fabric makes strapping easy to apply




perf	sheet	case qty
1% 1/8" x 12" x 18" blue/gray	24-5350-1	24-5350-8 8 ea
1% 1/8" x 12" x 18" beige	24-5351-1	24-5351-8 8 ea
1% 1/8" x 18" x 24" blue/gray	24-5340-1	24-5340-4 4 ea
1% 1/8" x 18" x 24" beige	24-5341-1	24-5341-4 4 ea
5% 1/8" x 12" x 18" blue/gray	24-5352-1	24-5352-8 8 ea
5% 1/8" x 12" x 18" beige	24-5353-1	24-5353-8 8 ea
5% 1/8" x 18" x 24" blue/gray	24-5342-1	24-5342-4 4 ea
5% 1/8" x 18" x 24" beige	24-5343-1	24-5343-4 4 ea
36% 1/16" x 12" x 18" blue/gray	24-5300-1	24-5300-8 8 ea
36% 1/16" x 12" x 18" beige	24-5301-1	24-5301-8 8 ea
36% 1/16" x 18" x 24" blue/gray	24-5310-1	24-5310-4 4 ea
36% 1/16" x 18" x 24" beige	24-5311-1	24-5311-4 4 ea
36% 3/32" x 12" x 18" blue/gray	24-5330-1	24-5330-8 8 ea
36% 3/32" x 12" x 18" beige	24-5331-1	24-5331-8 8 ea
36% 3/32" x 18" x 24" blue/gray	24-5320-1	24-5320-4 4 ea
36% 3/32" x 18" x 24" beige	24-5321-1	24-5321-4 4 ea

Manosplint™ Wisconsin non-stick sheets

- Coated material resists fingerprints
- Material will not adhere to itself unless the coating is removed
- Allows for easy modification and edge finishing
- Confirming and sculpted finish provides great support and function



PERFORATION PATTERNS



perf	sheet	case qty
solid 3/32" x 12" x 18" white	24-5540-1	24-5540-8 8 ea
solid 3/32" x 18" x 24" white	24-5530-1	24-5530-4 4 ea
solid 1/8" x 12" x 18" blue	24-5520-1	24-5520-8 8 ea
solid 1/8" x 12" x 18" beige	24-5521-1	24-5521-8 8 ea
solid 1/8" x 12" x 18" white	24-5522-1	24-5522-8 8 ea
solid 1/8" x 18" x 24" blue	24-5510-1	24-5510-4 4 ea
solid 1/8" x 18" x 24" beige	24-5511-1	24-5511-4 4 ea
solid 1/8" x 18" x 24" white	24-5512-1	24-5512-4 4 ea
solid 1/8" x 24" x 36" blue	24-5500-1	24-5500-2 2 ea
solid 1/8" x 24" x 36" beige	24-5501-1	24-5501-2 2 ea
solid 1/8" x 24" x 36" white	24-5502-1	24-5502-2 2 ea
1% 3/32" x 12" x 18" white	24-5541-1	24-5541-8 8 ea
1% 3/32" x 18" x 24" white	24-5531-1	24-5531-4 4 ea
1% 1/8" x 12" x 18" blue	24-5523-1	24-5523-8 8 ea
1% 1/8" x 12" x 18" beige	24-5524-1	24-5524-8 8 ea
1% 1/8" x 12" x 18" white	24-5525-1	24-5525-8 8 ea
1% 1/8" x 18" x 24" blue	24-5513-1	24-5513-4 4 ea
1% 1/8" x 18" x 24" beige	24-5514-1	24-5514-4 4 ea
1% 1/8" x 18" x 24" white	24-5515-1	24-5515-4 4 ea
1% 1/8" x 24" x 36" blue	24-5503-1	24-5503-2 2 ea
1% 1/8" x 24" x 36" beige	24-5504-1	24-5504-2 2 ea
1% 1/8" x 24" x 36" white	24-5505-1	24-5505-2 2 ea

CanDo® Eddy® heat gun

- Use heat gun to warm and mold thermoplastic material
- Variable temperature control
- Air concentrators allow heat gun to direct air flow to a specific area
- Air spreader allows airflow to be directed to a larger area
- Heat gun kit includes: heat gun, 3/8" air concentrator, 3" air spreader, and case



24-4070	heat gun kit
24-4071	heat gun only

accessories

24-4072	heat gun attachment, 7/8" Air Concentrator
24-4073	heat gun attachment, 3/8" Air Concentrator
24-4074	heat gun attachment, 2.25" Air Spreader
24-4075	heat gun attachment, 3" Air Spreader
24-4076	heat gun stand

Manosplint™ Carolina sticky sheets

- Use for large orthoses
- Designed for patients who are unable to assist with positioning
- Offers exceptional control during the fabrication process
- Non-coated material is slightly sticky and stays in place during application
- Resists fingerprints



PERFORATION PATTERNS



perf	sheet	case qty
solid 1/8" x 12" x 18" white	24-5120-1	24-5120-8 8 ea
solid 1/8" x 18" x 24" white	24-5110-1	24-5110-4 4 ea
solid 1/8" x 24" x 36" white	24-5100-1	24-5100-2 2 ea
1% 1/8" x 12" x 18" white	24-5121-1	24-5121-8 8 ea
1% 1/8" x 18" x 24" white	24-5111-1	24-5111-4 4 ea
1% 1/8" x 24" x 36" white	24-5101-1	24-5101-2 2 ea

Manosplint™ dry heat oven


- Use to heat thermoplastic splinting material
- Features water-free technology and digital temperature controls

24-5550 dry oven



Forma-Splint® splint pans

- Use to heat thermoplastic splinting material
- Large and medium pans have an insulated lid and dual thermostatic controls
- Small pan has durable stainless steel construction and a fast heat up time



			max sheet
24-4090	large	SP-1501	24" x 18"
24-4091	medium	SP-1502	19" x 11"
24-4094	small	LS-1	16" x 10"

Fiskars® splinting scissors

- Perfect for cutting thermoplastic materials
- Heavy duty scissors have knife edge and rounded tips for safety
- 8" length

24-6500	heavy duty
24-6501	left-handed
24-6502	right-handed



nylon spatula

- Use to safely lift warmed thermoplastic material out of heating pans
- Dishwasher safe; heat resistant up to 450°F

24-4060 spatula



splints and supports

Uriel® arm & shoulder support sleeve

- Provides compression of the common extensor muscles
- Aids with pain relief
- Fits left or right arm / shoulder



24-9062	medium
24-9063	large
24-9064	x-large
24-9065	xx-large

Uriel® elbow compression sleeve

- Provides uniform compression and support around the elbow
- Assists reducing swelling and pain
- Fits left or right elbow



24-9072	medium
24-9073	large
24-9074	x-large
24-9075	xx-large

Uriel® thumb splint



- Designed to support the thumb
- Material is breathable and comfortable
- Recommended for injuries or use post-cast removal
- Fits left or right hand

24-9002	medium
24-9003	large
24-9004	x-large
24-9005	xx-large

Uriel® rigid thumb splint



- Designed to provide firm support for the thumb
- Contains a flexible spiral rod that limits thumb movement
- Fits left or right hand

24-9019	universal size
---------	----------------

Uriel® wrist & thumb splint

- Holds the wrist in a functional position (35°)
- Restricts thumb movement and supports the wrist
- Fits left or right wrist / thumb



24-9029	universal size
---------	----------------

Uriel® neoprene max wrist support

- Provides firm wrist support while allowing for normal daily function
- Holds wrist at a functional 35° angle



24-9039L	universal, left
24-9039R	universal, right

Uriel® neoprene wrist support



- Use to treat wrist pain
- Allows for normal daily function
- Fits left or right wrist

24-9049	universal size
---------	----------------

Uriel® thermo wrist support

- Provides firm support to minimize strain on wrist
- Allows for normal daily function
- Fits left or right wrist



24-9059	universal size
---------	----------------

Comfy Splints™ hand/wrist/thumb splint

- Adjustable splint. Fits left or right hand / wrist / thumb
- **Hand/wrist**; side wings prevent ulna and radial deviation
- **Hand/thumb**; supports the thumb without stressing the web space between index finger and thumb



hand/
wrist



hand/
thumb

hand/wrist		hand/thumb	
24-3096	adult, small	24-3113	adult, small
24-3090	adult, medium	24-3110	adult, medium
24-3093	adult, large	24-3116	adult, large
24-3099	pediatric, small	24-3119	pediatric, small
24-3102	pediatric, medium	24-3122	pediatric, medium
24-3105	pediatric, large	24-3125	pediatric, large

Comfy Splints™ C-Grip hand splint



- Supports the wrist and uses soft finger rolls to gradually encourage finger extension
- Includes 2 graduated finger rolls

left hand		right hand	
24-3040L	adult, small	24-3040R	adult, small
24-3041L	adult, medium	24-3041R	adult, medium
24-3042L	adult, large	24-3042R	adult, large

Comfy Splints™ finger separator

- Soft cushion between fingers protects and prevents pressure areas between the PIP and DIP joints of adjoining digits
- Can be used faced downward to prevent hyperextension
- Fits left or right hand



24-3131	adult, small	24-3132	pediatric, small
24-3130	adult, medium	24-3133	pediatric, medium
		24-3134	pediatric, large

AFH wrist / thumb support



- Provides pleasant compression and support for the hand and thumb joint
- Fits left or right hand

24-6320	wrist and thumb support, velcro, deluxe ambidextrous
---------	--

kinesiology tape and IASTM tools

CanDo® kinesiology tape

- Use for kinesiology taping treatments of muscle fascia and tendinitis; improves circulation
- Made from cotton / rayon fiber; latex-free and hypoallergenic
- Wave-like acrylic adhesive is sweat resistant
- Tape is water-resistant and air-permeable allowing it to be worn for days, even in the shower, pool or sauna
- Can be stretched lengthways by 30-40%
- Available in 4 colors and two sizes
- Larger size includes easy-to-dispense box



2" x 16.5 ft roll

beige	24-4850	24-4850-10
blue	24-4851	24-4851-10
black	24-4852	24-4852-10
pink	24-4853	24-4853-10

case of 10

2" x 103 ft roll

beige	24-4860
blue	24-4861
black	24-4862
pink	24-4863



Free-Up® massage cream



- Formulated for healthcare professionals
- Noted for its glide, lubricity and tissue perception
- Slow absorption rate
- Bacteriostatic and hypo allergenic; does not contain beeswax

13-3244	8 oz jar
13-3241	16 oz jar

HawkGrips® soft tissue therapy tools

- Features convex and concave treatment edges
- Textured body prevents tool from slipping during use
- Use to treat a variety of body parts and sizes
- Perfect for clinicians on-the-go
- HawkGrips® Pro Multi-Instrument includes: drawstring pouch, 1 jar of regular emollient, and user manual
- Emollient helps HawkGrips® instruments glide smoothly over the skin while softening and protecting skin
- Sets include: a selection of HawkGrips® tools, roll up carrying case, 1 jar of regular emollient and user manual

Instrument

13-4023 HawkGrips® Pro Multi-Instrument

Sets

13-4000	introductory 3 instrument set
13-4001	silver 6 instrument set
13-4002	gold 7 instrument set
13-4003	platinum 9 instrument set

Emollients (Creams)

regular fragrance free

1 jar	13-4030	13-4031
5 jars	13-4030-5	13-4031-5
10 jars	13-4030-10	13-4031-10
20 jars	13-4030-20	13-4031-20



Introductory set

- Small Multi-Curve
- Scanner
- "Tongue Depressor"

Silver Set

- Small Multi-Curve
- Medium Multi-Curve
- Large Multi-Curve
- Scanner
- Boomerang
- "Tongue Depressor"

Gold Set

- Medium Handlebar
- Small Multi-Curve
- Medium Multi-Curve
- Large Multi-Curve
- Scanner
- Boomerang
- "Tongue Depressor"

Platinum Set

- Small Handlebar
- Medium Handlebar
- Large Handlebar
- Small Multi-Curve
- Medium Multi-Curve
- Large Multi-Curve
- Scanner
- Boomerang
- "Tongue Depressor"



massage furniture and percussor machines

blue
(B)

burgundy
(BUR)

black
(BLK)

green
(G)

massage furniture colors

deluxe table



- High quality portable massage table
- Vinyl covered 3" foam top
- Includes adjustable face cradle, hanging arm hammock and carrying case
- Rubber feet prevent scratching and sliding on floor
- 76"L x 30"W x 23" – 33"H; 35.5 lb
- 550 lb capacity

15-3732B blue
15-3732BUR burgundy
15-3732BLK black
15-3732G green



basic portable massage table



- Economical portable massage table made of birch wood
- 2" dense foam covered with oil- and water-proof vinyl
- Includes adjustable face cradle, side arm extenders, hanging arm hammock and carrying case
- Rubber feet on legs prevents scratching and sliding on floor
- 78"L x 32"W x 23" – 33"H; 35.5 lb
- 550 lb capacity

15-3731B blue
15-3731BUR burgundy
15-3731BLK black
15-3731G green



basic stationary table



- Flat top table
- Manually adjustable height from 25" to 32.7"
- 2.25" thick vinyl covered top offers patient comfort
- Features open shelf storage and hanging arm rests
- 74"L x 27"W x 25" - 32.7"H; 70.4 lbs
- 600 lb weight capacity

15-3740B blue
15-3740BLK black

table with adjustable back

- Multi-purpose portable table
- Tilt back adjusts to 10 different levels
- 73"L x 30"W x 24" – 32"H; 34 lb
- 450 lb capacity



15-3733B blue
15-3733BUR burgundy

aluminum massage table

- Lightweight, portable massage table
- Features aluminum headrest frame, armrest shelf and carrying case
- Height adjustable from 24" to 33"
- Foam top enhances patient comfort
- 73"L x 28"W x 24" - 33"H; 28 lbs
- 350 lb weight capacity



15-3742B blue
15-3742BLK black

massage chair

- Portable ergonomic design maintains proper body alignment
- Seat height, sternum and headrest are adjustable
- Includes carrying bag
- 500 lb capacity



15-3730B blue
15-3730BUR burgundy
15-3730BLK black

G5® TherAssist

- Ideal for stress reduction, myofascial release, sports rehabilitation, postural drainage
- Includes pro-pack applicators and trolley
- Can be used on table-top or wall-mounted (hardware sold separately)
- Speed 15-60 CPS



TherAssist 14-1503

G5® Flimm Fighter

- Ideal for COPD, Cystic Fibrosis, Chronic Bronchitis and Asthma
- Speed 15-30 CPS
- Tabletop mounted



for adult 14-1520
pediatric / 14-1521
geriatric

G5® hand-held portable massagers

- Speed 15-55 CPS
- Ideal for light duty clinics and homes
- Includes pro-port pack applicators
- GBM® measures: 10"L x 4"W x 4"H; 4 lb
- ProPower® measures: 12"L x 3"W x 3"H; 3 lb



GBM® 14-1501
ProPower® 14-1506

G5® Vibracare®

- Portable hand-held percussor
- Comfortable soft rubber hand grip
- Speed 15-55 CPS
- Adjunct for in-home postural drainage
- Includes therapy package



for adult 14-1522
pediatric / 14-1523
geriatric

massage tools and hook and loop

CanDo® massage balls and rollers


- Achieve muscle and tension release for all parts of the body
- “Bumpy” textured surface promotes sensory stimulation and is perfect for manual dexterity exercises
- Rollers are effective for feet and large areas



ball size			each	bulk (12 each)
2.8"	7 cm	green	30-1995	30-1995-12
3.2"	8 cm	yellow	30-1996	30-1996-12
3.6"	9 cm	red	30-1997	30-1997-12
4.0"	10 cm	blue	30-1998	30-1998-12
6.0"	15 cm	yellow	30-1999	30-1999-12
roller size			pair	bulk (12 each)
2.6" x 6.3"		blue	30-1994	30-1994-12

Bongers® percussion massager

- Long flexible handle facilitates self-massage to lower back, upper shoulders, arms and legs, or to a partner
- Proprioception tool that strengthens wrists, arms and shoulders
- Adjustable strap allows Bongers® to be hung or worn



14-1410	blue pair
14-1411	purple pair
14-1412	cosmic pair

Point Relief® mini massager


- Comfortable vibrating trigger-point massage
- Available with or without heat
- Compact, portable and small enough to fit in purse or pocket
- Battery operated (1 “C” battery, included)



includes 4 massage tips



massager with heat includes 5th tip for HEAT massage




	each	25 per case
massager	14-1050	14-1050-25
massager with heat	14-1051	14-1051-25



hook or loop rolls

- **Standard** woven-back tape is sewn or glued onto object
- **Self-adhesive** backing is peeled off and tape is pressed onto an object
- Hook and loop sold separately; specify white (W), black (BLK), tan (T) when ordering


	Standard		Self-Adhesive	
	hook	loop	hook	loop
½" x 10 yd	24-7066	24-7067	24-7056	24-7057
1" x 10 yd	24-7002	24-7003	24-7012	24-7013
1½" x 10 yd	24-7006	24-7007	24-7026	24-7027
2" x 10 yd	24-7022	24-7023	24-7032	24-7033
½" x 25 yd	24-7064	24-7065	24-7054	24-7055
1" x 25 yd	24-7000	24-7001	24-7010	24-7011
1½" x 25 yd	24-7004	24-7005	24-7024	24-7025
2" x 25 yd	24-7020	24-7021	24-7030	24-7031



hand therapy tools

scar tissue massage sticks

- Massage sticks are used for breaking apart scar tissue and releasing trigger points
- 4 different tip sizes allow for more concentrated or dispersed pressure
- Massage sticks are available in stainless steel or gold plated
- Sold individually or in 2 piece sets



14-1430	stainless steel, w/box, very-fine
14-1431	stainless steel, w/box, fine
14-1432	stainless steel, w/box, medium
14-1433	stainless steel, w/box, large
14-1440	gold plated, w/box, very-fine
14-1441	gold plated, w/box, fine
14-1442	gold plated, w/box, medium
14-1443	gold plated, w/box, large
14-1434	2 piece set, stainless steel, w/box, very-fine/fine
14-1435	2 piece set, stainless steel, w/box, medium/large
14-1444	2 piece set, gold plated, w/box, very-fine/fine
14-1445	2 piece set, gold plated, w/box, medium/large
14-1436	protective case for massage sticks

massage roller

- For desensitizing and massage
- 4 different degrees of hardness:

yellow = soft
red = medium
green = hard
blue = extra-hard

- A few minutes of invigorating massage to the problem areas improve the blood flow and the skin image considerably

14-1450	massage roller w/4-pads (Y,R,G,B)
14-1452	4 additional roller pads (Y,R,G,B)



mini cupping kit

- When the cup with vacuum is placed on target area, the skin is lifted off over the muscles and organs. Fresh blood rushes to fill the space bringing with it new oxygen and nutrients.
- The small size cups are ideal for placement on hands and fingers.
- The cups are transparent to monitor the condition of the skin.



16-1640	mini cupping kit, standard
16-1641	mini cupping kit, sensitive
16-1642	mini cupping kit with dynamic pump

swan neck ring splint

- Use swan neck ring splints for mild to moderate hyperextension of the PIP or DIP joint
- Worn on the fingers or the thumb the splint blocks hyperextension
- While wearing the splint you will be able to maintain full flexion of the finger

24-6300	circumference 41mm
24-6301	circumference 44mm
24-6302	circumference 47mm
24-6303	circumference 50mm
24-6304	circumference 53mm
24-6305	circumference 57mm
24-6306	circumference 60mm
24-6307	circumference 63mm
24-6308	circumference 66mm
24-6309	circumference 69mm
24-6310	circumference 72mm
24-6311	circumference 75mm
24-6312	12 piece ring set w/storage box and template
24-6313	storage box
24-6314	ring template, 1-36mm



thumb helper massagers

- Use to perform deep connective tissue and trigger point massage
- Use with oil, lotions or creams

14-1460Y	thumb saver massager, yellow
14-1460G	thumb saver massager, green
14-1460B	thumb saver massager, blue
14-1460R	thumb saver massager, red
14-1460LB	thumb saver massager, light blue
14-1460GRY	thumb saver massager, gray
14-1460BLK	thumb saver massager, black



Thumbsavers® massagers

- Provides support and reduces stress on joints and wrists for deep tissue massage work
- Advance model features an advanced sensitive tip, improved venting and a soft flexible opening

Thumbsavers® Classic

14-1470	small, red
14-1471	medium, white
14-1472	large, blue

Thumbsavers® Advance

14-1473	small, teal
14-1474	large, orange



Point Relief® topical pain relief and WaxWel® paraffin treatment

point relief® coldspot™ topical analgesic



- Relieves minor muscle and joint discomfort caused by backache, arthritis, strains and/or sprains
- Penetrates deep into the muscle
- All-natural ingredients, including menthol, MSM, and arnica
- Water-based, clear formula minimizes grease and residue
- Helps treat bruising



	each	bulk	
4 oz gel tube	11-0730-1	11-0730-12	12 ea
4 oz hands-free applicator gel tube	11-0731-1	11-0731-12	12 ea
3 oz roll-on	11-0720-1	11-0720-12	12 ea
4 oz spray	11-0701-1	11-0701-12	12 ea
16 oz gel dispenser bottle	11-0710-1	11-0710-24	24 ea
32 oz gel dispenser bottle	11-0711-1	11-0711-16	16 ea
1 gallon gel dispenser bottle	11-0712-1	11-0712-4	4 ea

display box set

multi-pack display includes 4 ea: 4 oz gel, 4 oz spray, 3 oz roll-on	11-0765-12	12 ea
spray display includes 12 ea: 4 oz spray bottles	11-0760-12	12 ea
gel display includes 12 ea: 4 oz gel tube	11-0761-12	12 ea
roll-on display includes 12 ea: 3 oz roll-on	11-0762-12	12 ea

point relief® lidospot™ topical anesthetic

- Each patch contains the highest permissible amount of lidocaine without a prescription (4%) and 1% menthol
- Patch can be worn for up to 12 hours
- can be easily cut to fit hard-to-fit areas, such as the elbow
- Strong adhesive keeps patch in place even during strenuous activity



13-5310-5	5/pack
13-5310-20	20/pack

WaxWel® paraffin bath unit and wax

- Use for hand and wrist, foot and ankle, or elbow
- 6 lb (4 liter) unit with oversized tank accommodates a man's size 12½ foot
- Lightweight and portable; perfect for high traffic clinics or for home use
- Operating temperature 125-134°F, 110 Volt
- Kit includes unit, washable terry cloth mitt and bootie, 65 hand/foot plastic liners and (6) 1 lb wax blocks



11-1601 WaxWel® unit only

11-1600	WaxWel® unit with 6 lb unscented kit
11-1603	WaxWel® unit with 6 lb wintergreen kit
11-1604	WaxWel® unit with 6 lb lavender kit
11-1605	WaxWel® unit with 6 lb citrus kit
11-1606	WaxWel® unit with 6 lb peach kit
11-1607	WaxWel® unit with 6 lb rose kit

WaxWel® accessories for all paraffin units

11-1711	terry mitt (6 each)
11-1712	terry bootie (6 each)
11-1713	plastic hand/foot liners (100 each)
11-1710	set: mitt (1 each), bootie (1 each), liners (65 each)

WaxWel® paraffin blocks

	6 lb box	36 lb case
unscented	11-1720-6	11-1720-36
wintergreen	11-1722-6	11-1722-36
lavender	11-1714-6	11-1714-36
citrus (orange)	11-1716-6	11-1716-36
peach	11-1718-6	11-1718-36
rose	11-1730-6	11-1730-36

WaxWel® paraffin beads

	3 lb bucket	6 lb box	(6 - 6 lb boxes)
unscented	11-1750-3	11-1750-6	11-1750-36
wintergreen	11-1752-3	11-1752-6	11-1752-36
lavender	11-1744-3	11-1744-6	11-1744-36
citrus (orange)	11-1746-3	11-1746-6	11-1746-36
peach	11-1748-3	11-1748-6	11-1748-36
rose	11-1760-3	11-1760-6	11-1760-36

Therabath® paraffin bath

- 6-9 lb (4-5½ liter) paraffin capacity
- Hospital grade paraffin warmer with lid, grill and guide
- Operating temperature of 130°F



11-1171	unit with 6 lb unscented wax - 110 V
11-1171-220V	unit with 6 lb unscented wax - 220 V
11-1170	unit with 6 lb wintergreen wax - 110 V
11-1170-220V	unit with 6 lb wintergreen wax - 220 V

heat therapy

Hydrocollator® moist heat packs

- Provides up to 30 minutes of soothing, moist heat
- Brown tabs on moist packs for rotation in clinics
- Cost-effective; provides hundreds of treatments



packs only		each	bulk (12 each)
A	standard 10" x 12"	00-1006	00-1006-12
B	neck 5" x 24"	00-1002	00-1002-12
C	half size 5" x 12"	00-1014	00-1014-12
D	oversize 15" x 24"	00-1004	00-1004-12
E	knee-shoulder 10" x 20"	00-1012	00-1012-12
F	spine small 10" x 18"	00-1008	00-1008-12
G	spine large 10" x 24"	00-1010	00-1010-12
H	hand 6.5" x 12.5"	00-1026	00-1026-12
pack and cover sets		each	bulk (12 each)
standard pack and cover set		00-1064	00-1064-12
neck pack and cover set		00-1062	00-1062-12

Hydrocollator® stationary heating units and carts

- Stationary heating units are small enough to be kept on the countertop or cart
- Water bath keeps packs at 168°F



		units with packs	L x W x H
N	00-2102	E-1 with 4 standard	12½" x 8" x 15"
	00-2102-2	E-1 with 3 standard, 1 neck	12½" x 8" x 15"
	00-2102-3	E-1 with 2 standard, 2 neck	12½" x 8" x 15"
O	00-2802	E-2 with 2 OS, 3 std., 1 neck	14" x 12" x 16½"
		carts	
P	11-1390	cart (holds E-1 unit)	
	11-1391	cart (holds E-2 unit)	
Q	00-4002	2 shelf MB	18¾" x 16¾" x 19½"
R	00-4004	3 shelf MB-T	18¾" x 16¾" x 29½"
		accessories	
S	11-1397	hot pack tongs	
	00-4228	dial thermometer	

Theratherm® digital moist heat pads

- Unique digital hand control allows for adjustable treatment temperature and time
- Safe and easy to use to relieve pain with heat therapy



00-1030	small (7 x 15")
00-1031	medium (14 x 14")
00-1032	large (14 x 27")
00-1033	shoulder / neck (23 x 20")

Hydrocollator® hot pack covers



all terry covers		each	bulk (12 each)
I	standard 20" x 24"	00-1118	00-1118-12
J	neck 9" x 24"	00-1104	00-1104-12
K	oversize 24" x 30"	00-1126	00-1126-12
foam-filled terry covers			
L	standard 20" x 24"	00-1102	00-1102-12
	neck 9" x 24"	00-1120	00-1120-12
M	oversize 24" x 30"	00-1124	00-1124-12
	standard with pocket	00-1108	00-1108-12
	oversize with pocket	00-1110	00-1110-12
00-4016 wall-mounted towel rack (33" x 2" with six 13" hooks)			

Hydrocollator® mobile heating units

- Units are large enough for a complete department
- No plumbing required; simply fill with water and plug in
- Thermostatically maintains temperature of 168°F



units with packs		L x W x H
00-2302	SS-2 with 8 standard	20¾" x 15½" x 33"
00-2402	M-2 with 12 standard	26½" x 15½" x 33½"
00-2502	M-4 with 24 standard	35" x 20" x 33"

moist heat pads - analog and digital

- Provides moist heat without adding water (fomentation)
- Analog controller adjusts temperature
- Digital controller adjusts temperature and treatment time
- Includes flannel cover



analog pads		
11-1120	king	26" x 14"
11-1121	medium	18" x 14"
11-1123	small	19" x 7"
digital pads		
11-1135	king	26" x 14"
11-1136	medium	18" x 14"

Relief Pak® moist heat packs & covers

- Reusable
- Easy to heat in hot pack heating unit
- Maintains therapeutic moist heat for up to 30 minutes
- Use with Relief Pak® cover or towel
- Washable foam-filled terry cloth or all terry microfiber covers available for each size



foam-filled hot pack covers

	size	each	bulk (12 each)
standard	20" x 24"	11-1360	11-1360-12
neck	9" x 24"	11-1361	11-1361-12
oversize	24" x 36"	11-1362	11-1362-12
half size	12" x 15"	11-1363	11-1363-12
circular	12" dia	11-1226	11-1226-12
standard with pocket		11-1364	11-1364-12
oversize with pocket		11-1365	11-1365-12

all-terry covers

standard	20" x 24"	11-1366	11-1366-12
neck	9" x 24"	11-1367	11-1367-12
oversize	24" x 36"	11-1368	11-1368-12
half size	12" x 15"	11-1369	11-1369-12
circular	12" dia	11-1227	11-1227-12

moist heat packs

	size	each	bulk (12 each)
standard	10" x 12"	11-1310	11-1310-12
neck	7" x 24"	11-1311	11-1311-12
oversize	15" x 24"	11-1312	11-1312-12
half size	5" x 12"	11-1313	11-1313-12
spine small	10" x 18"	11-1314	11-1314-12
spine large	10" x 24"	11-1315	11-1315-12
knee-shoulder	10" x 20"	11-1316	11-1316-12
circular	10" dia	11-1205	11-1205-12

hot water bottle


- Hospital grade rubber with screw stopper
- 2-quart capacity



11-1140 small
11-1140-12 case of 12

Relief Pak® Hot Button® heat compress

- Instant heat therapy without hot water or electricity
- Activate compress by gently flexing disk within the compress
- Instantly heats to approximately 130°F (54°C)
- Generates heat for at least 30 minutes
- Recharge compress by boiling pad
- Pad can be recharged at least 100 times



standard



small



oversize

		each	bulk (12 each)
small	3½" x 5½"	11-1025	11-1025-12
standard	5" x 15"	11-1026	11-1026-12
neck	11" x 17"	11-1027	11-1027-12
oversize	10" x 11"	11-1028	11-1028-12

Relief Pak® moist heat pack heating units

- All heating units are insulated
- Durable seamless stainless steel for easy cleaning
- Heats moist heat packs in a water bath up to 170° F
- On/off switch lights green when unit is on
- Heating light glows orange when unit is warming to desired temperature



- Comes with removable hospital-grade plug set
- Available in 110 volt and 220 volt units

heating units only (110 volt)

stationary	
4 Pack	11-1981
6 Pack	11-1982

mobile	
8 Pack	11-1984
12 Pack	11-1985
24 Pack	11-1986



heating units with packs (110 volt)

4 pack stationary units	13.5" x 10.5" x 18" (L x W x H)
4 std packs	11-1961
3 std, 1 neck pack	11-1961-2
2 std, 2 neck packs	11-1961-3

6 pack stationary units	19.5" x 17.5" x 20" (L x W x H)
3 std, 1 neck, 2 oversize packs	11-1962
6 std packs	11-1962-2

8 pack mobile units	19.5" x 17.5" x 29.5" (L x W x H)
8 std packs	11-1964
4 std, 4 neck packs	11-1964-2
2 std, 2 neck, 2 oversize packs	11-1964-3


12 pack mobile units	22.5" x 18.5" x 33.5" (L x W x H)
12 std packs	11-1965
3 std, 3 neck, 3 oversize packs	11-1965-2
4 std, 4 neck, 2 oversize packs	11-1965-3
4 std, 4 oversize packs	11-1965-4

24 pack mobile units	32.5" x 20.5" x 34.5" (L x W x H)
24 std packs	11-1966
8 oversize, 8 neck packs	11-1966-2
6 std, 6 neck, 6 oversize packs	11-1966-3

cold therapy

ColPac® reusable cold packs

- Reusable cold packs remain pliable and soft, even below freezing
- Place in a chilling unit or freezer
- Maintains therapeutic level for up to 30 minutes
- Use with hand towel between skin and pack
- **Black urethane packs** are more heavy duty than blue vinyl packs



Blue Vinyl ColPac®

each

bulk (12 each)

A	standard (11" x 14")	00-1500	00-1500-12
B	neck (6" x 23")	00-1508	00-1508-12
C	oversize (11" x 21")	00-1512	00-1512-12
D	half size (7" x 11")	00-1506	00-1506-12
E	quarter size (5.5" x 7.5")	00-1504	00-1504-12
F	throat (3" x 11")	00-1502	00-1502-12
G	eye (3" x 8")	00-1510	00-1510-12

Black Urethane ColPac®

each

bulk (12 each)

H	standard (10" x 13.5")	00-1552	00-1552-12
I	neck (6" x 21")	00-1554	00-1554-12
J	oversize (12.5" x 18.5")	00-1556	00-1556-12
K	half size (6.5" x 11")	00-1562	00-1562-12



ColPac® freezer



- 5.0 cubic foot capacity
- Stainless steel shelves, trays, and cup holders are designed to quickly freeze cold packs
- 21"L x 26"W x 33"H; 150 lb
- Comes with 12 standard ColPac® packs

09-0910K freezer w/ packs

ColPac® C-5 chilling unit



- Mobile unit has large storage and chilling capacity
- Requires no additional plumbing
- 27"L x 16"W x 33"H; 65 lb
- Comes with 6 standard and 6 half-size ColPac® packs

00-3102 C-5 unit w/ packs

Relief Pak® Cold n' Hot® Sensaflex® compress

- Cold n' Hot® Sensaflex® compress has a comfortable soft fabric that conforms to the curves of your body, is gentle on the skin, and washable
- Chill in freezer or heat in microwave
- Non-toxic, non-crystallizing gel allows pack to contour to body
- Standard hook and loop adheres to pack to hold in place

		each	bulk
circular	10" dia	11-1297	11-1297-12 (12 ea)
small	3" x 5"	11-1291	11-1291-12 (12 ea)
medium	7" x 12"	11-1292	11-1292-12 (12 ea)
large	10" x 13"	11-1293	11-1293-6 (6 ea)
neck	7" x 22"	11-1294	11-1294-12 (12 ea)
tri-sectional	8" x 16"	11-1298	11-1298-6 (6 ea)



Relief Pak® Cold n' Hot® Donut®

- Combine cold or hot treatment with compression for reduced swelling
- Roll-on pack is easy to apply and stays in place
- Chill in freezer or heat in microwave

	circ	each	bulk (case qty)
Donut® finger	1"	11-1530	11-1530-20 (20 ea)
Donut® small	4"-10"	11-1531	11-1531-10 (10 ea)
Donut® medium	10"-15"	11-1532	11-1532-10 (10 ea)
Donut® large	15"-21"	11-1533	11-1533-10 (10 ea)
Donut® x-large	21"-28"	11-1534	11-1534-10 (10 ea)



Nylatex® wraps



- Use to securely hold cold packs in place

00-1202	2½" x 18"; 3 ea
00-1264	2½" x 24"; 3 ea
00-1204	2½" x 36"; 3 ea
00-1206	2½" x 48"; 3 ea
00-1208	4" x 18"; 3 ea
00-1210	4" x 36"; 3 ea
00-1212	4" x 48"; 3 ea
00-1213	4" x 60"; 3 ea
00-1214	6" x 18"; 3 ea
00-1216	6" x 36"; 3 ea
00-1218	6" x 48"; 3 ea
00-1219	6" x 60"; 3 ea
00-1200	9 pc set (1 each, excluding 00-1213, 00-1264, 00-1219)

Relief Pak® cold pack front loading chilling unit



- 4.6 cubic foot capacity with 2 full-width door shelves and 2 full-width wire shelves
- Adjustable thermostat control
- Defrost drain and manual defrost

11-0501 front loading
11-0506 front loading with 12 Relief Pak® vinyl standard packs (11-1000)

Relief Pak® cold packs

- Effective cold therapy for sprains and strains, muscle soreness, bruises, swelling, and tension headache
- Maintains therapeutic cold for 30 minutes
- Remains pliable, soft even after freezing
- Uniform gel distribution seal prevents gel from dispensing
- **Black urethane packs** made of heavier material than blue vinyl; can chill to lower temperatures



standard



neck

blue vinyl each

standard (11" x 14")	11-1000
neck (6" x 23")	11-1001
oversize (11" x 21")	11-1002
half size (7" x 11")	11-1003
quarter size (5" x 7")	11-1004
throat (3" x 11")	11-1005
circular (10" diameter)	11-1007

bulk (12 each)

11-1000-12
11-1001-12
11-1002-12
11-1003-12
11-1004-12
11-1005-12
11-1007-12

black urethane each

standard (11" x 14")	11-1250
neck (6" x 23")	11-1251
oversize (11" x 21")	11-1252
half size (7" x 11")	11-1253
circular (10" diameter)	11-1257

bulk (12 each)

11-1250-12
11-1251-12
11-1252-12
11-1253-12
11-1257-12

Relief Pak® english-style ice cap

- Bags made of waterproof, rubberized fabric with plastic screw caps



diameter	each	bulk (12 each)
6"	11-1060	11-1060-12
9"	11-1061	11-1061-12
11"	11-1062	11-1062-12



Relief Pak® instant cold packs

- Pack goes from room temperature to 34°F in three seconds
- Squeeze the pack to activate
- Disposable; single use

11-1020	6" x 9"	12 ea
11-1021	4" x 6"	12 ea



Relief Pak® LiTE® Cold n' Hot® pack

- Economical, reusable pack
- Heat in microwave or chill in freezer
- Flexible gel packs conform to body contours

	each	bulk (case quantity)
5" x 6"	11-1054-1	11-1054-48 48 ea
5" x 10"	11-1055-1	11-1055-24 24 ea
8" x 14"	11-1056-1	11-1056-12 12 ea



Relief Pak® Cold n' Hot® elastomer pack

- Pack remains usable if punctured
- Designed for moist heat and soothing cold applications
- Tough, flexible gel, specially selected fabric and Velcro straps make for an exceptionally durable product
- Maintains flexibility at -20°F
- Stretch cover allows for conformity, heat transfer, and comfort
- Heats in microwave, chills in freezer



	each	bulk (12 each)
pad		
12½" x 12½"	11-1631	11-1631-12
wraps w/strap		
5" x 9½"	11-1633	11-1633-12
7¼" x 17¼"	11-1634	11-1634-12
10" x 24"	11-1635	11-1635-12
neck/shoulder w/strap		
11½" x 20"	11-1636	11-1636-12
head w/strap		
3" x 12½"	11-1637	11-1637-12
eye w/strap		
3" x 9"	11-1638	11-1638-12

Game Ready® and Squid cold compression therapy

Game Ready® cold compression

- **Game Ready® GRPro 2.1 control unit** uses active compression that squeezes greater benefits out of cold therapy alone. Circumferential wraps envelope the injury or surgery site and pneumatic compression works to conform the wrap to the contours of the body, allowing for better surface contact and more effective cooling



13-2500 **GRPro 2.1 control unit**
(includes AC adapter & 6 ft connector hose)

- **Game Ready® Med4 Elite™** integrates iceless cold (cryotherapy), controllable heat (thermotherapy), rapid contrast therapy (RCT) and intermittent pneumatic compression therapy (IPC) in one device



13-2501 **Med4 Elite™ control unit**

- **Game Ready ATX® anatomical wraps** feature two internal chambers—one for air, one for ice water. The patented construction enables the integrated, no-mess delivery of adjustable compression and cold therapies.

order Game Ready® control units with wraps of choice

	ATX® wraps	add'l sleeves
upper extremity	upper extremity	upper extremity
hand/wrist	13-2519	13-2561
elbow	13-2523	13-2564
flexed elbow	13-2524	13-2565
shoulder right, medium	13-2525	13-2566
shoulder right, large	13-2526	13-2567
shoulder left, medium	13-2527	13-2568
shoulder left, large	13-2528	13-2569
lower extremity	lower extremity	lower extremity
full leg boot, medium	13-2536	13-2533
full leg boot, large	13-2537	13-2534
ankle, large	13-2511	13-2553
ankle, x-large	13-2512	13-2554
half leg boot	13-2513	13-2555
knee, straight	13-2514	13-2556
knee, articulated	13-2515	13-2557
mid body	mid body	mid body
hip/groin right	13-2516	13-2558
hip/groin left	13-2517	13-2559
back	13-2518	13-2560
CT spine	13-2529	13-2547
cooling vest	cooling vest	cooling vest
cooling vest	13-2601	13-2606
cryo cap	cryo cap	cryo cap
cryo cap with ATX	13-2600	-
traumatic amputee	amputee	amputee
above knee, left	13-2508	13-2604
above knee, right	13-2509	13-2603
below knee, large	13-2510	13-2552
utility	13-2507	13-2605

Squid active cold compression

- Use for temporary relief of aches and pains
- Portable, air pressure pump is small and lightweight (under 1 lb / 0.45 kg)
- 4 preprogrammed levels of compression from low (30 mmHg) to high (85 mmHg) to
- Sets include pump, wrap, and gel pack



shown with Squid pump, ankle wrap and gel pack



- **Squid clinical units** have a rechargeable battery that allows for **50 uses** on one charge
- **Squid Go home units** have a rechargeable battery that allows for **15 uses** on one charge



	clinical set	home set	wrap/gel pack only
Upper Extremity Wrap			
wrist	13-2465	13-2480	13-2468
elbow	13-2425	13-2473	13-2428
shoulder, left, regular*	13-2450	13-2474	13-2452
shoulder, left, large*	13-2457	13-2475	13-2459
shoulder, right, regular*	13-2453	13-2476	13-2455
shoulder, right, large*	13-2460	13-2477	13-2462
*only offers intermittent compression			
Lower Extremity Wrap			
leg, regular	13-2435	13-2478	13-2438
leg, large	13-2441	13-2479	13-2442
ankle, regular	13-2410	13-2470	13-2413
ankle, large	13-2414	13-2471	13-2417
knee	13-2430		13-2433
leg circumferential, regular	13-2439		13-2447
leg circumferential, large	13-2445		13-2448
Mid Body Wrap			
back	13-2420	13-2472	13-2422

Ossur® Cold Rush® cold compression

- Insulated cooler keeps ice cold for up to 6 hours before unit needs to be refilled with ice / water; Use with anatomically designed pads sold with each unit in sets and separately

Cold Rush® standard unit

- Whisper-quiet (30 decibel) brushless motor provides a strong adjustable flow rate; ideal for patients to use while resting
- Locking lid prevents spills
- Lightweight design; 4.6 lbs

Units Only

Cold Rush® unit only 13-2700

Units with standard pad

universal 13-2702
rectangular, small 13-2708

Units with Upper Extremity pad

shoulder, small 13-2705
shoulder, large 13-2706

Units with Lower Extremity pad

knee 13-2701
hip, left 13-2703
hip, right 13-2704
foot and ankle 13-2707

Units with Mid Body pad

lumbar 13-2709

Cold Rush® compact unit

- Highly portable unit weighs 3.5 lbs
- Quiet motor has single flow rate

Compact

13-2710

pads only

13-2712 13-2723
13-2718 13-2724

13-2715 13-2721
13-2716 13-2722

13-2711 13-2720
13-2713 13-2728
13-2714 13-2727
13-2717 13-2726

13-2719 13-2725



Cold Rush® standard unit w/ universal pad



Cold Rush® compact unit w/ universal pad

Fluidotherapy® dry heat



single

double

- Increase circulation and ROM with the soothing dry heat
- Controls include: temperature (max 125°), airflow (velocity of particulate medium), pulse mode (pressure bursts)
- Includes 20 lb of Cellex™ medium



FLU110D FT110 single
FLU115D FT115 double
MED0001 Cellex™ 10 lbs

Richmar® TheraTouch® DX2 shortwave diathermy

- Shortwave diathermy device with scientific dosage controls
- Exclusive Scientific Dosage Control (SDC) provides an increased thermal effect with decreased risk
- Allows for 6 dosage levels:
dose I - subthermal;
dose II - mild heat;
dose III - low heat;
dose III SDC M - moderate heat;
dose III SDC H - high heat;
dose IV - vigorous heat
- Auto-tuning technology adjusts to ensure maximum energy transmission without resorting to a higher dosage level
- Use to treat muscle spasms, joint contractures and general pains
- Easy-to-use touchscreen interface
- Save up to 99 user-defined settings
- Dual treatment indicators on unit and inductive arm



13-3270 TheraTouch® diathermy

Mettler Auto*Therm® 390 shortwave diathermy

- Pulsed and continuous shortwave diathermy
- Truly portable allowing to move it easily between treatment rooms
- For a more versatile modality, the 391 package allows the clinician to use both capacitive pads and coil induction
- ME390 package:** includes two 4.7" x 7.1" soft rubber plate applicators with flexible cables, 6 felt spacers, 2 cloth covers and two 48" elastic velcro strap (no cart or arm)
- ME391 package:** ME390 PLUS instrument cart, four-joint folding arm, one 14 cm diameter coil applicator and one special field coil applicator



Auto*Therm® 390 packages

13-3062 ME390 package
13-3063 ME391 package

Chattanooga® Intellect® SWD100 shortwave diathermy

- Effective deep tissue heating therapy
 - Full-color touch-screen interface, onboard clinical resources, customizable treatment parameters and more
- Ideal for:**
- Pain relief
 - Contractures
 - Bursitis
 - Decreasing joint stiffness
 - Increased blood flow
 - Reduction of muscle spasm
 - Chronic inflammatory conditions



shortwave diathermy

00-1600 Intellect® SWD100

accessories

01-4790 electrode arm (left)
01-4791 electrode arm (right)
22-0000-2 monode
02-0969-553 diplode

capacitive electrodes

02-0453-216 80 mm
02-0453-217 120 mm
02-0453-218 165 mm

soft rubber electrodes

02-0453-266 180 mm x 120 mm
02-0453-267 250 mm x 145 mm

ultrasound

Sonicator® 718 portable ultrasound

- Small portable size
- Handheld operation
- Rechargeable battery and charger for portable power
- Full clinical power, 2 W/cm² continuous, 3 W/cm² pulsed portable ultrasound
- 5.5 cm² soundhead
- Package includes rechargeable battery (13-3408), battery charger (13-3410), power cable (13-3411) and battery case (13-3409)

13-3400 Sonicator 718® package



US Pro 2000™ portable ultrasound



- 3 intensity settings and 3 timer modes
- Perfect for clinical or home care use
- 1 MHz working frequency
- 8"L x 2"W x 2¾"D, 2 lb
- Includes: 2.8 oz conductive gel, AC adapter and carry case

13-1471 US Pro 2000™ ultrasound

US 1000™ portable ultrasound



- 3 intensity settings and 30 minute timer setting with auto off
- Perfect for clinical or home care use
- Includes: ultrasound gel, AC adapter and carry case

13-1470 US 1000™ ultrasound

Intellect® ultrasound gel



Intellect® ultrasound gel

00-4266-1 5 liter (1.3 gal)
00-4266-4 5 liter (4/cs)

Conductor™ gel

- Highly conductive gel
- 00-4248-1 8.5 oz
00-4248-12 8.5 oz (12/cs)
00-4248-24 8.5 oz (24/cs)
00-4238-1 5 liter (1.3 gal)
00-4238-4 5 liter (4/cs)

Aquasonic® 100 ultrasound gel



- Non-staining, hypoallergenic, and non-irritating

50-5810-1 0.25 liter
50-5810-12 0.25 liter (12/box)
50-5810-72 0.25 liter (72/box)
50-5812-1 5 liter
50-5812-4 5 liter (4/cs)

Polysonic® ultrasound gel



without aloe vera


50-6001-1 8.5 oz tube
50-6001-12 8.5 oz tube 12/box
50-6001-72 8.5 oz tube 72/box
50-6002-1 1-gallon dispenser
50-6002-4 1 gallon dispenser 4/cs

with aloe vera

50-6004-1 8.5 oz tube
50-6004-12 8.5 oz tube 12/box
50-6004-72 8.5 oz tube 72/box
50-6005-1 1 gallon
50-6005-4 1 gallon dispenser 4/cs

Ionto+® hi-performance electrodes and ibox2 unit

- Iontophoresis unit delivers both single and dual phase treatments
- Electrodes feature easy-fill window
- Large adhesive area prevents pad movement
- Silver chloride conductor for higher conduction
- Buffering agents for pH stabilization



iontophoresis unit


13-5201 ibox2 unit

12 treatment electrodes per box

- A 13-5250 wings
B 13-5251 small square
C 13-5252 medium square
D 13-5253 large square

IontoPatch® transdermal drug delivery system

- Self-contained wireless technology eliminates need for external batteries or wires
- Allows patient to return to daily activities while receiving time-released iontophoresis
- Made of flexible, breathable material
- Technician applies medication directly to the electrode
- Shuts off automatically when prescribed dose has been administered
- 6 treatment kits per box (patch, saline vial, alcohol prep)



	IontoPatch®	active area	dosage	fill volume	wear time
13-5220	80	15.5 cm ²	80 mA-min	1.2 - 1.6 cc	14 hours
13-5221	SP	8.0 cm ²	40 mA-min	1.2 - 1.6 cc	14 hours
13-5223	Stat	12.0 cm ²	80 mA-min	1.2 - 1.4 cc	4 hours
13-5224	Extra Strength™	12.0cm ²	120 mA-min	1.2 - 1.6 cc	8 hours

FabStim® TENS 2-channel

- 3 modes (Burst, Modulation, Constant) and timer
- Adjustable pulse rate and width
- Includes carry case, stimulator, lead wires, 4 self-stick electrodes and 9V battery



13-1323 analog TENS

FabStim® digital TENS 2-channel


- 5 modes: Burst, Modulation, Constant, MRW, SD1 and SD2
- Adjustable pulse rate and width
- Timer (5-60 mins or continuous) and compliance meter
- Includes carry case, stimulator, lead wires, 4 self stick electrodes and 9V battery



13-1308 digital TENS

FabStim® electrodes

- Self-adhesive (Tyco® gel) electrodes good for 8-12 uses
- Conveniently packed in re-sealable bags of 4 for individual patients; 10 bags per case (40 electrodes/case)
- Compatible with electrotherapy devices including TENS, EMS, IF, High Volt and Russian



	bag of 4	case (10 bags)
3" round	13-1288-1	13-1288-10 40/case
2" round	13-1287-1	13-1287-10 40/case
2" square	13-1291-1	13-1291-10 40/case
2" x 3.5" rectangle	13-1293-1	13-1293-10 40/case

Thermasonic® gel warmers

- Rapidly heat ultrasound gel and consistently monitor temperature
- Powerful heating element uniformly heats 8oz (250 ml) bottles of gel
- For wall mounting or countertop use



- E 50-5863 3-unit bottle warmer LCD
F 50-5862 3-unit bottle warmer LED
G 50-5861 single bottle warmer

3 bottle electric warmer

- Warm up to three 8 oz (250 ml) bottles (included)
- Specify blue (B), white (W) or gray (GRY)



50-6068 8 oz warmer (B/W/GRY*)

*Note: if color is not specified, blue is shipped

VitalStim® PLUS 4-channel stim

- Rehab or strengthen swallowing muscles
- Preset or customize programs
- Screen Mirroring helps clinicians guide patients through personalized, targeted swallowing exercises
- sEMG Biofeedback System projects progress onto computer or tablet screen



unit

68-0066 VitalStim® Plus electrotherapy system

accessories

- 68-0067 VitalStim® Plus snap lead wires
68-0068 VitalStim® Plus reference Emg Lead wire
68-0069 VitalStim® Plus hand switch
68-0070 VitalStim® Plus rubber sleeve
68-0071 VitalStim® Plus Stylus
68-0072 VitalStim® Plus Stand
68-0073 VitalStim® Plus Battery Door
68-0074 VitalStim® Plus reference emg electrode

Richmar stimulation and ultrasound units

SoundCare™ Plus ultrasound

- Ultrasound device operates in 1MHz and 3MHz frequencies
- Includes 2 sound heads, 1 cm² and 5 cm²
- Easy push button switches between applicator heads
- 20 preset programs

13-3380 SoundCare™ Plus



Quattro™ 2.5 4-channel stim

- Waveforms include: TENS, EMS, 2-pole interferential, 4-pole interferential and Russian
- Includes: 8 sets of electrodes, 8 sets of rubber electrodes with sponges, 8 straps, 4 sets of lead wires and user manual

13-3370 Quattro™ 2.5 4-channel stim



ComboCare™ stim / ultrasound combo

- 2 channel stim and ultrasound device
- 5 treatment waveforms: IF 4-pole, IF 2-pole, Russian, EMS and TENS
- Preset programs: 50 wave, 20 ultrasound

13-3390 ComboCare™ stim and ultrasound



TheraSound EVO ultrasound

- Dual frequency (1 MHz and 3 MHz) ultrasound therapy device
- Comes with the 2-in-1 Therapy Hammer applicator, providing both 5 cm² and 2 cm² sound heads

13-3385 TheraSound EVO ultrasound



Winner EVO system

- Two or four channel electrostimulation (ST2 & ST4) and combination two or four channel electrostimulation and ultrasound (CM2 & CM4) clinical devices
- Light Therapy and Cold Laser Therapy can be added to combo units with addition of Laser Module

stim / ultrasound combo units

- 69-0101 EVO CM2, 2-channel
- 69-0102 EVO CM4, 4-channel
- 69-0103 EVO CM4 combo unit with cart

stim units

- 69-0222 EVO ST2, 2 channel
- 69-0223 EVO ST4, 4-channel

accessories

- 69-0226 AutoSound hands-free ultrasound applicator
- 69-0227 light cluster probe (500 mW - 870 nm LED / 633 nm SLD)
- 13-4550 laser module
- 13-4551 LaserPrism tethered laser (200 mW / 785 nm)
- 13-4552 LaserPrism cordless laser (450 mW / 808nm)



stim unit only



stim / ultrasound system with cart

TheraTouch™ professional 4-channel system

- Four channel stim (EX4) and combination stim and ultrasound (CX4) clinical devices
- Menu-driven color screen interface intuitively groups and displays comprehensive clinical protocols as well as user-defined manual operation



CX4 unit



EX4 unit

CX4 4 channel stim / ultrasound combo

- 13-3391 CX4 unit
- 13-3392 CX4 unit w/ cart

EX4 4-channel stim unit

- 13-3371 EX4 unit
- 13-3372 EX4 unit, w/ cart

professional series cart

- 13-3399 mobile cart

Hivamat 200 oscillation



- Provides deep oscillation therapy for pain relief, muscle relaxation, and increasing circulation
- **Portable model** uses predefined programs from 3 practice-specific treatment protocol cards: athletic, rehab, and chiropractic

Hivamat 200 portable

- 70-0371 athletic
- 70-0372 rehab
- 70-0373 chiropractic

Chattanooga stimulation and ultrasound units

TranSport® 2-channel

- Lightweight, portable stim / ultrasound units can be used in the clinic on a cart or taken on the go
- 2-channel stim
- 4 waveforms: Interferential, Premodulated, High Volt, and Russian
- Dual frequency ultrasound 1 / 3.3 MHz; 5cm² head standard; 1, 2, and 10 cm² heads available



combo stim / ultrasound units

- 00-2738K 2-channel system
- 00-2738KC w/ mobile cart
- 00-2738KB w/ bag and battery pack

stim units

- 00-2783K 2-channel system
- 00-2894K w/ mobile cart
- 00-2783B w/ bag and battery pack

ultrasound units

- 00-2782K ultrasound
- 00-2893K ultrasound with cart

Intellect® LEGEND XT™ combo

- 6 waveforms: Interferential, Premodulated, High Volt, Russian, Symmetrical Biphasic and Microcurrent
- Dual frequency ultrasound 1 / 3.3 MHz; 5cm² head standard; 1, 2, and 10 cm² heads available



combo stim / ultrasound units

- 00-2760 2-channel system
- 00-2791 2-channel system w/cart
- 00-2788 4-channel system
- 00-2795 4-channel system w/ cart

stim units

- 00-2763 2-channel system
- 00-2793 2-channel system w/cart
- 00-2786 4-channel system
- 00-2797 4-channel system w/ cart

Vectra® Neo clinical system

- 5 modules: 2 and 4 channel stim, laser, Ultrasound, EMG Biofeedback, EMG-Triggered Stim
- Dual frequency ultrasound 1 / 3.3 MHz; 5cm² head standard; 1, 2, and 10 cm² heads available
- Touchscreen offers 100 protocol options



Vectra® Neo unit

00-6000 Unit Only

modules for Vectra® Neo

- 07-0000 Stim Channels 1 & 2
- 07-0004 Stim Channels 1 & 2 + Emg
- 07-0003 Stim Channels 3 & 4
- 07-0002 Ultrasound
- 07-0005 Laser

accessories

- 07-0001 Cart
- 07-0008 Patient Remote

Vectra® Genisys®

- Modular therapy system offers ultrasound, stim, ultrasound / stim combo, laser, sEMG, and sEMG + stim
- Stim available in 2 or 4 channels, with or without sEMG
- Dual frequency ultrasound 1 / 3.3 MHz; 5cm² head standard; 1, 2, and 10 cm² heads available



combo stim / ultrasound with sEMG

- 00-2761 2-channel system
- 00-2792 2-channel system w/cart
- 00-2789K 4-channel system
- 00-2796 4-channel system w/ cart

combo stim / ultrasound w/o sEMG

- 00-2872 2-channel system
- 00-2876 2-channel system w/cart
- 00-2874 4-channel system
- 00-2878 4-channel system w/ cart

stim with sEMG

- 00-2764K 2-channel system
- 00-2787K 4-channel system
- 00-2798 4-channel system w/ cart

stim without sEMG

- 00-2871 2-channel system
- 00-2873 4-channel system

Cube Laser

- Laser / biostimulation pain relief system
- Stimulates tissue growth and repair at cellular level

- 13-4600 CUBE 3, 3 wavelengths
- 13-4601 CUBE 4, 4 wavelengths
- 13-4602 CUBE 4 Plus, w/ optical handpiece



Intellect® mobile 2 RPW radial pressure wave

- Powerful, yet portable, radial pressure wave device
- Quiet, built-in compressor
- Pressure, frequency and pulses can be monitored on LED display
- Pulses can be initiated via a trigger on the handpiece



- 13-4700 mobile 2 RPW

Focus Shockwave

- Electromagnetically generates an extracorporeal shockwave in a short amount of time
- Allows for deep, focused shockwave treatment at a depth of 1 inch
- LCD touchscreen display



- 13-4750 focus shockwave

Mettler stimulation and ultrasound units

Sys*Stim® 208 and 208A low-volt stims

- Knobs control stimulation intensity and frequency, digital timer and remote patient cut-off switch
- 208 single-channel stimulator** has Surge, Pulse and Tetanize modes
- 208A dual-channel stimulator** has Surge, Pulse, Tetanize and Reciprocation modes



13-3015 Sys*Stim® 208, 1-channel
13-3016 Sys*Stim® 208A, 2-channel

Sys*Stim® 206 stim

- Single channel stimulator with 4 waveforms
- Narrow pulse is suited for TENS while wide and AC pulses can be used for muscle stimulation
- Direct current can be used to stimulate denervated muscles or to perform iontophoresis



13-3021 Sys*Stim® 206, 1-channel

Sonicator® Plus 992 and 994 combo

- Microprocessor-controlled combo
- 6 waveforms: Interferential, Premodulated, Medium (Russian), Biphasic, Hi-Volt and Micro-current
- Includes 5 cm² applicator

Sonicator® Plus 992

- 2-channel stimulation
- 3 protocols can be set-up using different timers

Sonicator® Plus 994

- 4-channel stimulation
- 5 protocols can be set-up using different timers



13-3001 Sonicator® Plus 992, 2-channel
13-3000 Sonicator® Plus 994, 4-channel

Each Sonicator® Plus comes standard with a 5 cm² dual frequency 1 and 3 MHz applicator

Sonicator® 740, 740x ultrasound

- 2-frequency (1, 3.3 MHz) delivered in continuous or 10, 20 or 50% pulsed modes
- Applicator is contoured for easy application, has an LED in the handle that lights up when adequate coupling achieved
- 740** has a 5 cm² applicator
- 740x** has 3 applicators (1, 5, 10 cm²)



13-3300 Sonicator® 740 w/ 1 head
13-3301 Sonicator® 740x w/ 3 heads
13-3320 battery (optional)
13-3034 universal applicator cable
13-3064 carry bag

Sonicator® Plus 921 and 941 combo

- Intensity control knobs for stimulation and ultrasound precisely set treatment intensity
- Treatment parameters can be changed during treatment allowing patients to feel **real time** change ensuring they feel comfortable and safe
- Vivid 7" color touch screen allows quick and easy setup
- Store up to 200 treatment setups
- 8 waveforms: Interferential (IF-4), Premodulated (IF-2), EMS, Russian, TENS, Hi-Volt, Microcurrent, and DC Galvanic
- Use with small (2cm², not included) and large (11 cm², included) ultrasound applicators; each has 1MHz and 3MHz frequencies



13-4800 Sonicator® Plus 921, 2-channel
13-4801 Sonicator® Plus 941, 4-channel

Sonicator® Plus 920 combo

- Touch screen combination offers Interferential (4-pole), Premodulated (2-pole), Medium frequency, Biphasic, Hi-Volt, TENS, Micro-current and DC waveforms
- Allows clinician to treat 3 patients at the same time
- 61 pre-set protocols with additional space for your own protocols
- 11 languages and internal safety features
- Additional 0.9 cm² applicator available



13-3450 Sonicator® Plus 920, 2-channel

upholstery
color



sand



camel



slate



wedgewood



black

hand therapy tables and stools

mobile revolving stool

- Height adjustable (18"-24") "screw-type" revolving stool with round footrest
- 2" hooded carpet casters
- 14" diameter, 4" thick padded seat
- 250 lb capacity



16-1102 stool, black
16-1100 specify color from above chart

treatment (P.T.) stool

- 3" padded seat
- 18" fixed height
- Available with or without 2" hooded carpet casters
- 250 lb capacity



mobile
16-1602 mobile, gray
16-1600 specify color from chart
stationary (not shown)
16-1601 specify color from chart

pneumatic stools

- Therapy stool with round, 3" upholstered swivel seat allows easy access to the patient
- Gas spring for easy raising and lowering of stool (18"-22")



without back



with back

without back
07-7063 black
07-7060 blue
07-7061 burgundy
07-7062 gray

with back
07-7067 black
07-7064 blue
07-7065 burgundy
07-7066 gray

hand therapy table

- Hand therapy table for therapist and patient
- Patient can use while sitting, standing, ambulating, or in a wheelchair. 13" height range
- Lightweight, easily maneuverable for use anywhere



no drawer
15-4305 37" x 19" x 29"-42"
with drawer
15-4306 37" x 19" x 29"-42"
with front cut out
15-4307 37" x 24" x 29"-42"

4-cutout hydraulic table



- Hydraulic foot activated mechanism raises table height 28½" to 34"
- Laminate top swivels 360° and locks in place
- Table accommodate up to 4 wheelchair patients

66-0019 48" x 48" x 28.5"-34"

horseshoe therapy table



- Multi-purpose wheelchair accessible table is ideal for work, play or any other activity
- U-shape top accommodates up to 4 wheelchair patients
- Top is 18" deep for easy access
- 5 adjustable metal legs

15-4300 72" x 48" x 26"-34"

mobile cabinet for splinting and supplies

- Grab-bar handle for mobility
- Door cabinet with lock
- High impact polystyrene drawer
- 3 vertical cubby sections for splinting material behind hinged door
- Durable dowel and glue construction
- 4" twin-wheel locking casters



shown with splinting pan
(see page 45)

15-4301 42"x 24" x 31" cabinet

Cubex™ therapy system on wheels

- Reversible activity board in top has pegboard pad on one side and smooth easel on the other side
- Foldaway bolt board on right side and pegboard on left side
- Includes 100 plastic stackable pegs, 24 assorted bolts, 2 screwdrivers, 2 wrenches and Allen wrench
- 2 door cabinet with lock and 4 tubs
- Table mounted finger ladder offers progressive mobility of arm at shoulder, elbow and wrist joints



15-4302 Cubex™ therapy system
15-4304 Cubex™ therapy system w/ finger ladder
15-4303 finger ladder only

fixed height treatment tables

- 30" height
- 2" upholstered top
- Steel frame
- 4 corner stability
- 10 year limited warranty
- 500 lb weight limit



L x W	basic	H-brace	shelf
72" x 27"	15-5095	15-5085	15-5090
72" x 30"	15-5096	15-5086	15-5091
78" x 27"	15-5097	15-5087	15-5092
78" x 30"	15-5098	15-5088	15-5093

classroom/trainer tables

- 30" height
- 2" removable mat top (optional)
- Plastic laminated surface
- 10 year limited warranty
- 400 lb weight limit



L x W	table w/o mat	removable mat
72" x 27"	15-5100	15-5105
72" x 30"	15-5101	15-5106
78" x 27"	15-5102	15-5107
78" x 30"	15-5103	15-5108

motorized hi-lo treatment tables

- Electric hi-lo / back / leg 23.25" to 34.75"
- 2" upholstered top
- Generous storage area
- Steel frame
- 10 year limited warranty
- 500 lb weight limit



L x W	1 section	2 section	3 section
80" x 28"	15-5080	15-5082	15-5081

bariatric motorized hi-lo treatment tables

- Electric hi-lo / back / leg 17.5" to 36.75"
- 2" removable mat top
- Locking casters
- Hand control
- Steel frame
- 10 year limited warranty
- 750 lb weight limit

L x W	1 section
78" x 32"	15-5110
78" x 36"	15-5111
78" x 42"	15-5112
78" x 48"	15-5113
L x W	2 section
78" x 32"	15-5115
78" x 36"	15-5116
78" x 42"	15-5117
78" x 48"	15-5118
L x W	3 section
78" x 32"	15-5120
78" x 36"	15-5121
78" x 42"	15-5122
78" x 48"	15-5123




treatment table accessories

15-5070	caster (4 ea)		
15-5071	paper holder	15-5073	enclosed linen shelf
15-5072	siderails	15-5074	facial slot

motorized parallel bars

- Electric adjustable height: 25.5" to 39.75"; individual or dual operation
- Manual adjustable width: **Standard** 15.5" to 24.5"; **Bariatric** 21.5" to 30.5"



- Steel base with 1ft tapered end and/or 2 ft platform extension (optional)
- Mat for active surface
- Aluminum handrail bar diameter 1.6"; Child's low bar diameter 1.3"
- 10 year limited warranty
- Weight limit: **Standard** 600 lb; **Bariatric** 1,000 lb; **child's low bar** 250 lb

	standard				
	Length	parallel bars	platform width	active width	child's low bar
	6 ft	15-5015	48"	31.5"	15-5000
	8 ft	15-5016	48"	31.5"	15-5001
	10 ft	15-5165	48"	31.5"	15-5002
	12 ft	15-5017	48"	31.5"	15-5003
	16 ft	15-5166	48"	31.5"	15-5004
	18 ft	15-5018	48"	31.5"	15-5005
	tapered extensions				
	1 ft	15-5012			
	3 ft	15-5010			

	bariatric			
	Length	parallel bars	platform width	active width
	6 ft	15-5160	54"	37.5"
	8 ft	15-5161	54"	37.5"
	10 ft	15-5162	54"	37.5"
	12 ft	15-5019	54"	37.5"
	16 ft	15-5163	54"	37.5"
	tapered extensions			
	1 ft	15-5013		
	3 ft	15-5011		

motorized hi-lo mat tables

600 lb weight limit

- Electric adjustable height
 - 2" upholstered top
 - Foot switch
 - Steel frame with rounded corners
 - 10 year limited warranty
- | L x W | height range | |
|-------------|------------------|---------|
| 7 ft x 4 ft | 18.5" to 29.5" | 15-5020 |
| 7 ft x 5 ft | 18.5" to 29.5" | 15-5021 |
| 8 ft x 6 ft | 18.5" to 29.5" | 15-5022 |
| 7 ft x 4 ft | 23.25" to 34.75" | 15-5025 |
| 7 ft x 5 ft | 23.25" to 34.75" | 15-5026 |
| 8 ft x 6 ft | 23.25" to 34.75" | 15-5027 |



motorized hi-lo deluxe mat tables

750 lb weight limit

- Electric adjustable height: 17" to 39"
 - 2" upholstered top
 - Hand switch
 - Locking casters
 - Steel frame w/ rounded corners
 - 10 year limited warranty
- | L x W | |
|-------------|---------|
| 7 ft x 4 ft | 15-5035 |
| 7 ft x 5 ft | 15-5036 |
| 7 ft x 6 ft | 15-5037 |



bariatric motorized hi-lo mat tables

1,000 lb weight limit

- Electric height / back / leg section 18" to 38"
- 2" upholstered (removable) top
- Hand switch
- Steel frame with rounded corners
- Locking and/or leveling casters
- 10 year limited warranty



L x W	1 section	2 section	3 section
7 ft x 5 ft	15-5040	15-5045	15-5050
7 ft x 6 ft	15-5041	15-5046	15-5051

fixed height mat tables

1,000 lb weight limit

- 18" fixed height
- 2" upholstered top
- Steel frame with rounded corners
- 4 corner stability
- 10 year limited warranty



L x W	
7 ft x 4 ft	15-5030
7 ft x 5 ft	15-5031
7 ft x 6 ft	15-5032

bariatric mat table accessories

- 15-5056 sit to stand bars (2 ea)
- 15-5057 guard rail
- 15-5058 T-bar
- 15-5059 bracket

resistance band poles

- 15-5060 48" pole
- 15-5061 72" pole



motorized tilt tables



	L x W	height	weight limit	
standard	78" x 28"	33"	350 lb	15-5152
wide footboard	78" x 28"	33"	350 lb	15-5150
heavy duty	81.5" x 31"	34.5"	500 lb	15-5153
bariatric	81.5" x 31"	34.5"	750 lb	15-5154

- Electric tilt 0 to 90°
- Angle indicator
- 2" upholstered top
- 5" locking casters
- Hand control
- Steel frame
- Safety belts
- 10 year limited warranty

motorized hi-lo tilt table

- Electric adjustable height / tilt: 24" to 35.5" height, 0 to 86° tilt
- Angle indicator
- 2" upholstered top
- 4" locking swivel casters
- Hand control
- Fold down foot platform
- Safety belts
- 10 year limited warranty
- 350 lb weight limit



L x W	
81.5" x 30"	15-5151


motorized Trendelenberg / tilt table

- Electric tilt / Trendelenberg
- Height: 33"
- 2" upholstered top
- 3" tubular steel
- 5" casters
- Safety belts
- Removable foot board
- 10 year limited warranty
- 350 lb weight limit



range	L x W	
0 to -12° Trendelenberg	78" x 28"	15-5155
0 to 86° tilt		

choose vinyl
top color



motorized hi-lo treatment tables

- Electric hi-lo / back / leg 18" to 37"
- 1.5" foam top
- Footswitch control
- Power assisted caster system
- Dimensions: 76"L x 27"W
- 400 lb weight limit



1 section	
standard	15-1700
2 section	
w/ face opening	15-1702
w/ no face opening	15-1815
w/ 3 pc head	15-1705
3 section	
w/ fixed center	15-1709
w/ fixed center, 3 pc head	15-1716
w/ elevated center	15-1706
w/ elevated center, 3 pc head	15-1715
4 section	
standard	15-1712

motorized single pedestal hi-lo treatment tables

- Electric hi-lo / back / leg 24" to 35.5"
- 1.5" foam top
- Footswitch control
- One-step caster system
- 400 lb weight limit



1 section	L x W	
standard	76" x 27"	15-1736
2 section		
w/ face opening	76" x 27"	15-1737
w/ prenatal top, 3" soft foam	76" x 30"	15-1738
3 section		
fixed center	76" x 27"	15-1740
w/ elevated center	76" x 27"	15-1739
w/ elevated center, 3 pc head	76" x 27"	15-1743
4 section		
traction table	76" x 27"	15-1741
traction table w/ 3 pc head	76" x 27"	15-1742
6 section		
w/ 4 locking swivel casters	77" x 25"	15-1744

motorized bariatric hi-lo treatment tables

- Electric hi-lo / back / leg 18" to 37"
- 1.75" foam top
- Footswitch control
- Power assisted caster system
- 500 lb weight limit



1 section	L x W	
w/ 1.5" foam top	76" x 34"	15-1816
w/ 1.5" foam top	76" x 40"	15-1817
2 section		
w/ face opening	76" x 34"	15-1703
w/ no face opening	76" x 40"	15-1704
3 section		
w/ fixed center	76" x 34"	15-1707
w/ fixed center	76" x 40"	15-1708
3 section, 800lb weight limit, hi-lo 22" to 38", casters		
w/ fixed center	76" x 36"	15-1711


motorized SX hi-lo treatment tables

- Electric hi-lo / back / leg 18" to 37"
- 1.5" foam top
- Footswitch control
- 4 locking swivel casters
- 400 lb weight limit



1 section	L x W	
standard, no side rails	72" x 27"	15-1746
changing table w/ side rails	60" x 25"	15-1747
changing table w/ side rails	72" x 25"	15-1753
2 section		
w/ face opening	72" x 27"	15-1748
3 section		
fixed center	72" x 27"	15-1750
w/ elevated ctr	72" x 27"	15-1749
w/ elevated ctr, 3 pc head	72" x 27"	15-1814
4 section		
standard	76" x 27"	15-1751

choose vinyl
top color



Armedica
Manufacturing Corporation

fixed height steel mat table

- 20" fixed height
- 2" firm density foam top
- Steel frame

L x W	weight limit	1 section
7 ft x 4 ft	1,000 lb	15-1762
7 ft x 5 ft	1,000 lb	15-1764
8 ft x 6 ft	1,000 lb	15-1765
2 section		
7 ft x 4 ft	900 lb	15-1763



manual hi-lo steel mat table

- Manual hand crank adjustable height 20" to 30"
- 2" firm density foam top
- Steel frame
- 900 lb weight limit

L x W	1 section
7 ft x 4 ft	15-1758
2 section	
7 ft x 4 ft	15-1759



motorized hi-lo steel mat table


- Electric adjustable height 20" to 30"
- 2" firm density foam top
- Steel frame
- 900 lb weight limit

L x W	1 section
7 ft x 4 ft	15-1754
7 ft x 5 ft	15-1756
8 ft x 6 ft	15-1757
2 section	
7 ft x 4 ft	15-1755



motorized platform mounted parallel bars

- **Motorized platform model** electronically adjusts bar height from 26" to 48" and shows exact height on digital display
- Width manually adjusts 15" to 29";
- Platform measures 10 ft L x 43" W



Length

7 ft	15-1780
10 ft	15-1779

manual platform and floor mounted parallel bars



- **Platform and floor mounted models:** height manually adjusts 26" to 41"; width manually adjusts 19" to 26"

platform mounted

15-1781	7 ft length
15-1782	10 ft length

floor mounted

15-1783	7 ft length
15-1784	10 ft length

CANDo®
treatment furniture

upholstery
color options

black

slate
blue

royal
blue

warm
gray

desert
tan

finish options

natural

dark
cherry

CanDo® Staircases

- Steel fasteners and metal braces provide support and stability
- Molded, solid hardwood handrails offer patient comfort
- Safety treads on each step
- Two-sided staircase patients can ascend an descend without turning
- One-sided staircase saves space
- 500 lb weight capacity



15-4230	two-sided staircase	54"L x 34"W x 48"H
15-4238	one-sided staircase	56"L x 34"W x 60"H

CanDo® wall mounted mat platform tables

- Space saving, folding mat platform table
- 3 front hardwood legs lock in place when folded
- 2" firm foam padding
- 500 lb weight capacity
- Choose upholstery color and finish from table above
- 21" fixed height

15-2062	7 ft L x 4 ft W
15-2063	7 ft L x 5 ft W




MEDCOR PROFESSIONALS • 152 US Route One, Scarborough, ME 04074 • phone: 207-222-2828

Contact Info: Jacob Currier phone: 207-949-2215 email: jcurrier@medcorpro.com

Jason Finoia phone: 207-844-4203 email: jfinoia@medcorpro.com • Kevin Melod phone: 207-838-1172 email: kmelody@medcorpro.com

CanDo® treatment table with shelf

- Features hardwood legs, laminate shelf
- Scratch resistant topcoat finish
- Seamless round corner top
- 2" firm foam padding
- 400 lb weight capacity
- Available with adjustable back or fixed flat top
- Choose upholstery color and finish from table above



flat top



adjustable back

flat top


15-4249 72"L x 27"W x 31"H
15-4250 72"L x 30"W x 31"H

adjustable back

15-4231 72"L x 27"W x 31"H
15-4232 72"L x 30"W x 31"H
15-4233 78"L x 30"W x 31"H

CanDo® H-Frame treatment table

- Hardwood legs
- Seamless round corner top
- 2" firm foam padding
- 400 lb weight capacity
- Available with adjustable back or fixed flat top
- Choose upholstery color and finish from table above



flat top



adjustable back

flat top

15-4247 72"L x 27"W x 31"H
15-4248 72"L x 30"W x 31"H

adjustable back

15-4245 72"L x 27"W x 31"H
15-4246 72"L x 30"W x 31"H

CanDo® upholstered mat table

- Features solid hardwood legs and plywood frame
- Heavy duty, upholstered mat platform
- 6 leg design with center support adds stability
- 2" ultra-firm, foam padding
- Gray edge bumper guard helps prevent tears
- 700 lb weight capacity
- Choose upholstery color from table above



15-4235 84"L x 48"W x 18"H
15-4236 84"L x 60"W x 18"H
15-4237 96"L x 72"W x 18"H

CanDo® hi-lo upholstered mat table

- Electric foot-controlled height adjusts from 22" to 34"
- Ultra-strong uniframe base; powder-coated steel frame
- Custom bumper strip protects top's edges from tearing
- 2" foam padding
- 600 lb weight capacity
- Choose upholstery color from table above




15-4251 84"L x 48"W
15-4252 84"L x 60"W
15-4253 96"L x 72"W

Skillbuilders® rolls, half rolls and wedges

- Skillbuilders® premium rolls and wedges have integral waterproof polyurethane coating that is washable
- Use for exercise therapy and general positioning



rolls



wedges

Skillbuilders® Wedges

	L x W x H	incline
30-1010	22" x 20" x 4"	11°
30-1011	22" x 20" x 6"	16°
30-1012	22" x 20" x 8"	20°
30-1013	22" x 20" x 10"	25°
30-1014	26" x 24" x 6"	9°
30-1015	26" x 24" x 8"	13°
30-1016	26" x 24" x 10"	17°
30-1017	26" x 24" x 12"	22°



half rolls

Skillbuilders® Half Rolls

	L x W	Elevation	Height
30-1230	18" x 3"	0"	2"
30-1231	18" x 4"	1"	3"
30-1232	24" x 6"	1½"	4½"
30-1233	8" x 30"	2"	6"

Skillbuilders® Rolls

L x Dia

30-1000	24" x 4"
30-1001	24" x 6"
30-1002	24" x 8"
30-1003	36" x 8"
30-1004	36" x 10"
30-1005	36" x 12"
30-1006	48" x 12"
30-1007	48" x 14"
30-1008	48" x 16"

privacy screens

- Lightweight, foldable screens
- Available in 6-mil white (W), 4-mil yellow (Y), and 6-mil green (G) vinyl curtains; with or without casters



without
casters



with casters

vinyl color options

DEFAULT COLOR



white



yellow



green

without caster

		extended size
65-0100	2-panel	35"W x 68"H
65-0101	3-panel	48"W x 68"H
65-0102	4-panel	62"W x 68"H
65-0103	5-panel	84"W x 68"H

with casters

		extended size
65-0110	3-panel	48"W x 68"H
65-0111	4-panel	62"W x 68"H

original Dycem® non-slip rolls

Grips on both sides to prevent movement. Ideal for providing support, stability and confidence during everyday activities and exercise

- Effective anti-slip material grips on both sides without the need for adhesive
- Can be cut to any size
- Malleable so it can be fitted around various objects
- Yellow (Y), red (R), blue (B), black (BLK), white (W), forest green (G), lime (LIM), pink (PNK), silver (S)

8" width		16" width	
50-1503	8" x 16 yard roll	50-1508	16" x 16 yard roll
50-1500	8" x 10 yard roll	50-1505	16" x 10 yard roll
50-1518	8" x 5½ yard roll*	50-1519	16" x 5½ yard roll*
50-1501	8" x 6½ foot roll	50-1506	16" x 6½ foot roll
50-1502	8" x 3¼ foot roll	50-1507	16" x 3¼ foot roll
50-1509	roll rack dispenser**		

*available in yellow (Y), red (R) and blue (B) only

**fits up to 16" x 16 yards of any Dycem® roll



Dycem® non-slip table mats / pads

- Grips at angles up to 45 degrees
- Mats for healthcare, home and leisure provide a secure surface to anchor items such as cups, plates, tools, puzzles and note pads to table
- Extra stability during therapy or exercise can facilitate a variety of coordination and spatial awareness exercises
- Yellow (Y), red (R), blue (B), black (BLK), forest green (G), lime (LIM), pink (PNK), silver (S)



circular pad

50-1595	5½" diameter
50-1596	7½" diameter
50-1597	8½" diameter*
50-1598	10" diameter*

*8½" and 10" diameter circular pads available in yellow (Y), red (R) and blue (B) only



rectangular pad

50-1590	10" x 7¼"
50-1591	10" x 14"
50-1592	15" x 18"

Dycem® non-slip floor mat



- Use where there might be a risk of slipping
- Use for increased stability

50-1642	14" x 20", black
50-1641	18" x 24", gray (GRY) or blue (B)

Dycem® netting on roll

- Soft mesh non-slip material
- Improves grip and provides padding
- Can be easily cut to size

50-1550 24" x 6½ ft roll



Dycem® jar opener

- Effective on screw-top jars including child-resistant and tamper-proof lids
- Blue (B), yellow (Y), silver (S)

50-1650	jar opener
50-1651	bottle opener



Dycem® self-adhesive strips

- Adhesive backing permanently secures the strip
- Blue (B), yellow (Y)

50-1560 16" x 11/8" strip, set of 3



FabLife® hip kits

- Perfect for those recovering from hip, knee or back surgery

Kits include:

- reacher
- sock aid
- shoehorn
- sponge
- dressing stick



deluxe kit

standard kits

item	reacher length	sock aid	shoehorn	sponge
86-0070	26"	formed, 2 handles	plastic, 18"	contoured
86-0071	32"	formed, 2 handles	plastic, 18"	contoured
86-0072	26"	flexible, loop handle	metal, 24"	contoured
86-0073	32"	flexible, loop handle	metal, 24"	contoured

deluxe kits

item	reacher length	sock aid	shoehorn	sponge	dressing stick
86-0074	26"	formed, 2 handles	plastic, 24"	contoured	24"
86-0075	32"	formed, 2 handles	plastic, 18"	contoured	24"
86-0084	26"	formed, 2 handles	-	contoured	24"
86-0085	32"	formed, 2 handles	-	contoured	24"

FabLife® open-jaw reachers

- Durable metal frame
- Aids patient with limited reach and hand strength
- Reachers have magnetic tip to help retrieve objects
- Drag lug helps position products for easy pickup
- Open-jaw (closes when squeezed)

non-folding

- 50-1108 20" pistol grip
- A 50-1100 26", 6 oz pistol grip
- 50-1113 26", red pistol grip
- B 50-1101 32", 8 oz pistol grip

deluxe non-folding

- C 50-1110 20" full hand grip
- 50-1111 26" full hand grip
- 50-1112 32" full hand grip

folding

- 50-1102 26", folds to 14" pistol grip
- D 50-1106 32", folds to 16" pistol grip



FabLife® shoehorns

- Shoehorn is economical and convenient to remove shoes and socks



plastic shoehorn

- 86-0391 18"L
- 86-0391-25 18"L, 25 ea
- 86-0392 24"L
- 86-0392-25 24"L, 25 ea

metal shoehorn

- 86-0360 24"L
- 86-0360-25 24"L, 25 ea

FabLife® dressing stick

- Combination hook/pusher assists in putting on shirt and jacket, pulling up slacks or removing socks



- 86-0030 dressing stick, 24"L
- 86-0030-25 25 each
- 86-0031 deluxe dressing stick w/ foam grip
- 86-0031-25 25 each

FabLife® straight arm sponges

- Extends reach to wash lower back and extremities



- 45-2380 round sponge
- 45-2381 contoured sponge

FabLife® sock aids

- Designed to guide sock over heel
- **Formed sock aids** made of rigid, plastic material and have a foam patch to help keep sock in place
- **Flexible sock aids** are made of washable terry cloth and have a nylon inside to reduce friction



Formed Sock Aids

- 86-0041 2 handles
- 86-0040 continuous loop handle
- 86-0045 continuous loop handle, x-wide

Flexible Sock Aids

- 86-0001 sock aid with ridges, loop strap
- 86-0002 2 handles

FabLife® transfer boards



- Tapered boards for those with lower extremity disabilities
- 350 lb capacity

- solid wood**
- 50-3000 8 x 24" no cutouts
- 50-3001 8 x 30" no cutouts
- 50-3002 8 x 24" 1 cutout
- 50-3003 8 x 30" 1 cutout
- 50-3004 8 x 24" 2 cutouts
- 50-3005 8 x 30" 2 cutouts
- plastic**
- 50-3007 8 x 28" no cutouts

7-Day MediPlanner



- Organize pills by day and time of day (Morning, Noon, Evening, Bedtime)
- 28 total compartments
- See-through lids allow pills to be easily seen

- 85-0101 MediPlanner

Saunders® HomeTrac® cervical traction

- Directs traction force towards occiput
- Applies up to 50 lb of traction
- User has total control of force at all times
- Provides a comfortable and cost-effective option to continuing clinical traction treatments at home
- Includes user's guide and deluxe carrying case



00-7044 Saunders® cervical traction

Saunders® HomeTrac® lumbar traction

- Actively moving treatment surface ensures efficient application of traction
- Effectively applies up to 200 lb of traction
- Friction-free surface ensures smooth motion



50-1005 Saunders® lumbar traction

ComforTrac® home cervical traction

- Individual adjustable neck cushions slide in and out for custom fit
- Quick release allows user to immediately release all pressure; increasing safety and convenience
- Allows three angles of traction at 10, 15 and 20 degrees
- Comes with carrying case



50-1070 ComforTrac® cervical traction

ComforTrac® home lumbar traction

- Automatic carriage return increases ease of use, and aids in intermittent traction applications
- Effectively applies up to 200 lb of traction
- Quick release device on pump allows user to immediately release all pressure; increasing safety and convenience
- Hand pump displays both pounds in kilograms



50-1075 ComforTrac® lumbar traction

FabTrac® overdoor traction

- Mounts over door
- Economical for continued traction at home
- Includes head halter, 20 lb water bag and hardware
- Halter is machine washable



50-1000 1 ea
50-1000-25 25 ea

TracCollar® cervical traction

- Traction collar gently stretches and relaxes muscles for pain relief
- Portable, lightweight collar fastens easily with buckles
- For professional and home use



50-1095 small / medium neck, 14-16", red
50-1096 large / x-large neck, 16-18", blue

NeckPro II® overdoor traction

- Provides easy-to-use cervical traction without the need for water bags or heavy weights
- Includes NeckPro II® compression spring and ratchet assembly, head halter and overdoor bracket (fits doors measuring up to 1½" (3.8 cm) in thickness)



00-4395 NeckPro II®

HipTrac™ hip traction

- For hip osteoarthritis, FAI/labral pathology, joint stiffness, pain or muscle tightness
- Use at the clinic and at home to reduce pain and increase mobility to improve exercise compliance and performance
- Adjustable ankle binding



50-0920 HipTrac™

traction therapy

Saunders® cervical clinic traction

- Pulls at base of the occiput for cervical traction
- Performs horizontal traction
- Clevis and horizontal space bar included



00-7040K traction device

ComforTrac® cervical clinic traction

- Pulls at base of the occiput for cervical traction
- Performs horizontal traction
- Clevis and horizontal space bar included



50-1071 traction device

Galaxy TTET400, 4 section hi-lo traction table w/ foot bar lift

- 4 section traction table with scissor frame and 4 multi-directional wheels
- Head section has a removable face hole cover
- Adjustable height 22" to 40"
- 500 lb lifting capacity



35-2171 Galaxy TTET400, 4 section hi-lo traction table

TX traction unit

- Digital touch screen interface
- Intermittent, static and cyclic traction programs
- User-defined hold, rest, treatment times
- Traction tension parameters 0 - 200 lb
- Store up to 10 user protocols and up to 14 patient treatment sessions on patient data card

TX® traction unit

00-4759 TX® traction unit

TX® traction stand

00-3600 mobile traction stand



TX® with mobile traction stand

TX deluxe universal belt traction package

- For pelvic or lumbar traction application
- Has thoracic restraint, pelvic belt, 2 restraining straps, 2 cervical pillows and spreader bar



complete accessory package

00-1440 deluxe package

components

00-1433 HD pelvic package

00-1301 TX pillow (ea)

06-0030 17" spreader bar

sku index

sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg		
00-1002	52	00-1213	54	00-2782	61	00-4759	72	10-0205	10	10-0356	10	10-0634	20	10-0841	13	10-0883	13	10-0958	14	10-1163	37	10-1434	14
00-1004	52	00-1214	54	00-2783	61	00-6000	61	10-0206	10	10-0357	10	10-0638	20	10-0842	13	10-0884	13	10-0959	14	10-1164	37	10-1435	14
00-1006	52	00-1216	54	00-2786	61	00-7040	72	10-0207	10	10-0358	10	10-0690	12	10-0843	13	10-0885	13	10-0960	17	10-1178	8	10-1436	14
00-1008	52	00-1218	54	00-2787	61	00-7044	71	10-0208	10	10-0359	10	10-0710	9	10-0844	13	10-0886	13	10-0962	17	10-1179	8	10-1440	14
00-1010	52	00-1219	54	00-2788	61	01-4790	57	10-0209	10	10-0360	10	10-0711	9	10-0845	13	10-0900	14	10-0963	17	10-1385	8	10-1441	14
00-1012	52	00-1264	54	00-2789	61	01-8025	9	10-0210	10	10-0361	10	10-0712	9	10-0846	13	10-0901	14	10-0964	17	10-1400	14	10-1442	14
00-1014	52	00-1301	72	00-2791	61	01-8030	9	10-0211	10	10-0362	10	10-0717	9	10-0850	13	10-0902	14	10-0967	14	10-1401	14	10-1443	14
00-1026	52	00-1433	72	00-2793	61	06-0030	72	10-0212	10	10-0363	10	10-0718	9	10-0851	13	10-0903	14	10-0968	14	10-1402	14	10-1444	14
00-1030	52	00-1440	72	00-2795	61	07-0000	61	10-0213	10	10-0440	11	10-0740	15	10-0852	13	10-0905	14	10-0969	14	10-1403	14	10-1445	14
00-1031	52	00-1500	54	00-2796	61	07-0001	61	10-0214	10	10-0441	11	10-0741	15	10-0853	13	10-0906	14	10-0970	14	10-1404	14	10-1446	14
00-1032	52	00-1502	54	00-2797	61	07-0002	61	10-0215	10	10-0442	11	10-0742	15	10-0854	13	10-0907	14	10-0994	14	10-1405	14	10-1460	14
00-1033	52	00-1504	54	00-2802	52	07-0003	61	10-0216	10	10-0443	11	10-0743	15	10-0855	13	10-0908	14	10-0995	14	10-1406	14	10-1461	14
00-1062	52	00-1506	54	00-2871	61	07-0004	61	10-0217	10	10-0444	11	10-0744	15	10-0856	13	10-0910	14	10-1074	6	10-1410	14	10-1462	14
00-1064	52	00-1508	54	00-2872	61	07-0005	61	10-0218	10	10-0445	11	10-0745	15	10-0860	13	10-0911	14	10-1075	6	10-1411	14	10-1463	14
00-1102	52	00-1510	54	00-2873	61	07-0008	61	10-0219	10	10-0550	11	10-0746	15	10-0861	13	10-0912	14	10-1076	6	10-1412	14	10-1464	14
00-1104	52	00-1512	54	00-2874	61	07-7060	63	10-0250	10	10-0551	11	10-0747	15	10-0862	13	10-0913	14	10-1077	6	10-1413	14	10-1466	14
00-1108	52	00-1552	54	00-2876	61	07-7061	63	10-0251	10	10-0552	11	10-0749	15	10-0863	13	10-0918	14	10-1078	6	10-1414	14	10-1467	14
00-1110	52	00-1554	54	00-2878	61	07-7062	63	10-0252	10	10-0553	11	10-0750	15	10-0864	13	10-0919	14	10-1079	6	10-1415	14	10-1468	14
00-1118	52	00-1556	54	00-2893	61	07-7063	63	10-0253	10	10-0554	11	10-0751	15	10-0865	13	10-0920	14	10-1089	6	10-1416	14	10-1469	14
00-1120	52	00-1562	54	00-2894	61	07-7064	63	10-0254	10	10-0555	11	10-0752	15	10-0866	13	10-0921	14	10-1090	6	10-1420	14	10-1470	14
00-1124	52	00-1600	57	00-3102	54	07-7065	63	10-0255	10	10-0556	11	10-0753	15	10-0870	13	10-0923	14	10-1091	6	10-1421	14	10-1471	14
00-1126	52	00-2102	52	00-3600	72	07-7066	63	10-0256	10	10-0557	11	10-0756	17	10-0871	13	10-0924	14	10-1092	6	10-1422	14	10-1473	14
00-1132	52	00-2302	52	00-4002	52	07-7067	63	10-0257	10	10-0558	11	10-0757	17	10-0872	13	10-0925	14	10-1093	6	10-1423	14	10-1480	14
00-1200	54	00-2402	52	00-4004	52	09-0910	54	10-0258	10	10-0559	11	10-0758	15	10-0873	13	10-0926	14	10-1095	6	10-1424	14	10-1481	14
00-1202	54	00-2502	52	00-4016	52	10-0150	10	10-0350	10	10-0560	11	10-0759	15	10-0874	13	10-0940	14	10-1120	13	10-1425	14	10-1482	14
00-1204	54	00-2738	61	00-4228	52	10-0200	10	10-0351	10	10-0561	11	10-0775	17	10-0875	13	10-0941	14	10-1130	13	10-1426	14	10-1483	14
00-1206	54	00-2760	61	00-4238	58	10-0201	10	10-0352	10	10-0563	11	10-0784	17	10-0876	13	10-0943	14	10-1131	13	10-1430	14	10-1484	14
00-1208	54	00-2761	61	00-4248	58	10-0202	10	10-0353	10	10-0564	11	10-0785	17	10-0880	13	10-0944	14	10-1150	13	10-1431	14	10-1485	14
00-1210	54	00-2763	61	00-4266	58	10-0203	10	10-0354	10	10-0575	11	10-0786	17	10-0881	13	10-0956	14	10-1160	13	10-1432	14	10-1490	16
00-1212	54	00-2764	61	00-4395	71	10-0204	10	10-0355	10	10-0578	11	10-0787	17	10-0882	13	10-0957	14	10-1162	37	10-1433	14	10-1491	16

sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg		
10-1492	16	10-1896	16	10-2632	14	10-3382	10	10-5264	6	10-5564	6	10-5758	3	10-6176	5	10-6457	3	11-1531	54	12-0286	24	12-1205	29
10-1493	16	10-1904	20	10-2633	14	10-3383	10	10-5265	6	10-5565	6	10-5760	6	10-6177	5	10-6458	3	11-1532	54	12-0290	24	12-1210	29
10-1494	16	10-1910	6	10-2634	14	10-3384	10	10-5266	6	10-5566	6	10-5761	6	10-6178	5	10-6600	8	11-1533	54	12-0291	24	12-1211	29
10-1495	16	10-1911	6	10-2635	14	10-3738	15	10-5267	6	10-5569	6	10-5762	6	10-6210	5	10-6601	8	11-1534	54	12-0292	24	12-1212	29
10-1496	16	10-1912	6	10-2636	14	10-3760	15	10-5269	6	10-5610	2	10-5763	6	10-6211	5	10-6602	8	11-1600	51	12-0293	24	12-1220	29
10-1510	16	10-1913	6	10-2640	14	10-3761	15	10-5270	2	10-5611	2	10-5764	6	10-6212	5	10-6603	8	11-1601	51	12-0297	24	12-1221	29
10-1511	16	10-1914	6	10-2641	14	10-3762	15	10-5271	2	10-5612	2	10-5765	6	10-6213	5	10-6604	8	11-1603	51	12-0300	27	12-1222	29
10-1512	16	10-1915	6	10-2642	14	10-3763	15	10-5272	2	10-5613	2	10-5766	6	10-6214	5	10-6605	8	11-1604	51	12-0301	27	12-1360	32
10-1513	16	10-1970	17	10-2643	14	10-3764	15	10-5273	2	10-5614	2	10-5767	6	10-6215	5	10-6606	8	11-1605	51	12-0302	27	12-1361	32
10-1514	16	10-1971	17	10-2644	14	10-3765	15	10-5274	2	10-5615	2	10-5769	6	10-6216	5	10-6607	8	11-1606	51	12-0303	27	12-1362	32
10-1515	16	10-1972	17	10-2645	14	10-3766	15	10-5275	2	10-5616	2	10-5820	6	10-6217	5	10-6610	8	11-1607	51	12-0304	27	12-1370	33
10-1516	16	10-1973	17	10-2646	14	10-3767	15	10-5276	2	10-5617	2	10-5821	6	10-6218	5	10-6611	8	11-1631	55	12-0340	27	12-1390	33
10-1520	16	10-1975	17	10-2650	14	10-3768	15	10-5277	2	10-5618	2	10-5822	6	10-6220	5	10-6612	8	11-1633	55	12-0341	27	12-1391	33
10-1521	16	10-1980	17	10-2651	14	10-3769	15	10-5278	2	10-5620	2	10-5823	6	10-6221	5	10-6613	8	11-1634	55	12-0342	27	12-1392	33
10-1522	16	10-1981	17	10-2652	14	10-3770	15	10-5279	2	10-5621	2	10-5824	6	10-6222	5	10-7900	8	11-1635	55	12-0343	27	12-1393	33
10-1523	16	10-1982	17	10-2653	14	10-3771	15	10-5280	3	10-5622	2	10-5825	6	10-6223	5	11-0501	55	11-1636	55	12-0380	28	12-1395	32
10-1524	16	10-1983	17	10-2654	14	10-3772	15	10-5282	3	10-5623	2	10-5826	6	10-6224	5	11-0506	55	11-1637	55	12-0381	28	12-1396	32
10-1525	16	10-1985	17	10-2655	14	10-3773	15	10-5284	3	10-5624	2	10-5829	6	10-6225	5	11-0701	51	11-1638	55	12-0387	27	12-1397	32
10-1526	16	10-1990	17	10-2656	14	10-3774	15	10-5291	6	10-5625	2	10-5841	6	10-6226	5	11-0710	51	11-1710	51	12-0388	27	12-1398	32
10-1530	16	10-1991	17	10-2660	14	10-3775	15	10-5292	6	10-5626	2	10-5842	6	10-6227	5	11-0711	51	11-1711	51	12-0392	27	12-1406	33
10-1531	16	10-1992	17	10-2661	14	10-3776	15	10-5293	6	10-5627	2	10-5843	6	10-6228	5	11-0712	51	11-1712	51	12-0393	27	12-1407	33
10-1532	16	10-1993	17	10-2662	14	10-3777	15	10-5294	6	10-5628	2	10-5844	6	10-6229	5	11-0720	51	11-1713	51	12-0394	27	12-1408	33
10-1533	16	10-1995	17	10-2663	14	10-3778	15	10-5295	6	10-5629	2	10-5845	6	10-6230	5	11-0730	51	11-1714	51	12-0397	27	12-1440	33
10-1534	16	10-2270	17	10-2664	14	10-3779	15	10-5296	6	10-5630	2	10-5846	6	10-6231	5	11-0731	51	11-1715	51	12-0398	27	12-1441	33
10-1535	16	10-2271	17	10-2665	14	10-3785	15	10-5297	6	10-5631	2	10-5847	6	10-6232	5	11-0760	51	11-1718	51	12-0399	27	12-1442	33
10-1536	16	10-2272	17	10-2670	14	10-3786	15	10-5299	6	10-5632	2	10-5849	6	10-6233	5	11-0761	51	11-1720	51	12-0400	28	12-1443	33
10-1600	11	10-2273	17	10-2671	14	10-3787	15	10-5310	7	10-5633	2	10-5860	4	10-6234	5	11-0762	51	11-1722	51	12-0401	28	12-1444	33
10-1601	11	10-2275	17	10-2672	14	10-3797	15	10-5320	7	10-5634	2	10-5861	4	10-6235	5	11-0765	51	11-1730	51	12-0402	28	12-1445	34
10-1602	11	10-2280	17	10-2680	14	10-3798	15	10-5325	7	10-5635	2	10-5862	4	10-6236	5	11-1000	55	11-1744	51	12-0403	28	12-1446	34
10-1603	11	10-2281	17	10-2681	14	10-3799	15	10-5330	7	10-5636	2	10-5863	4	10-6237	5	11-1001	55	11-1746	51	12-0455	24	12-1450	34
10-1604	11	10-2282	17	10-2682	14	10-3820	15	10-5334	7	10-5637	2	10-5864	4	10-6238	5	11-1002	55	11-1748	51	12-0460	27	12-1455	34
10-1605	11	10-2283	17	10-2690	14	10-3821	15	10-5336	7	10-5638	2	10-5865	4	10-6260	5	11-1003	55	11-1750	51	12-0461	28	12-1460	34
10-1606	11	10-2285	17	10-2810	14	10-3840	15	10-5340	7	10-5639	2	10-5866	4	10-6261	5	11-1004	55	11-1752	51	12-0475	25	12-1461	34
10-1607	11	10-2290	17	10-2811	14	10-3841	15	10-5380	3	10-5650	2	10-5867	4	10-6262	5	11-1005	55	11-1760	51	12-0480	25	12-1462	34
10-1608	11	10-2291	17	10-2812	14	10-4000	16	10-5382	3	10-5651	2	10-5869	4	10-6263	5	11-1007	55	11-1961	53	12-0502	35	12-1464	34
10-1609	11	10-2292	17	10-2813	14	10-4001	16	10-5384	3	10-5652	2	10-5870	4	10-6264	5	11-1020	55	11-1962	53	12-0503	35	12-1465	34
10-1610	11	10-2293	17	10-2814	14	10-4002	16	10-5450	7	10-5653	2	10-5871	4	10-6265	5	11-1021	55	11-1964	53	12-0580	25	12-1466	34
10-1611	11	10-2295	17	10-2818	14	10-4003	16	10-5480	7	10-5654	2	10-5872	4	10-6266	5	11-1025	53	11-1965	53	12-0604	24	12-1467	34
10-1612	11	10-2500	10	10-2819	14	10-4004	16	10-5481	7	10-5655	2	10-5873	4	10-6267	5	11-1026	53	11-1966	53	12-0608	24	12-1468	34
10-1614	11	10-2501	10	10-2820	14	10-4005	16	10-5482	7	10-5656	2	10-5874	4	10-6268	5	11-1027	53	11-1981	53	12-1000	31	12-1469	34
10-1615	11	10-2502	10	10-2821	14	10-5104	7	10-5483	7	10-5657	2	10-5875	4	10-6270	5	11-1028	53	11-1982	53	12-1001	31	12-1470	34
10-1616	11	10-2503	10	10-2822	14	10-5105	7	10-5485	7	10-5658	2	10-5876	4	10-6271	5	11-1054	55	11-1984	53	12-1002	31	12-1471	34
10-1617	11	10-2504	10	10-2824	14	10-5106	7	10-5488	6	10-5660	3	10-5877	4	10-6272	5	11-1055	55	11-1985	53	12-1005	31	12-1472	34
10-1618	11	10-2505	10	10-2826	14	10-5111	7	10-5490	2	10-5682	3	10-5879	4	10-6273	5	11-1056	55	11-1986	53	12-1006	31	12-1475	34
10-1619	11	10-2506	10	10-2828	14	10-5180	3	10-5491	2	10-5684	3	10-5880	4	10-6274	5	11-1060	55	12-0070	22	12-1010	30	12-1476	34
10-1620	11	10-2507	10	10-2830	14	10-5181	3	10-5492	2	10-5687	3	10-5882	4	10-6275	5	11-1061	55	12-0071	22	12-1011	30	12-1477	34
10-1623	11	10-2508	10	10-2832	14	10-5182	3	10-5493	2	10-5688	3	10-5884	4	10-6276	5	11-1062	55	12-0072	22	12-1012	31	12-1480	33
10-1624	11	10-2509	10	10-2834	14	10-5183	3	10-5494	2	10-5689	3	10-6110	5	10-6277	5	11-1120	52	12-0080	22	12-1014	30	12-1481	33
10-1636	11	10-2510	10	10-2836	14	10-5184	3	10-5495	2	10-5691	2	10-6111	5	10-6278	5	11-1121	52	12-0081	22	12-1015	30	12-1490	33
10-1637	11	10-2511	10	10-2838	14	10-5185	3	10-5496	2	10-5692	2	10-6112	5	10-6301	4	11-1123	52	12-0082	22	12-1016	30	12-1492	33
10-1638	11	10-2512	10	10-2840	14	10-5186	3	10-5497	2	10-5693	2	10-6113	5	10-6302	4	11-1135	52	12-0090	22	12-1025	31	12-1495	33
10-1639	11	10-2513	10	10-2890	9	10-5187	3	10-5498	2	10-5694	2	10-6114	5	10-6303	4	11-1136	52	12-0091	22	12-1026	31	12-1500	34
10-1650	11	10-2514	10	10-2900	13	10-5188	3	10-5510	3	10-5695	2	10-6115	5	10-6304	4	11-1140	53	12-0092	22	12-1027	31	12-1501	34
10-1651	11	10-2515	10	10-2902	13	10-5191	2	10-															

sku index

sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg	sku	pg
12-1666	32	12-4521	39	13-2552	56	13-4801	62	15-5050	65	24-5260	44	24-5653	40	24-5835	39	24-7022	49	30-1998	49
12-1668	32	12-4522	39	13-2553	56	13-5220	59	15-5080	64	24-5261	44	24-5654	40	24-5842	39	24-7023	49	30-1999	49
12-1670	32	12-4523	39	13-2554	56	13-5221	59	15-5081	64	24-5270	44	24-5655	40	24-5843	39	24-7024	49	30-2100	19
12-1671	32	12-4524	39	13-2555	56	13-5223	59	15-5082	64	24-5271	44	24-5656	40	24-5844	39	24-7025	49	30-2101	19
12-1672	32	12-4525	39	13-2556	56	13-5224	59	15-5085	64	24-5272	44	24-5657	40	24-5845	39	24-7026	49	30-2110	19
12-1720	32	12-4526	39	13-2557	56	13-5250	59	15-5090	64	24-5273	44	24-5658	40	24-5846	39	24-7027	49	30-2111	19
12-1721	32	12-4527	39	13-2558	56	13-5251	59	15-5095	64	24-5274	44	24-5659	40	24-5847	39	24-7030	49	30-2200	19
12-1722	32	12-4528	39	13-2559	56	13-5252	59	15-5100	64	24-5280	44	24-5660	40	24-5850	43	24-7031	49	30-2201	19
12-1723	32	12-4556	39	13-2560	56	13-5253	59	15-5105	64	24-5281	44	24-5661	40	24-5851	43	24-7032	49	30-2210	19
12-1724	32	12-4558	39	13-2561	56	13-5310	51	15-5110	64	24-5282	44	24-5666	40	24-5852	43	24-7033	49	30-2211	19
12-1725	32	12-4609	39	13-2564	56	14-1050	49	15-5115	64	24-5283	44	24-5667	40	24-5853	43	24-7054	49	30-2280	19
12-1726	32	12-4613	39	13-2565	56	14-1051	49	15-5120	64	24-5284	44	24-5668	40	24-5854	43	24-7055	49	30-2281	19
12-1727	32	12-4614	39	13-2566	56	14-1371	49	15-5151	65	24-5300	44	24-5669	40	24-5855	43	24-7056	49	30-2290	19
12-1728	32	12-4619	39	13-2567	56	14-1410	49	15-5152	65	24-5301	44	24-5670	40	24-5870	42	24-7057	49	30-2291	19
12-1729	32	12-4620	39	13-2568	56	14-1411	49	15-5155	65	24-5310	44	24-5671	40	24-5871	42	24-7064	49	30-2340	19
12-1730	32	12-4628	39	13-2569	56	14-1412	49	15-5160	64	24-5311	44	24-5672	40	24-5872	42	24-7065	49	30-2341	19
12-1731	32	13-1287	59	13-2600	56	14-1430	50	16-1100	63	24-5320	44	24-5673	40	24-5873	42	24-7066	49	30-2400	21
12-1732	32	13-1288	59	13-2601	56	14-1431	50	16-1102	63	24-5321	44	24-5674	40	24-5874	42	24-7067	49	30-2401	21
12-1733	32	13-1291	59	13-2603	56	14-1432	50	16-1600	63	24-5330	44	24-5675	40	24-5875	42	24-9002	46	30-2500	21
12-1734	32	13-1293	59	13-2604	56	14-1433	50	16-1601	63	24-5331	44	24-5676	40	24-5880	43	24-9003	46	30-2502	21
12-1735	32	13-1308	59	13-2605	56	14-1434	50	16-1602	63	24-5340	44	24-5677	40	24-5881	43	24-9004	46	30-4000	18
12-1736	32	13-1323	59	13-2606	56	14-1435	50	16-1640	50	24-5341	44	24-5678	40	24-5900	42	24-9019	46	30-4001	18
12-1737	32	13-1470	58	13-2700	57	14-1436	50	16-1641	50	24-5342	44	24-5680	40	24-5901	42	24-9029	46	30-4002	18
12-1738	32	13-1471	58	13-2701	57	14-1440	50	16-1642	50	24-5343	44	24-5681	40	24-5902	42	24-9039	46	30-4003	18
12-1739	32	13-2410	56	13-2702	57	14-1441	50	24-3040	46	24-5350	44	24-5682	40	24-5903	42	24-9049	46	30-4010	18
12-1740	32	13-2413	56	13-2703	57	14-1442	50	24-3041	46	24-5351	44	24-5683	40	24-5904	42	24-9059	46	30-4011	18
12-1741	32	13-2414	56	13-2704	57	14-1443	50	24-3042	46	24-5352	44	24-5684	40	24-5905	42	24-9062	46	30-4012	18
12-1743	32	13-2417	56	13-2705	57	14-1444	50	24-3090	46	24-5353	44	24-5686	40	24-5910	43	24-9063	46	30-4013	18
12-1745	32	13-2420	56	13-2706	57	14-1445	50	24-3093	46	24-5400	44	24-5687	40	24-5911	43	24-9064	46	30-4050	19
12-1746	32	13-2422	56	13-2707	57	14-1450	50	24-3096	46	24-5401	44	24-5688	40	24-5912	43	30-1000	68	30-4051	19
12-1926	28	13-2425	56	13-2708	57	14-1452	50	24-3099	46	24-5410	44	24-5689	40	24-5914	43	30-1001	68	30-4052	19
12-1927	28	13-2428	56	13-2709	57	14-1460	50	24-3102	46	24-5411	44	24-5690	40	24-5915	43	30-1002	68	32-1230	21
12-1935	28	13-2430	56	13-2710	57	14-1470	50	24-3105	46	24-5420	44	24-5710	42	24-5916	43	30-1003	68	32-1232	21
12-1954	38	13-2433	56	13-2711	57	14-1471	50	24-3110	46	24-5421	44	24-5711	42	24-5917	43	30-1004	68	32-1234	21
12-1955	38	13-2435	56	13-2712	57	14-1472	50	24-3113	46	24-5430	44	24-5712	42	24-5918	43	30-1005	68	32-1236	21
12-1956	38	13-2438	56	13-2713	57	14-1473	50	24-3116	46	24-5431	44	24-5713	42	24-5919	43	30-1006	68	32-1238	21
12-2004	38	13-2439	56	13-2714	57	14-1474	50	24-3119	46	24-5440	44	24-5714	42	24-5920	43	30-1007	68	32-1240	21
12-2100	38	13-2441	56	13-2715	57	14-1501	48	24-3122	46	24-5441	44	24-5715	42	24-5921	43	30-1008	68	32-1247	21
12-2200	38	13-2442	56	13-2716	57	14-1503	48	24-3125	46	24-5450	44	24-5716	42	24-5922	43	30-1010	68	32-1248	21
12-2210	38	13-2445	56	13-2717	57	14-1520	48	24-3130	46	24-5451	44	24-5720	41	24-5923	43	30-1011	68	32-1257	21
12-2211	38	13-2447	56	13-2718	57	14-1522	48	24-3131	46	24-5500	45	24-5721	41	24-5924	43	30-1012	68	32-1258	21
12-2220	38	13-2448	56	13-2719	57	14-1630	50	24-3132	46	24-5501	45	24-5722	41	24-5925	43	30-1013	68	32-1500	21
12-2240	38	13-2450	56	13-2720	57	15-1171	8	24-3133	46	24-5502	45	24-5723	41	24-5926	43	30-1014	68	45-2380	70
12-2241	38	13-2452	56	13-2721	57	15-1700	66	24-3134	46	24-5503	45	24-5740	41	24-5927	43	30-1015	68	45-2381	70
12-2242	38	13-2453	56	13-2722	57	15-1702	66	24-4060	45	24-5504	45	24-5741	41	24-5928	43	30-1016	68	50-0920	71
12-2243	38	13-2455	56	13-2723	57	15-1711	66	24-4070	45	24-5505	45	24-5742	41	24-5929	43	30-1017	68	50-1000	71
12-2244	38	13-2457	56	13-2724	57	15-1736	66	24-4071	45	24-5510	45	24-5743	41	24-5930	43	30-1230	68	50-1005	71
12-2250	38	13-2459	56	13-2725	57	15-1737	66	24-4072	45	24-5511	45	24-5760	41	24-5931	43	30-1231	68	50-1010	8
12-2261	38	13-2460	56	13-2726	57	15-1754	67	24-4073	45	24-5512	45	24-5761	41	24-5932	43	30-1232	68	50-1023	8
12-2262	38	13-2462	56	13-2727	57	15-1758	67	24-4074	45	24-5513	45	24-5762	41	24-5933	43	30-1233	68	50-1024	8
12-3000	37	13-2465	56	13-2728	57	15-1816	66	24-4075	45	24-5514	45	24-5763	41	24-5934	43	30-1789	18	50-1025	8
12-3001	37	13-2468	56	13-3000	62	15-2062	67	24-4076	45	24-5515	45	24-5770	41	24-5940	41	30-1790	18	50-1034	8
12-3002	37	13-2470	56	13-3001	62	15-3730	48	24-4090	45	24-5520	45	24-5771	39	24-5941	41	30-1791	18	50-1035	8
12-3005	37	13-2471	56	13-3015	62	15-3731	48	24-4091	45	24-5521	45	24-5772	39	24-5942	41	30-1792	18	50-1036	8
12-3010	37	13-2472	56	13-3016	62	15-3732	48	24-4094	45	24-5522	45	24-5773	39	24-5945	42	30-1795	18	50-1037	8
12-3011	37	13-2473	56	13-3021	62	15-3733	48	24-4850	47	24-5523	45	24-5774	39	24-5946	42	30-1797	18	50-1070	71
12-3020	37	13-2474	56	13-3034	62	15-3740	48	24-4851	47	24-5524	45	24-5775	39	24-5947	42	30-1798	18	50-1075	71
12-3021	37	13-2475	56	13-3062	57	15-3742	48	24-4852	47	24-5525	45	24-5776	39	24-6300	50	30-1800	18	50-1095	71
12-3030	37	13-2476	56	13-3063	57	15-4230	67	24-4853	47	24-5530	45	24-5777	39	24-6301	50	30-1801	18	50-1096	71
12-3031	37	13-2477	56	13-3064	62	15-4231	68	24-4860	47	24-5531	45	24-5778	39	24-6302	50	30-1802	18	50-1100	70
12-3050	37	13-2478	56	13-3241	47	15-4232	68	24-4861	47	24-5540	45	24-5779	39	24-6303	50	30-1803	18	50-1101	70
12-3051	37	13-2479	56	13-3270	57	15-4233	68	24-4862	47	24-5541	45	24-5780	39	24-6304	50	30-1804	18	50-1102	70
12-3052	37	13-2480	56	13-3300	62	15-4235	68	24-4863	47	24-5550	45	24-5781	39	24-6305	50	30-1805	18	50-110	

product description / brand name index

item	pg	item	pg	item	pg	item	pg
3B charts and models	39	dynamometer, back-leg-chest	28	Hausmann	63	Orfit®	40 - 43
Absolute+Axis®	31	dynamometer, digital push/pull	27	HawkGrips®	47	Orfitube™	43
AcuAngle®	35	dynamometer, electronic hand	22	HeartSmart	38	Ossur®	57
adjustable hand grip exercisers	16	dynamometer, electronic hand smedley	24	heating unit	52 - 53	paraffin therapy	51
AFH	46, 50	dynamometer, electronic push	28	hip kits	70	paraffin wax, blocks and pellets	51
Airex®	21	dynamometer, electronic push/pull	28	HipTrac™	71	parallel bars	64, 67
Allen diagnostic tools	36	dynamometer, hydraulic hand	23	hook and loop	49	pedal exercisers	9
anatomical charts	37	dynamometer, hydraulic push/pull	27	Hot Button®	53	pedometers	38
anatomical models	37	dynamometer, hydraulic wrist	27	hot water bottle	53	pegboard	37
aneroid sphygmomanometers	38	dynamometer, mechanical push/pull	27	Hydrocollator®	52	Performance Attainment	35
Aquasonic®	58	dynamometer, smedley	24	IASTM	47	PhysioTrainer	9
Armedica treatment furniture	66 - 67	dynamometer, squeeze bulb	24	inclinometer	35 - 36	pill dispensers	70
Babinski hammers	34	Dynex™	24	inclinometer, digital	35	pinch gauge, electronic	25
balance pads	21	electric massagers	48 - 49	Intelect®	58, 61	pinch gauge, hydraulic	25
balance trainers	18 - 22	electrodes	59	Ionto+®	59	pinch gauge, mechanical	25
ball chairs	18	ErgoFET™	27 - 28	IontoPatch®	59	pinwheel	34
band and tubing handles and anchors	7	exercise balls	18	iontophoresis	59	Point Relief®	49, 51
band and tubing storage	7	exercise band	2 - 7	Jamar®	24	portable electrotherapy units	59 - 62
Baseline®	22 - 38	exercise band loops	6	Jumper® stability dome	19	portable ultrasound	58, 61
BIMS™	22	exercise band, latex-free	2 - 6	kettlebells	11	Power Hand Grip	16
body fat scale	38	exercise skate	13	kinesiology tape	47	privacy screens	68
Bongers®	49	exercise tubing	3 - 4	Lafayette®	28, 31, 37	Puff LITE®	14
Bosu®	19	exercise tubing with cuffs	6	large joint protractor	31	pulse oximeters	28
bottle warmers	59	exercise tubing with handles	6	lidocaine patches	51	putty, antimicrobial	14
braces	46	exercise tubing, latex-free	3 - 4	LidoSpot™	51	putty, LITE®	14
CanDo® 2 - 3, 6 - 21, 45, 47, 49, 67 - 68		FabLife®	70	magneciser exercisers	9	putty, standard	14
Chattanooga	52, 54, 57 - 59, 61	FabStim®	59	Magneciser®	9	putty, tools	14
chilling units	54 - 55	FabStretch®	8	Manosplint™	44 - 45	putty, variable strength	14
cold compression	56 - 57	FabTrac®	71	manual massagers	49 - 50	Puttycise®	14
Cold n' Hot®	54 - 55	finger circumference gauge	29	massage balls	49	reachers	70
cold packs	54 - 55	finger flexion gloves	16	massage furniture	48	rebounders	12
ColPac®	54	finger ladder	13	mat tables	65, 67	Relief Pak®	53 - 55
ComfortTrac®	71 - 72	finger pinch exerciser	13	mats	21	REP Band®	6
Comfy Splints™	46	Fiskars	45	measurement tapes	29	Richmar®	57, 60
Compass Health	58	fixed hand grip exercisers	16	medicine balls	12	Rickshaw™ rehab exerciser	12
Corbett	37	fluidotherapy	57	Mettler	57 - 58, 62	ring splint	50
core trainer dome	19	foam rollers	19	MicroFET™	24, 28	Saunders®	71 - 72
cuff weights	10	Free-Up®	47	mini cupping	50	shoehorns	70
cuff weights, adjustable	10	G5® percussors	48	moist heat packs	52 - 53	shortwave diathermy	57
Dellon®	33	Game Ready	56	moist heat pads	52	shoulder exerciser	9
dexterity test	37	gel ball exerciser	16	monofilaments, disposable	32	shoulder pulley	8
Digi-Extend n' Squeeze®	17	goniometers, digital	31	monofilaments, fold-up	32	Shuttle® MiniPress®	20
Digi-Extend®	17	goniometers, extendable	31	monofilaments, multi-filament	33	Sit 'n Reach® inclinometer	36
Digi-Flex LITE®	15	goniometers, finger	30	monofilaments, personal	32	Skillbuilders®	68
Digi-Flex Multi®	15	goniometers, metal	30	monofilaments, retracting	33	skinfold calipers	38
Digi-Flex Thumb®	15	goniometers, plastic	30 - 31	monofilaments, tactile	32 - 33	sock aids	70
Digi-Flex®	15	goniometers, x-ray	30	Multi-Axial platform system	20	SoftGrip®	10
Digi-Squeeze®	17	hand evaluation sets	26	multi-grip exerciser	8	Sonicator®	58, 62
discriminator	33	hand exercisers	13 - 17	MVP®	20	spinal assessment devices	35
dolorimeter	33	hand webs	13	NeckPro II®	63	splint pans	45
donut balls	18	hand weights	10 - 12	non-slip material	69	splinting accessories	45
dressing stick	70	HandMaster Plus™	17	Nylatex®	54	splinting material	40 - 45
dumbbells	11	hand-wrist exerciser	13	Orficast™	42	sponges	70
Dycem	69	Handy Grip™	11	Orflight®	41	Squid	56
						stethoscopes	38
						stim	59 - 62
						stim / ultrasound combo	60 - 62
						stools	63
						stopwatch and timer	38
						stretch straps	8
						stretching	8
						Sup-R Band®	4
						Sup-R Tubing®	4
						Sys*Stim®	62
						Taylor hammers	34
						the Cuff®	10
						Therabath®	51
						TheraPutty®	14
						TheraPutty® Plus	14
						therapy furniture	63 - 68
						Theratherm®	52
						Thermasonic	59
						tilt tables	65
						TOGU®	18 - 22
						topical analgesic	51
						Touch-Test®	33
						TracCollar®	71
						traction	71 - 72
						transfer boards	70
						treatment tables	64, 66
						TRI W-G furniture	64 - 65
						tuning forks	34
						Twist-n' Bend®	16
						TX traction	72
						ultrasound	58, 60 - 62
						ultrasound gel	58
						Uriel®	46
						Val-u-Band®	5
						vestibular discs / wedges	18
						vestibular wedges	19
						VitalStim®	59
						volumetric measuring	29
						WaTE™ balls	10
						WaTE™	10 - 11
						WaxWel®	51
						weight bars	11
						weight rack	12
						West	33
						Whitehall	45
						wooden staircase	67
						work bench	37
						wrist wheel exerciser	13
						Xtensor™	17



CAN DO[®] MGE[™] exerciser

- Easy-to-use grips allow for unlimited exercise options
- Each exerciser contains both large (6") and small (2") grips. Use small grips as anchor or for finger exercises
- Wrap around hands and feet to perform high impact upper and lower body exercises
- **Cloth exerciser contains no latex, no scent, and no powder coating**
- Can be machine washed and dried
- Easily change resistance by moving 1 grip up or down the exerciser, or by using a different resistance
- Progressively exercise with 8 different resistance levels
- Available in 9-grip (6 ft. exercisers) and 135-grip (30 yd. rolls) 15 exercisers



9 grips, 6 foot exerciser		9 grips, 6 foot exerciser (24 ea)	135 grips, 30 yard roll, 15 exercisers
tan	xx-light	10-6600	10-6610
yellow	x-light	10-6601	10-6611
red	light	10-6602	10-6612
green	medium	10-6603	10-6613
blue	heavy	10-6604	10-6614
black	x-heavy	10-6605	10-6615
silver	xx-heavy	10-6606	10-6616
gold	xxx-heavy	10-6607	10-6617



9 grip (6 foot)
exerciser



135 grips,
30 yard roll,
15 exercisers

Cut exerciser from roll
to any length (number
of grips) desired.



CAN DO[®] Band



- Available in the Preferred Colors[®] sequence
- **CanDo[®] Band** is made with Tear-Stop[®] technology to reduce punctures and tear "snap" breakage
- **Sup-R Band[®]** latex-free bands eliminate the need to stock both latex and latex-free bands

Twin-Pak[®] 100-yard box

Latex	Latex-Free
yellow x-light 10-5491	10-5651
red light 10-5492	10-5652
green medium 10-5493	10-5653
blue heavy 10-5494	10-5654
black x-heavy 10-5495	10-5655
yel red grn blu blk set (1 ea) 10-5498	10-5658

50-yard dispenser box

Latex	Latex-Free
yellow x-light 10-5221	10-5621
red light 10-5222	10-5622
green medium 10-5223	10-5623
blue heavy 10-5224	10-5624
black x-heavy 10-5225	10-5625
yel red grn blu blk set (1 ea) 10-5228	10-5628

6-yard dispenser box

Latex	Latex-Free
yellow x-light 10-5211	10-5611
red light 10-5212	10-5612
green medium 10-5213	10-5613
blue heavy 10-5214	10-5614
black x-heavy 10-5215	10-5615
yel red grn blu blk set (1 ea) 10-5218	10-5618

(30) 5-foot Singles[®] in dispenser box

Latex	Latex-Free
yellow x-light 10-5181	10-5741
red light 10-5182	10-5742
green medium 10-5183	10-5743
blue heavy 10-5184	10-5744
black x-heavy 10-5185	10-5745
yel red grn blu blk set (1 ea) 10-5188	10-5748



Twin-Pak[®] 100-yard box

Latex-Free
yellow x-light 10-6331
red light 10-6332
green medium 10-6333
blue heavy 10-6334
black x-heavy 10-6335
yel red grn blu blk set (1 ea) 10-6338

50-yard dispenser box

Latex-Free
yellow x-light 10-6321
red light 10-6322
green medium 10-6323
blue heavy 10-6324
black x-heavy 10-6325
yel red grn blu blk set (1 ea) 10-6328

6-yard dispenser box

Latex-Free
yellow x-light 10-6311
red light 10-6312
green medium 10-6313
blue heavy 10-6314
black x-heavy 10-6315
yel red grn blu blk set (1 ea) 10-6318

(30) 5-foot Singles[®] in dispenser box

Latex-Free
yellow x-light 10-6371
red light 10-6372
green medium 10-6373
blue heavy 10-6374
black x-heavy 10-6375
yel red grn blu blk set (1 ea) 10-6378



MedCOR Professionals
152 US Route One,
Scarborough, ME 04074
phone: 207-222-2828

Contact Info:

Jacob Currier 207-949-2215 jcurrier@medcorpro.com
Jason Finoia 207-844-4203 jfinonia@medcorpro.com
Kevin Melod 207-838-1172 kmelody@medcorpro.com